



LAWRENCE
LIVERMORE
NATIONAL
LABORATORY

Lipid Membrane Composition Analyzed by Multi-isotope Imaging Mass Spectrometry

C. G. Marxer, E. S. Park, P. K. Weber, I. D.
Hutcheon, S. G. Boxer

February 9, 2004

Biophysical Society Annual Meeting
Baltimore, MD, United States
February 14, 2004 through February 18, 2004

Disclaimer

This document was prepared as an account of work sponsored by an agency of the United States Government. Neither the United States Government nor the University of California nor any of their employees, makes any warranty, express or implied, or assumes any legal liability or responsibility for the accuracy, completeness, or usefulness of any information, apparatus, product, or process disclosed, or represents that its use would not infringe privately owned rights. Reference herein to any specific commercial product, process, or service by trade name, trademark, manufacturer, or otherwise, does not necessarily constitute or imply its endorsement, recommendation, or favoring by the United States Government or the University of California. The views and opinions of authors expressed herein do not necessarily state or reflect those of the United States Government or the University of California, and shall not be used for advertising or product endorsement purposes.

Lipid membrane composition analyzed by Multi-isotope Imaging Mass Spectrometry

The lateral organization of lipids and membrane-associated proteins in biological membranes is often detected by fluorescence microscopy. Although extremely sensitive, fluorescent labels, particularly those attached to lipid molecules, may alter their physical properties.

Multi-isotope imaging mass spectrometry (MIMS) using a NanoSIMS50 (Cameca) offers the opportunity to determine the lateral composition with very high lateral resolution (50nm), with component-specific information encoded by elemental and isotopic composition. Data will be presented on proteins deposited on SiO₂ (40 nm thick) by microcontact printing. Patterns of uniformly ¹⁵N labeled proteins can be readily distinguished from natural abundance (mostly ¹⁴N) proteins by detection of ¹²C¹⁵N⁻ vs. ¹²C¹⁴N⁻ fragments (see figure).

Microcontact printed proteins can be used to pattern supported bilayers [1]. Methods have been developed to freeze dry the bilayer largely retaining its lateral organization. Phospholipids can be imaged using CN⁻ or P⁻ and distinguished from protein barriers. By manipulating the isotopic composition of different lipid and/or membrane, spontaneous lateral organization or reorganization in response to an electric field parallel to the bilayer surface [1] can be probed.

[1] Acc. of Chem. Res., 35, 149 (2002)




Figure: *Fibronectin* (¹⁴N) and ¹⁵N-*Acp2* microcontact printed on SiO₂.

Authors:

Carine Galli Marxer, Chemistry Department, Stanford University, Stanford, CA 94305

Eun Sun Park, Chemistry Department, Stanford University, Stanford, CA 94305

Peter K. Weber, Lawrence Livermore National Laboratory, Livermore, CA 94551

Ian D. Hutcheon, Lawrence Livermore National Laboratory, Livermore, CA 94551

Steven G. Boxer, Chemistry Department, Stanford University, Stanford, CA 94305

This work was performed under the auspices of the U.S. Department of Energy by the University of California, Lawrence Livermore National Laboratory under Contract W-7405-Eng-48.