



LAWRENCE
LIVERMORE
NATIONAL
LABORATORY

Autonomous Pathogen Detection System

Benjamin Hindson, Mary McBride, Anthony Makarewicz, Bruce Henderer, Ujwal Sathyam, Shanavaz Nasarabadi, Kodumudi Venkateswaran, Dora Gutierrez, Sally Smith, Tom Metz, Bill Colston

February 20, 2004


2nd Joint Conference on Point Detection for Chemical and Biological Defense
Williamsburg, VA, United States
March 3, 2004 through March 3, 2004

Disclaimer

This document was prepared as an account of work sponsored by an agency of the United States Government. Neither the United States Government nor the University of California nor any of their employees, makes any warranty, express or implied, or assumes any legal liability or responsibility for the accuracy, completeness, or usefulness of any information, apparatus, product, or process disclosed, or represents that its use would not infringe privately owned rights. Reference herein to any specific commercial product, process, or service by trade name, trademark, manufacturer, or otherwise, does not necessarily constitute or imply its endorsement, recommendation, or favoring by the United States Government or the University of California. The views and opinions of authors expressed herein do not necessarily state or reflect those of the United States Government or the University of California, and shall not be used for advertising or product endorsement purposes.



Autonomous Pathogen Detection System



APDS team


- Assays
 - Mary McBride
 - Shanavaz Nasarabadi
 - Sally Smith
 - Venkat Venkateswaran
- Software, control, comm.
 - Bruce Henderer
 - Ujwal Sathyam
 - Dean Hadley
 - Robert Johnson
 - Paul Sargis
- Hardware
 - Ben Hindson
 - Tony Makarewicz
 - Bill Benett
 - Dora Gutierrez
 - Tom Metz
 - George Dougherty
- Project
 - John Dzenitis
 - Bill Colston

This work was performed under the auspices of the U. S. Department of Energy by the University of California, Lawrence Livermore National Laboratory under Contract No. W-7405-Eng-48.

Funding from DOE, DHS, and DoD.



BASIS and BioWatch: Centralized testing of air filters for biological agents



Autonomous Pathogen Detection System: Analysis at collector, networked reporting of results



- Aerosol collection
 - Up to 3,000 Lpm
 - Particle size selection
 - Samples are archived, can be cultured
- Sample preparation
 - Sequential injection analysis platform
 - Flexible and expandable
- Multiplexed detection and identification
 - Bead-based, Luminex™ immunoassay panels
 - bacteria, viruses, protein toxins
 - 11-plex + 4 controls
 - PCR confirmation of DNA sequences
 - Any antibody or sequence can be incorporated
- Data acquisition and control
 - Custom acquisition and analysis software
 - Wireless, Cellular, & Ethernet networking

Civilian & base protection differs from the battlefield

- Threat less known
 - Must test for more agents
- Operation is never-ending
 - Operating cost must be lower
- Different impact of alarms
 - Much less tolerance for false positives
- Treating victims vs. force protection
 - Some speed can be sacrificed for certainty




The APDS has advantages over the state-of-the-art


- APDS advantages vs. manual filter collection
 - Lower operating cost
 - Shorter time from collection to answer
 - Higher frequency reporting
- APDS advantages vs. triggered strip tests
 - Lower operating cost
 - Higher-quality antibody assays
 - PCR assay included
 - Assay upgrade/expansion is simple
- Issues for military applications
 - Ruggedization, equipment cost, speed vs. strip test




APDS gives lab-quality answers before a sample could get back to a lab



APDS software and communications allow remote command and control




High flow-rate aerosol collection



- 200 – 3,000 Lpm air sample in
- 4 mL liquid sample out
- Multistage
 - Prefractionator cap
 - Virtual impactor
 - Wetted-wall cyclone



Deeply multiplexed immunoassays with Luminex™



Data-rich signals from flow cytometer are monitored to detect and identify biological agents


- Fluorescent intensities, numbers of beads, statistics
- Multiplex signals have extra information
- Internal controls are important for confidence
 - Instrument control (detector OK)
 - Fluorescent control (label OK)
 - Antibody control (labeling antibody OK)
 - Negative control (no nonspecific binding)



11



Example of multiplex immunoassay signals




12



Orthogonal identification using PCR

- Uses DNA instead of protein recognition
 - Looking for different signature, so “orthogonal”
 - Tremendous amplification gives great sensitivity
- TaqMan used for confirmatory PCR




13



APDS automated PCR is shockingly repeatable

24 consecutive PCR runs & 1 negative control



Proven in chamber testing at Dugway Proving Grounds




- September 2002
 - Multiplex immunoassay
 - Live-agent releases
 - *B. anthracis* (anthrax)
 - *Y. pestis* (plague)
- September 2003
 - Multiplex immuno. + PCR
 - Killed-agent releases
 - *B. anthracis* (anthrax)
 - *Y. pestis* (plague)
 - Simulant releases
 - Botulinum toxoid, *B. globigii*

15




Identification and confirmation of a *Ba* release

Immunoassay




PCR




Identification and confirmation of a *Yp* release

Immunoassay




PCR




Identification and confirmation of a *Bg* release


Immunoassay



PCR



Identification of a botulinum toxoid release



Proven in fully autonomous testing



- Washington DC subway
 - June 2003
 - 1 unit, 7 days
- Albuquerque airport
 - December 2002
 - 2 units, 4 days
- Laboratory runs
 - 24×7 for 3 weeks
 - Many shorter runs
- Continuing tests in field

Current work

- Department of Homeland Security
 - Field operation
 - Commercialization
- Department of Defense (Tech. Transition Program)
 - Triggered by early-warning detector (BAWS)

