



LAWRENCE  
LIVERMORE  
NATIONAL  
LABORATORY

# Results of LLNL's Participation in the 16th OPCW Proficiency Test

H. R. Gregg, A. Alcaraz

January 27, 2005

## Disclaimer

---

This document was prepared as an account of work sponsored by an agency of the United States Government. Neither the United States Government nor the University of California nor any of their employees, makes any warranty, express or implied, or assumes any legal liability or responsibility for the accuracy, completeness, or usefulness of any information, apparatus, product, or process disclosed, or represents that its use would not infringe privately owned rights. Reference herein to any specific commercial product, process, or service by trade name, trademark, manufacturer, or otherwise, does not necessarily constitute or imply its endorsement, recommendation, or favoring by the United States Government or the University of California. The views and opinions of authors expressed herein do not necessarily state or reflect those of the United States Government or the University of California, and shall not be used for advertising or product endorsement purposes.

This work was performed under the auspices of the U.S. Department of Energy by University of California, Lawrence Livermore National Laboratory under Contract W-7405-Eng-48.

# Results of LLNL's Participation in the 16<sup>th</sup> OPCW Proficiency Test

Hugh R. Gregg<sup>1</sup>, Armando Alcaraz<sup>2</sup>

<sup>1</sup>Chemistry & Chemical Engineering, Chemistry & Material Sciences, LLNL

<sup>2</sup>Forensic Science Center, Chemistry & Material Sciences, LLNL

In the summer and fall of 2004, LLNL scientists prepared samples for the 16<sup>th</sup> official Organisation for the Prohibition of Chemical Weapons (OPCW) Proficiency Test.

Our report of this effort is attached.



**ORGANISATION FOR THE PROHIBITION  
OF CHEMICAL WEAPONS**

**Report of the  
Sixteenth Official OPCW  
Proficiency Test**

***Part I: Sample Preparation***

Laboratory code: Sample Prep Lab

Total number of pages: 128<sup>1</sup>

---

<sup>1</sup> Total number of pages including cover page and all attachments

|                          |
|--------------------------|
| <b>TABLE OF CONTENTS</b> |
|--------------------------|

|   | <b>Page no.</b> |
|---|-----------------|
| <b>Cover Page</b> .....                     | <b>1</b>        |
| <b>Table of Contents</b> .....              | <b>2</b>        |
| <b>Contributing Scientists</b> .....        | <b>3</b>        |
| <b>1 Introduction</b> .....                 | <b>4</b>        |
| <b>2 Test Scenario</b> .....                | <b>4</b>        |
| <b>3 List of Scheduled Chemicals</b> .....  | <b>4</b>        |
| <b>4 Sample Preparation</b> .....           | <b>6</b>        |
| <b>5 Packaging and Transportation</b> ..... | <b>9</b>        |
| <b>6 Analysis of Samples</b> .....          | <b>10</b>       |
| <b>Annexes</b>                              |                 |
| <b>A: Dispatch Announcement</b> .....       | <b>15</b>       |
| <b>B: Letter Accompanying Package</b> ..... | <b>16</b>       |
| <b>C: Dispatch Notification</b> .....       | <b>17</b>       |
| <b>D: Sample Receipt</b> .....              | <b>18</b>       |

## **Contributing Scientists from Lawrence Livermore National Laboratory**

### **Management, Reporting & Quality Assurance**

Armando Alcaraz, Principle Investigator

Hugh Gregg, co-PI, reporting

Tuijauna Mitchell-Hall, QA

### **Synthesis**

Phil Pagoria

Andy Vance

### **Sample Prep and Analytical Analysis**

Armando Alcaraz, sample prep, GC/MS analysis

Sarah Chinn, NMR

Hugh Gregg, GC/AED & LC/MS analysis

Carolyn Koester, sample prep, GC/MS analysis, stability studies

Robert Maxwell, NMR

### **Packaging and Shipping**

Marina Chiarappa-Zucca

Sarah Chinn

Peter Nunes

Phil Pagoria

Tom Tillotson

Rich Whipple

LLNL shipping staff

## 1 Introduction

The Sixteenth Official OPCW Proficiency Test started in October 2004. The samples were prepared by scientists affiliated with the Forensic Science Center at the Lawrence Livermore National Laboratory in Livermore, California, USA. The work was funded by the US Department of Energy.

The test scenario and the spiking and background chemicals were discussed and agreed in advance with the OPCW. The samples were prepared in accordance with “Work Instruction for the Preparation of Test Samples for OPCW Proficiency Tests” (Document No.: QDOC/LAB/WI/PT2).

The preparation of the test samples and their analysis are described in this report.

## 2 Test Scenario

A State Party has presented in accordance with Article IX of the Convention an inspection request for a challenge inspection. The Executive Council did not decide against this request, and the Director-General, in accordance with Part X of the Verification Annex, sent a challenge inspection team to a certain facility that has been accused of producing chemical warfare agents.

The facility is described by the inspection team as being a chemical research plant with well-equipped analytical and synthesis laboratories. The inspection team collected one sample inside the analytical laboratory from a container labeled “Organic Waste”. Another sample was collected outside a synthesis laboratory from a large unlabeled plastic container, “Liquid”. The inspection team collected a third sample, “Soil”, outside the facility’s chemical loading/unloading area. The inspection team did not analyze the samples on-site and forwarded them directly for off-site analysis. In addition, the inspection team could not collect corresponding blanks and used pure organic solvents and pre-cleaned soil similar to that found at the facility.










Please analyze the samples for the presence of **any Scheduled chemicals and/or their degradation/reaction products**, taking into account the characteristics of the samples.

## 3 List of Scheduled Chemicals

The majority of the chemicals used in this Proficiency Test were obtained commercially. Several spiked chemicals (both scheduled and unscheduled) were synthesized in-house. The purities of the scheduled chemicals were determined by NMR ( $^1\text{H}$  and  $^{13}\text{C}\{^1\text{H}\}$ ) and those results are presented in Part II, *Purity Checks*, of this report.

The Scheduled chemicals spiked into the three matrices are presented in the table on the next page.

|   |
|---|
| <b>List of Spiking Chemicals for the Sixteenth Official OPCW Proficiency Test</b> |
|---|

| Sample code | No | Compound name                          | CAS number | Spiking Chemical Structure  | Schedule number |
|-------------|----|--|------------|---|-----------------|
| O           | A  | Trichloronitromethane                  | 76-06-2    |    | 3.A.04          |
| O           | B  | 2-(N-Ethyl-N-propylamino)ethylchloride |            |    | 2.B.10          |
| O           | C  | 2-(N,N-Diisopropylamino)ethylchloride  | 96-79-7    |    | 2.B.10          |
| O           | D  | 2-(N,N-Diethylamino)ethylchloride      | 100-35-6   |    | 2.B.10          |
| L           | E  | Triethanolamine                        | 102-71-6   |    | 3.B.17          |
| L           | G  | 2-(N,N-Diisopropylamino)ethanol        | 96-80-0    |   | 2.B.11          |
| S           | H  | Ethyldiethanolamine                    | 139-87-7   |  | 3.B.15          |
| S           | I  | Methyldiethanolamine                   | 105-59-9   |  | 3.B.16          |
| S           | J  | Triethanolamine                        | 102-71-6   |  | 3.B.17          |



## 4 Sample Preparation

Several preliminary stability studies were carried out in the months preceding the 16<sup>th</sup> Official Proficiency Test. These studies were used to down-select the list of scheduled chemicals used for spiking and helped to determine stable amounts of materials and containers. Based on these studies, and in consultation with OPCW, three sample matrices were spiked as shown below.

### 4.1 Organic sample (O)

The following chemicals constituted the Organic sample.

| Component                                     | Origin/purity                   | Concentration (mg/L) |
|---|---------------------------------|----------------------|
| <b>Trichloronitromethane</b>                  | In-house synthesis              | 30                   |
| <b>2-(N-Ethyl-N-propylamino)ethylchloride</b> | In-house synthesis <sup>1</sup> | 10                   |
| <b>2-(N,N-Diisopropylamino)ethylchloride</b>  | Aldrich <sup>1</sup> , 97%      | 6                    |
| <b>2-(N,N-Diethylamino)ethylchloride</b>      | Aldrich <sup>1</sup> , 99%      | 30                   |
| Diesel fuel                                   | Local gas station <sup>1</sup>  | 200                  |
| Carbon tetrachloride                          | J.T. Baker, 100%                | 20                   |
| 2,4-Dinitro-t-butylbenzene                    | In-house synthesis, 97%         | 40                   |
| 2,4-Dinitrotoluene                            | Eastman Kodak                   | 20                   |
| 2-Nitroaniline                                | Eastman Kodak                   | 20                   |
| Hexane  | Fluka, min. 99.5%               | Solvent              |

<sup>1</sup>The preparation of the chlorides and diesel fuel cleanup are presented in section 4.4.

### 4.2 Liquid sample (L)

The following chemicals constituted the Liquid sample.

| Component                                | Origin/purity                            | Concentration (mg/L) |
|--|--|----------------------|
| <b>Triethanolamine</b>                   | Aldrich, 98%                             | 500                  |
| <b>2-(N,N-Diisopropylamino)ethanol</b>   | Aldrich, 99%                             | 30                   |
| 2-(N,N-Diethylamino)ethanol              | Aldrich, 98%                             | 20                   |
| 2-(N,N-Dimethylamino)ethanol             | Aldrich, 99%                             | 20                   |
| 2-(N-Butyl-N-methylamino)ethanol         | In house synthesis                       | 30                   |
| 2,6-Dinitroaniline                       | Aldrich, 97%                             | 20                   |
| 2-Nitroaniline                           | Eastman Kodak                            | 20                   |
| Pyridine                                 | Sigma Aldrich, 99.9%                     | 20                   |
| Poly(ethyleneglycol), M <sub>n</sub> 200 | Aldrich                                  | 200                  |
| Methanol                                 | Riedel-deHaën,<br>Chromasolv®, min 99.9% | Solvent              |

The following chemicals constituted the Soil sample.

| Component            | Origin/purity                  | Concentration (mg/kg) |
|----------------------|--------------------------------|-----------------------|
| Ethyldiethanolamine  | Aldrich, 98%                   | 20                    |
| Methyldiethanolamine | Aldrich, 99.9+%                | 20                    |
| Triethanolamine      | Aldrich, 98%                   | 20                    |
| Diesel fuel          | Local gas station <sup>1</sup> | 200                   |
| Caffeine             | Aldrich, 99%                   | 40                    |
| Nicotine             | Sigma, 99%                     | 40                    |
| Malation             | Chem Service, 99.2%            | 40                    |
| Tributylamine        | EM, 98%                        | 40                    |
| Tripentylamine       | EM, 98%                        | 40                    |
| 1,3-Dinitrobenzene   | Aldrich, 97%                   | 40                    |
| 2,4-Dinitrotoluene   | Eastman Kodak                  | 20                    |
| Dichloromethane      | Aldrich, 99.9%                 | (dried)               |
| Sand, purified       | J.T. Baker                     | Matrix                |

<sup>1</sup>The preparation of the chlorides and diesel fuel cleanup are presented in section 4.4.

#### 4.4 Chemical preparation

The three scheduled chlorine-containing amines were purchased or synthesized as their hydrochloride salts. These salts were converted to the chloride form by adding 1.0 mL of aqueous sodium hydrogen carbonate to a known amount of material and extracting the formed chloride into three serial aliquots of 1.0 mL of hexane. The resulting hexane extracts were dried using a column of anhydrous sodium sulfate. To estimate the resulting solution concentrations, it was assumed that 100% conversions and recoveries of the chlorides into the hexane were accomplished.

Diesel fuel was obtained from a local gas station and washed prior to use. To wash the fuel, 100 mL of fuel was placed in a 500 mL round-bottomed flask and 100 mL of concentrated sulfuric acid was added. The mixture was stirred vigorously for 2 hours. The two phases were placed in a separatory funnel and the sulfuric acid layer was discarded. The fuel layer was washed sequentially with 20 mL water, 10 mL of a 5% solution of sodium hydrogen carbonate, and 20 mL of water. The washed fuel was dried over anhydrous sodium sulfate and finally dried over anhydrous sodium carbonate, yielding a clear, colorless liquid.

The remaining chemicals were obtained or synthesized neat and ready to use.

#### 4.5 Protocols for preparing test samples

The soil samples and blank soil samples were prepared on October 6, 2004 and the two liquid samples and their corresponding blanks were prepared on October 7, 2004. Two staff members of the OPCW laboratory witnessed the sample preparation process. A total of 40 sets of samples and their corresponding blanks were made. Seventeen sets of samples and their corresponding blanks

were sent to laboratories that were planning to participate in the 16<sup>th</sup> OPCW Proficiency Test, two sets of samples and blanks were sent to the evaluation laboratory, and two sets of samples and blanks were sent to the OPCW laboratory. Three sets of samples and blanks, chosen randomly by the OPCW staff members, were used by the preparation laboratory for stability testing. An additional set was used by the preparation laboratory for qualitative analysis. The remaining sample sets were retained by the sample preparation laboratory for contingencies.

#### **4.6 Organic sample (Code O) and Organic sample blank (Code OB) preparation**

Hexane was added to a 500 mL, glass, solvent-rinsed, volumetric flask. Next, 125 microliters of diesel fuel was added to the hexane. Stock solutions of individual test chemicals were made in hexane and ranged in concentration from 1-30 mg/mL. A sufficient volume of each stock solution was added to the hexane solution to yield the desired analyte concentrations, listed in the previous section. The hexane solution was brought to a final volume of 500 mL and analyzed by GC/MS to verify that the solution was of the expected composition.

Neat hexane was used as the organic solvent blank. The hexane was analyzed by GC/MS to verify that it contained none of the added test chemicals.

#### **4.7 Liquid sample (Code L) and Liquid sample blank (Code LB) preparation**

Methanol was added to a glass, solvent-rinsed, 2.5 L jug. Next, 250 microliters of poly(ethyleneglycol) was added to the methanol. Stock solutions of individual test chemicals were made in methanol and ranged in concentration from 5-50 mg/mL. A sufficient volume of each stock solution was added to the methanol solution to yield the desired analyte concentrations, listed in the previous section. The methanol solution was brought to a final volume of 1.0 L and analyzed by GC/MS to verify that the solution was of the expected composition.

Neat methanol was used as the organic solvent blank. The methanol was analyzed by GC/MS to verify that it contained none of the added test chemicals.

#### **4.8 Soil sample (Code S) and Soil sample blank (Code SB) preparation**

To make the soil samples, the test chemicals were first added to dichloromethane. The resulting dichloromethane solution was added to the sand and allowed to evaporate, leaving the test chemicals on the soil.

Dichloromethane was added into a glass, solvent-rinsed, 2.5 L jug. Next, 625 microliters of diesel fuel was added to the dichloromethane. Stock solutions of individual test chemicals were made in dichloromethane and ranged in concentration from 10-20 mg/mL. A sufficient volume of each stock solution was added to the dichloromethane solution to yield the desired analyte concentrations, listed in the previous section. The dichloromethane solution was brought to a final volume of 1.25 L and analyzed by GC/MS to verify that the solution was of the expected composition.

Into each of five, solvent-rinsed, 1-L beakers were placed 500 g sand. Into each beaker containing sand, 250 mL of the above dichloromethane solution were added. The resulting mixture was stirred with a metal spatula and placed in a chemical fume hood. The dichloromethane was allowed to evaporate overnight, leaving the test chemicals on the sand. When the sand was dry, it was again

mixed in the beaker and transferred into an empty bottle, which had previously contained clean sand. The sand contained in all of the beakers (2500 g total) was transferred to the bottle. The bottle containing spiked sand was then placed on a roller and homogenized for approximately 5 minutes. The resulting spiked sand was then sampled, extracted, derivitized and analyzed by GC/MS to verify that it was of the expected composition.

Blank samples were prepared by exposing 2500g of sand to 1.25 L of dichloromethane. The blank sand was prepared in a manner identical to that of the spiked soil, with the exception that no test chemicals were added. The resulting blank sand was sampled, extracted, derivitized and analyzed by GC/MS to verify that it was free of test chemicals.

#### **4.9 Standards**

Two sets of standards were prepared, one each for the OPCW laboratory and the evaluation laboratory. These standards were prepared by dilution from the stock solutions in the range of 1 to 2 mg/mL. A set of vials of each standard, along with a description, was shipped a few days after the samples were shipped.

### **5 Packaging and Transportation**

The sample preparation laboratory received the list of participating laboratories on 28 September 2004 and contacted them on several occasions to confirm delivery addresses and preferred shippers. It was gratifying that most laboratories responded promptly, however there were significant communication issues with several laboratories.

#### **5.1 Announcement of Sample Dispatch**

Annex A contains a copy of the fax and/or e-mail detailing the sample dispatch information.

#### **5.2 Packaging**

10 mL of the organic solution (hexane) and corresponding blank were pipetted into 20-mL vials (pre-cleaned borosilicate glass) and capped (Teflon-lined, silicon-backed seal). The caps were tightened, taped with a thin strip of duct tape followed by the application of a tamper indicating seal. These bottles were labeled and placed inside of plastic bottles with absorbent material. These secondary bottles were sealed with a heat-shrink plastic seal.

18 mL of the liquid solution (methanol) and corresponding blank were pipetted into 20-mL vials (pre-cleaned borosilicate glass) and capped (Teflon-lined, silicon-backed seal). The caps were tightened, taped with a thin strip of duct tape followed by the application of a tamper indicating seal. These bottles were labeled and placed inside of plastic bottles with absorbent material. These secondary bottles were sealed with a heat-shrink plastic seal.

50 grams of the soil and the soil blank were weighed into 2-oz. clear glass jar and sealed with caps (Teflon-bonded to polypropylene foam). The caps were tightened, taped with a thin strip of duct tape followed by the application of a tamper indicating seal. These bottles were labeled and sealed inside of a heat-sealable plastic bag.

All six containers (O, OB, L, LB, S, SB) were placed inside of a foam-lined cardboard box. Each box was accompanied by a letter (Annex B). The shipping information was affixed to the exterior of the box.

### 5.3 Confirmation of Dispatch

The LLNL shipping department took custody of the samples and provided tracking numbers. Participating laboratories were notified of the carrier and tracking number (see Annex C for an example). Due to multiple transportation restrictions, multiple carriers had to be used.

### 5.4 Sample Receipt

At least 10 of the 17 laboratories received the samples within one week of shipping. A summary of the shipping and receipt information is provided in Annex D.

## 6 Analysis of Samples

The three sample matrices were analyzed both quantitatively and qualitatively.

### 6.1 Quantitative Analysis

Each of the three matrices required different sample work up procedures.

#### 6.1.1 Analytes in Organic Solution

Analysis of spiked chemicals in the organic solution was performed by direct injection of one microliter of sample solution into the GC/MS (EI). Quantitation was performed using extracted ion chromatograms and using heptadecane (present as a component of the diesel fuel background) as an internal standard. A three point calibration curve (10 to 40 mg/L) was prepared fresh for each sampling day.

**m/z EICs used for quantitation**

| Cmpd # | Compound name                          | EIC m/z |
|--------|--|---------|
| A      | Trichloronitromethane                  | 117     |
| B      | 2-(N-Ethyl-N-propylamino)ethylchloride | 120     |
| C      | 2-(N,N-Diisopropylamino)ethylchloride  | 106     |
| D      | 2-(N,N-Diethylamino)ethylchloride      | 86      |

#### 6.1.2 Analytes in Liquid Solution

Analysis of spiked chemicals in the organic solution was performed by direct injection of one microliter of sample solution (compound E was diluted with methanol before injection) into the GC/MS (EI). Quantitation was performed using extracted ion chromatograms and using PEG200 (present as a component added to the background) as an internal standard. A three point calibration

curve (for compound E: 50 to 200 mg/L; for compound G: 10 to 40 mg/L) was prepared fresh for each sampling day.

**m/z EICs used for quantitation**

| Cmpd # | Compound name                         | EIC m/z |
|--------|---------------------------------------|---------|
| E      | Triethanolamine                       | 118     |
| G      | 2-(N,N-Diisopropylamino)ethanol amine | 114     |

### 6.1.3 Analytes in Soil Sample

Analysis of spiked chemicals in the soil sample was performed by extraction, derivatization, and direct injection of one microliter of sample solution into the GC/MS (EI). Quantitation was performed using extracted ion chromatograms using heptadecane (present as a component of the diesel fuel background) as an internal standard. A three point calibration curve (10 to 40 mg/L) was prepared fresh for each sampling day and was derivatized in the same manner and at the same time as the sample extracts.

The extraction of the sand was accomplished by sonicating 2.0 g of sand with 2.0 mL of methanol for 15 minutes. The resulting solution was centrifuged for 3-5 minutes and 500 µL of the clarified extract were transferred to a 4-mL vial. The extract was evaporated to dryness with a gentle stream of nitrogen and reconstituted immediately with 200 µL of N, O-Bis(trimethylsilyl)trifluoroacetamide (BSTFA). The resulting solution was sealed and heated at 60°C for 30 minutes to allow derivatization of the three amines.

**m/z EICs used for quantitation**

| Cmpd # | Compound name        | EIC m/z |
|--------|----------------------|---------|
| H      | Ethyldiethanolamine  | 174     |
| I      | Methyldiethanolamine | 160     |
| J      | Triethanolamine      | 262     |

### 6.1.4 Quantitative Analysis Results

**Chemical A: Trichloronitromethane (mg/L)**

| Sample         | Day 0 | Day 7 | Day 11 | Day 15 | Day 21 | Day 29 | Day 40 |
|----------------|-------|-------|--------|--------|--------|--------|--------|
| O/25           | 30    | 25    | 26     | 32     | 26     | 26     | 26     |
| O/35           | 33    | 30    | 24     | 30     | 26     | 25     | 27     |
| O/40           | 33    | 28    | 29     | 30     | 26     | 26     | 31     |
|                |       |       |        |        |        |        |        |
| <b>Average</b> | 32    | 28    | 26     | 31     | 26     | 26     | 28     |
| <b>%RSD</b>    | 5.4%  | 9.1%  | 9.6%   | 3.8%   | 1.0%   | 0.9%   | 8.0%   |

**Chemical B: 2-(N-Ethyl-N-propylamino)ethylchloride (mg/L)**

| Sample  | Day 0 | Day 7 | Day 11 | Day 15 | Day 21 | Day 29 | Day 40 |
|---------|-------|-------|--------|--------|--------|--------|--------|
| O/25    | 10.4  | 10.6  | 10.7   | 11.0   | 10.4   | 10.7   | 10.2   |
| O/35    | 10.9  | 11.5  | 10.2   | 10.4   | 10.9   | 9.5    | 10.7   |
| O/40    | 11.3  | 11.6  | 11.4   | 10.7   | 11.4   | 10.1   | 11.4   |
|         |       |       |        |        |        |        |        |
| Average | 10.9  | 11.2  | 10.8   | 10.7   | 10.9   | 10.1   | 10.8   |
| %RSD    | 4.2%  | 4.7%  | 5.4%   | 2.9%   | 4.5%   | 6.0%   | 5.7%   |

**Chemical C: 2-(N,N-Diisopropylamino)ethylchloride (mg/L)**

| Sample  | Day 0 | Day 7 | Day 11 | Day 15 | Day 21 | Day 29 | Day 40 |
|---------|-------|-------|--------|--------|--------|--------|--------|
| O/25    | 5.3   | 5.2   | 5.4    | 5.2    | 4.8    | 5.8    | 5.2    |
| O/35    | 6.0   | 5.6   | 5.2    | 5.0    | 5.1    | 5.1    | 5.4    |
| O/40    | 6.5   | 5.7   | 5.8    | 5.2    | 5.2    | 5.5    | 5.7    |
|         |       |       |        |        |        |        |        |
| Average | 5.9   | 5.5   | 5.5    | 5.1    | 5.1    | 5.5    | 5.4    |
| %RSD    | 9.9%  | 4.2%  | 5.3%   | 2.1%   | 3.9%   | 5.8%   | 4.8%   |

**Chemical D: 2-(N,N-Diethylamino)ethylchloride (mg/L)**

| Sample  | Day 0 | Day 7 | Day 11 | Day 15 | Day 21 | Day 29 | Day 40 |
|---------|-------|-------|--------|--------|--------|--------|--------|
| O/25    | 24    | 23    | 23     | 25     | 25     | 23     | 24     |
| O/35    | 25    | 25    | 22     | 24     | 26     | 21     | 25     |
| O/40    | 25    | 25    | 24     | 24     | 27     | 22     | 26     |
|         |       |       |        |        |        |        |        |
| Average | 25    | 24    | 23     | 24     | 26     | 22     | 25     |
| %RSD    | 3.8%  | 5.9%  | 4.9%   | 3.8%   | 4.6%   | 4.7%   | 4.6%   |

**Chemical E: Triethanolamine (mg/L)**

| Sample  | Day 0 | Day 7 | Day 11 | Day 15 | Day 21 | Day 29 | Day 40 |
|---------|-------|-------|--------|--------|--------|--------|--------|
| O/25    | 433   | 374   | 384    | 418    | 426    | 516    | 500    |
| O/35    | 487   | 490   | 477    | 409    | 476    | 564    | 557    |
| O/40    | 530   | 516   | 461    | 482    | 506    | 512    | 579    |
|         |       |       |        |        |        |        |        |
| Average | 483   | 460   | 441    | 436    | 469    | 531    | 545    |
| %RSD    | 10%   | 16%   | 11%    | 9.0%   | 8.6%   | 5.5%   | 7.5%   |

**Chemical G: 2-(N,N-Diisopropylamino)ethanol (mg/L)**

| Sample  | Day 0 | Day 7 | Day 11 | Day 15 | Day 21 | Day 29 | Day 40 |
|---------|-------|-------|--------|--------|--------|--------|--------|
| O/25    | 30    | 28    | 29     | 33     | 33     | 28     | 26     |
| O/35    | 32    | 29    | 25     | 32     | 31     | 26     | 30     |
| O/40    | 33    | 30    | 27     | 31     | 31     | 28     | 29     |
|         |       |       |        |        |        |        |        |
| Average | 32    | 29    | 27     | 32     | 32     | 28     | 28     |
| %RSD    | 4.6%  | 3.2%  | 7.6%   | 4.0%   | 4.3%   | 4.5%   | 7.7%   |

**Chemical H: Ethyldiethanolamine (mg/kg)**

| Sample  | Day 0 | Day 8 | Day 12 | Day 16 | Day 22 | Day 30 | Day 40 |
|---------|-------|-------|--------|--------|--------|--------|--------|
| O/25    | 10.3  | 8.6   | 7.5    | 10.9   | 7.9    | 10.0   | 9.3    |
| O/35    | 11.6  | 9.6   | 10.9   | 9.7    | 8.1    | 10.0   | 10.2   |
| O/40    | 10.2  | 10.3  | 10.8   | 11.8   | 7.3    | 10.4   | 10.8   |
|         |       |       |        |        |        |        |        |
| Average | 10.7  | 9.5   | 9.7    | 10.8   | 7.8    | 10.2   | 10.1   |
| %RSD    | 7.2%  | 8.8%  | 20%    | 9.7%   | 5.2%   | 2.3%   | 7.6%   |

**Chemical I: Methyldiethanolamine (mg/kg)**

| Sample  | Day 0 | Day 8 | Day 12 | Day 16 | Day 22 | Day 30 | Day 40 |
|---------|-------|-------|--------|--------|--------|--------|--------|
| O/25    | 7.5   | 7.6   | 7.4    | 6.6    | 5.3    | 7.0    | 7.3    |
| O/35    | 9.3   | 7.9   | 8.2    | 6.2    | 5.5    | 6.8    | 7.7    |
| O/40    | 8.3   | 7.7   | 7.7    | 7.8    | 4.5    | 7.2    | 7.7    |
|         |       |       |        |        |        |        |        |
| Average | 8.4   | 7.7   | 7.8    | 6.9    | 5.1    | 7.0    | 7.6    |
| %RSD    | 11%   | 2.1%  | 5.4%   | 12%    | 10%    | 3.1%   | 3.0%   |

**Chemical J: Triethanolamine (mg/kg)**

| Sample  | Day 0 | Day 8 | Day 12 | Day 16 | Day 22 | Day 30 | Day 40 |
|---------|-------|-------|--------|--------|--------|--------|--------|
| O/25    | 8.7   | 10.4  | 9.2    | 11.3   | 8.1    | 12.5   | 12.7   |
| O/35    | 9.6   | 10.9  | 9.5    | 8.7    | 8.0    | 12.3   | 13.0   |
| O/40    | 6.6   | 13.0  | 8.9    | 12.1   | 6.9    | 12.9   | 12.8   |
|         |       |       |        |        |        |        |        |
| Average | 8.3   | 11.4  | 9.2    | 10.7   | 7.7    | 12.6   | 12.8   |
| %RSD    | 19%   | 12%   | 3.3%   | 17%    | 8.2%   | 2.2%   | 1.1%   |



## 6.2 Qualitative Analysis

The results of the qualitative analysis are presented in Part III, *Qualitative Analysis*, of this report. That section of the report is documented using the standard forms for reporting results of an Official OPCW Proficiency Test.

The Organic samples were analyzed directly by GC/MS EI & CI. The Liquid samples were analyzed after derivitization with BSTFA to form the trimethylsilyl derivative and then analyzed by GC/MS EI & CI. The Soil samples were first extracted with methanol, derivitized with BSTFA to form the trimethylsilyl derivative and then analyzed by GC/MS EI & CI.

Detailed procedures are described in Part III of this report.

**Annex A: Dispatch Announcement**

October 6, 2004

University of California  
Lawrence Livermore National Laboratory  
7000 East Ave  
Livermore, CA 94550

**Date of Sample Dispatch for the Official 16<sup>th</sup> OPCW Proficiency Test**

Dear Colleagues,

The samples for the 16<sup>th</sup> Official OPCW Proficiency Test will be dispatched from our laboratory on **8 October 2004**. FedEx or other air freight carrier will deliver the packages. The airway bill numbers will be released on the day of sample dispatch. The test samples contain an organic sample, liquid sample, and soil with their corresponding blanks.

Please acknowledge the receipt of the fax or email to:

Fax: (925) 423-9014  
Email: [alcaraz1@llnl.gov](mailto:alcaraz1@llnl.gov)

Good luck with the analysis.

*Armando Alcaraz*

---

Sincerely,  
Armando Alcaraz  
Forensic Science Center  
Lawrence Livermore National Laboratory  
Voice: (925) 423-6889  
Fax: (925) 423-9014

**Annex B: Letter Accompanying Package**

October 8, 2004

University of California  
Lawrence Livermore National Laboratory  
7000 East Ave  
Livermore, CA 94550

Participating Laboratories  
16<sup>th</sup> Official OPCW Proficiency Test

**TEST SAMPLES FOR THE 16<sup>TH</sup> OFFICIAL OPCW PROFICIENCY TEST**

Dear Colleagues,

1. Please find enclosed samples for chemical analysis according to the information which you received from the OPCW regarding the 16<sup>th</sup> OPCW test scenario.
2. The package contains four vials sealed in plastic containers and two jars sealed in plastic wrap. The samples are an organic solvent, liquid, and sand with their corresponding blanks. They are labeled:

**Samples:**

The three samples coded are as follows:

- O: Sample from container "Organic Waste"
- L: Sample appearing to be an organic liquid
- S: Soil from the facility's chemical loading/unloading area

with their corresponding blanks

- OB
- LB
- SB

Your laboratory code is a two-digit number, which has been selected at random.

3. Please confirm the arrival date of the samples at your laboratory and their condition, by fax or email to:

Mr. Stefan Mogl, OPCW Laboratory (Fax:+31 15 2840679, email:  
OPCWrijj@worldonline.nl)

Regards,

---

Armando Alcaraz

**Annex C: Dispatch Notification**

The following is an example of e-mail or fax notification.

Dear Prof. Col József Fűrész and Lt Gellért Karvaly,

Below is the web address and FedEx tracking number for your OPCW test samples:

[http://www.fedex.com/cgi-bin/tracking?action=track&language=english&cntry\\_code=us&initial=x&tracknumbers=40090164874](http://www.fedex.com/cgi-bin/tracking?action=track&language=english&cntry_code=us&initial=x&tracknumbers=40090164874)

**FedEx Track numbers = 40090164874**

**Please acknowledge the receipt of the fax or email to:**

Fax: (925) 423-9014

Email: [alcaraz1@llnl.gov](mailto:alcaraz1@llnl.gov)

Best regards,  
Armando Alcaraz  
Program Element Leader  
Lawrence Livermore National Laboratory  
Forensic Science Center, L-178  
7000 East Ave  
Livermore, CA 94550-9234  
(925) 423-6889 phone  
(925) 423-9014 fax

**Annex D: Sample Receipt**

| <b>Country</b>        | <b>Name of laboratory and address</b>  | <b>Air Freight Carrier</b> | <b>Airway bill number</b> | <b>Receipt package date</b> |
|-----------------------|--|----------------------------|---------------------------|-----------------------------|
| <b>Austria</b>        | RD-ARWT/ABCUT<br>Vorgartenstr. 225<br>A-1024 Wien<br>Austria   | FedEx                      | 40090164815               | 13/10/2004                  |
| <b>Belgium</b>        | DLD (Departement Laboratoria van Defensie)<br>Kwartier Majoor Housiau<br>Martelarenstraat 181<br>B-1800 Vilvoorde (Peutie)<br>Belgium  | FedEx                      | 40090164841               | 12/10/2004                  |
| <b>Czech Republic</b> | Research Institute for Organic Syntheses<br>Rybitvi 296<br>532 18 Pardubice 20<br>Czech Republic   | FedEx                      | 40090164852               | 12/10/2004                  |
| <b>Denmark</b>        | Ministry of Defence,<br>Emergency Management Agency, Chemical Division<br>Universitetsparken 2,<br>Nørre Alle 67, 7 <sup>th</sup> floor<br>DK-2100 Copenhagen<br>Denmark             | FedEx                      | 4009016505025             | 12/10/2004                  |
| <b>Finland</b>        | Finnish Institute for Verification of the Chemical Weapons Convention (VERIFIN)<br>P.O. Box 55<br>00014 University of Helsinki<br>A.I.Virtasen aukio 1,<br>00560 Helsinki<br>Finland | Air France                 | 057-4146 8000             | 21/10/2004                  |
| <b>France</b>         | Centre d'Etudes du Bouchet (CEB)<br>Section Analyses Chimiques<br>PO Box 3<br>91710 Vert-le-Petit<br>France  | FedEx                      | 40090164863               | 11/10/2004                  |

|                           |   |                              |                |             |
|---------------------------|---|------------------------------|----------------|-------------|
| <b>Hungary</b>            | Department of Toxicology<br>Institute of Health Protection<br>Hungarian Defence Forces<br>44 Róbert károly körút,<br>Budapest, H-1134<br>P.O.Box 68, Budapest, H1555<br>Hungary | FedEx                        | 40090164874    | 12/10/2004  |
| <b>India-IICT</b>         | Indian Institute of Chemical<br>Technology<br>Analytical Division<br>Tarnaka<br>Hyderabad-500 007<br>Andhra Pradesh<br>India  | Emirates<br>Airlines         | 176-7000 2936  | 19/10/2004  |
| <b>India -<br/>VERTOX</b> | Defence Research &<br>Development Establishment<br>VERTOX Laboratory<br>Jhansi Road<br>Gwalior 474002<br>India  | Cathay<br>Pacific<br>Airways | 160-3462 0622  | 19/10/2004  |
| <b>Malaysia</b>           | Department of chemistry<br>Malaysia<br>Ministry of Science,<br>Technology and Innovation<br>Jalan Sultan<br>46661 Petaling Jaya<br>Malaysia                                     | Cathay<br>Pacific<br>Airways | 160-3462 0611  | 15/10/2004  |
| <b>Morocco</b>            | Laboratoire Officiel d'Analyses<br>et de Recherches Chimiques<br>(LOARC)<br>25, Rue Nichakra Rahal<br>20 000 Casablanca<br>Morocco  | Air France                   | 057-4146 79996 | 10/10/2004* |
| <b>Netherlands</b>        | TNO Prins Mauris Laboratory<br>Lange Kleiweg 137<br>2288 GJ Rijswijk<br>The Netherlands   | FedEx                        | 40090164911    | 12/10/2004  |
| <b>Romania</b>            | NBC & Ecological Defence<br>Scientific Research Centre<br>Sos. Oltenitei 225, sector 4<br>RO-75 6872 BUCHAREST<br>Romania   | Air France                   | 057-4146 8011  | 01/11/2004  |

|                            |   |                    |                            |                          |
|----------------------------|---|--------------------|----------------------------|--------------------------|
| <b>Russia CAL</b>          | Central Chemical Weapons<br>Destruction Analytical<br>Laboratory GosNIIOKhT<br>(CAL)<br>Shosse Entusiastov, 23<br>111024 Moscow<br>Russian Federation | Air France         | 057-4146 8022              | 27/11/2004*              |
| <b>Spain</b>               | Fábrica Nacional “La<br>Marañosa”,<br>Carretera San Martin de la<br>Vega. Km. 10.5<br>San Martin de la Vega<br>Madrid 28330<br>Spain                  | FedEx              | 40090164966                | 15/10/2004               |
| <b>Sweden</b>              | Swedish Defence Research<br>Agency (FOI)<br>Division of NBC Defence<br>Cementvägen 20<br>SE-901 82 UMEÅ<br>Sweden                                     | FedEx              | 40090164970                | 14/10/2004               |
| <b>Ukraine</b>             | State Analytical Laboratory<br>61024 Kharkiv<br>Petrovskogo St. 28<br>Ukraine   | Lufthanza<br>Cargo | 020-7274 8712              | 01/11/2004*              |
| <b>UK-DSTL</b>             | Defence Science and<br>Technology Laboratory,<br>Porton Down, Salisbury<br>Wiltshire, England<br>SP4 0JQ  | FedEx              | 40090164981<br>40090164756 | 11/10/2004<br>11/10/2004 |
| <b>OPCW<br/>Laboratory</b> | OPCW Laboratory<br>Heulweg 28-30<br>2288 GN Rijswijk<br>The Netherlands   | FedEx              | 40090165003<br>40090164760 | 12/11/2004<br>12/11/2004 |

\* This date represents the date when the Proficiency Test samples arrived at the point-of-entry only. These laboratories did not respond to the preparation laboratory, so the date (if any) of receipt of the samples at their facility is unknown.



**ORGANISATION FOR THE PROHIBITION  
OF CHEMICAL WEAPONS**

**Report of the  
Sixteenth Official OPCW  
Proficiency Test**

***Part II: Purity Checks***

Laboratory code: Sample Prep Lab



|                          |
|--------------------------|
| <b>TABLE OF CONTENTS</b> |
|--------------------------|

|   | <b>Page no.</b> |
|---|-----------------|
| <b>Cover page</b> .....   | 21              |
| <b>Table of contents</b> .....  | 22              |
| Compound number A   |                 |
| Analytical technique 1: NMR $^1\text{H}$ Observation .....                  | 23              |
| Analytical technique 1: NMR $^1\text{H}$ Observation (solvent).....         | 25              |
| Analytical technique 2: NMR $^{13}\text{C}\{^1\text{H}\}$ Observation ..... | 27              |
| Compound number B   |                 |
| Analytical technique 1: NMR $^1\text{H}$ Observation .....                  | 29              |
| Analytical technique 2: NMR $^{13}\text{C}\{^1\text{H}\}$ Observation ..... | 31              |
| Compound number C   |                 |
| Analytical technique 1: NMR $^1\text{H}$ Observation .....                  | 33              |
| Analytical technique 2: NMR $^{13}\text{C}\{^1\text{H}\}$ Observation ..... | 35              |
| Compound number D   |                 |
| Analytical technique 1: NMR $^1\text{H}$ Observation .....                  | 37              |
| Analytical technique 2: NMR $^{13}\text{C}\{^1\text{H}\}$ Observation ..... | 39              |
| Compound number E & J   |                 |
| Analytical technique 1: NMR $^1\text{H}$ Observation .....                  | 41              |
| Analytical technique 2: NMR $^{13}\text{C}\{^1\text{H}\}$ Observation ..... | 43              |
| Compound number G   |                 |
| Analytical technique 1: NMR $^1\text{H}$ Observation .....                  | 45              |
| Analytical technique 2: NMR $^{13}\text{C}\{^1\text{H}\}$ Observation ..... | 47              |
| Compound number H   |                 |
| Analytical technique 1: NMR $^1\text{H}$ Observation .....                  | 49              |
| Analytical technique 2: NMR $^{13}\text{C}\{^1\text{H}\}$ Observation ..... | 51              |
| Compound number I   |                 |
| Analytical technique 1: NMR $^1\text{H}$ Observation .....                  | 53              |
| Analytical technique 2: NMR $^{13}\text{C}\{^1\text{H}\}$ Observation ..... | 55              |


## NMR METHOD AND ANALYSIS RESULTS

Laboratory code: Prep Lab    Sample code(s): Purity Check    Compound number: A

**Aliquot codes:** CW-CK-1-132-2

**Sample:** Trichloronitromethane in CD<sub>2</sub>Cl<sub>2</sub>

**Compound structure:**



**NMR Method name:** <sup>1</sup>H NMR

### METHOD DESCRIPTION

|  |  |                                  |   |
|--|--|----------------------------------|---|
| <b>Instrument</b>                        | Bruker DRX 500   |                                  |   |
| <b>Manufacturer and Type:</b>            |  |                                  |   |
| <b>Frequency:</b>                        | 500.09 MHz   | <b>Temperature control unit:</b> | <input checked="" type="checkbox"/> Yes <input type="checkbox"/> No |
| <b>Probe head:</b>                       | TBI  | <b>Temperature:</b>              | 26.5 °C   |
| <b>Sample tube diameter:</b>             | 5 mm   | <b>Solvent:</b>                  | CD <sub>2</sub> Cl <sub>2</sub>                                     |
| <b>pH of:</b>                            | <b>Sample =</b>  | <b>Blank =</b>                   | <b>Reference =</b>  |
| <b>δ reference reagent in Sample:</b>    | <input checked="" type="checkbox"/> Internal = CD <sub>2</sub> Cl <sub>2</sub> <input type="checkbox"/> External = |                                  |   |
| <b>δ reference reagent in Reference:</b> | <input type="checkbox"/> Internal = <input type="checkbox"/> External =  |                                  |   |

### ANALYSIS

|  |   |
|--|---|
| <input type="checkbox"/> Compared to reference chemical:     | Source : <input type="checkbox"/> Own Synthesis <input type="checkbox"/> Commercial   |
| <input type="checkbox"/> Compared to library spectrum:       | Source : <input type="checkbox"/> OCAD (Code      ) <input type="checkbox"/> Other:   |
| <input type="checkbox"/> Standard addition:                  | Source :  |
| <input checked="" type="checkbox"/> Spectral interpretation: | E.Pretsch, P. Buhlmann, C. Affolter, <u>Structure Determination of Organic Compounds: Tables of Spectral Data</u> , (Berlin: Springer-Verlag) 2000. |

|                     |   |
|---------------------|---|
| <b>Experiments:</b> | <input checked="" type="checkbox"/> <sup>1</sup> H <input type="checkbox"/> <sup>13</sup> C{ <sup>1</sup> H} <input type="checkbox"/> <sup>19</sup> F <input type="checkbox"/> <sup>31</sup> P{ <sup>1</sup> H} <input type="checkbox"/> <sup>31</sup> P <input type="checkbox"/> COSY<br><input type="checkbox"/> Other: |
|---------------------|---|

| Assignment(s):        | Chemical shift(s)   | Coupling constant(s)  |                          |  |  |  |  |  |  |  |  |   |                      |                         |  |  |  |  |  |  |  |  |
|-----------------------|---|-----------------------|--------------------------|--|--|--|--|--|--|--|--|---|----------------------|-------------------------|--|--|--|--|--|--|--|--|
|                       | <table style="width: 100%;"> <tr> <th style="width: 50%;">Sample spectrum [ppm]</th> <th style="width: 50%;">Reference spectrum [ppm]</th> </tr> <tr><td> </td><td> </td></tr> <tr><td> </td><td> </td></tr> <tr><td> </td><td> </td></tr> <tr><td> </td><td> </td></tr> </table> | Sample spectrum [ppm] | Reference spectrum [ppm] |  |  |  |  |  |  |  |  | <table style="width: 100%;"> <tr> <th style="width: 50%;">Sample spectrum [Hz]</th> <th style="width: 50%;">Reference spectrum [Hz]</th> </tr> <tr><td> </td><td> </td></tr> <tr><td> </td><td> </td></tr> <tr><td> </td><td> </td></tr> <tr><td> </td><td> </td></tr> </table> | Sample spectrum [Hz] | Reference spectrum [Hz] |  |  |  |  |  |  |  |  |
| Sample spectrum [ppm] | Reference spectrum [ppm]  |                       |                          |  |  |  |  |  |  |  |  |   |                      |                         |  |  |  |  |  |  |  |  |
|                       |   |                       |                          |  |  |  |  |  |  |  |  |   |                      |                         |  |  |  |  |  |  |  |  |
|                       |   |                       |                          |  |  |  |  |  |  |  |  |   |                      |                         |  |  |  |  |  |  |  |  |
|                       |   |                       |                          |  |  |  |  |  |  |  |  |   |                      |                         |  |  |  |  |  |  |  |  |
|                       |   |                       |                          |  |  |  |  |  |  |  |  |   |                      |                         |  |  |  |  |  |  |  |  |
| Sample spectrum [Hz]  | Reference spectrum [Hz]   |                       |                          |  |  |  |  |  |  |  |  |   |                      |                         |  |  |  |  |  |  |  |  |
|                       |   |                       |                          |  |  |  |  |  |  |  |  |   |                      |                         |  |  |  |  |  |  |  |  |
|                       |   |                       |                          |  |  |  |  |  |  |  |  |   |                      |                         |  |  |  |  |  |  |  |  |
|                       |   |                       |                          |  |  |  |  |  |  |  |  |   |                      |                         |  |  |  |  |  |  |  |  |
|                       |   |                       |                          |  |  |  |  |  |  |  |  |   |                      |                         |  |  |  |  |  |  |  |  |

|                                |   |
|--------------------------------|---|
| <b>Interpretative comments</b> | There are no protons in this sample. All peaks in <sup>1</sup> H NMR spectrum are due to the solvent and its impurities (confirmed by background NMR analysis of the solvent). Purity checks were performed with <sup>13</sup> C{ <sup>1</sup> H} analysis. |
|--------------------------------|---|

Current Data Parameters  
NAME CWCK11322  
EXPNO 1  
PROCNO 1



F2 - Acquisition Parameters  
Date\_ 20041123  
Time 15:37  
INSTRUM spect  
PROBHD 5 mm TBI 1H/13  
PULPROG zg  
TD 131072  
SOLVENT CD2CL2  
NS 4  
DS 0  
SWH 7507.507 Hz  
FIDRES 0.057278 Hz  
AQ 8.7295122 sec  
RG 100  
DW 66.600 usec  
DE 6.00 usec  
TE 297.5 K  
D1 1.00000000 sec  
MCREST 0.00000000 sec  
MCWRK 0.01500000 sec

===== CHANNEL f1 =====  
NUC1 1H  
P1 4.00 usec  
PL1 0.00 dB  
SFO1 500.0894638 MHz

F2 - Processing parameters  
SI 65536  
SF 500.0870229 MHz  
WDW no  
SSB 0  
LB 0.00 Hz  
GB 0  
PC 1.00

1D NMR plot parameters  
CX 20.00 cm  
CY 12.50 cm  
FIP 6.934 ppm  
FI 3467.71 Hz  
F2P 0.609 ppm  
F2 304.50 Hz  
PPMCM 0.31627 ppm/cm  
HZCM 158.16031 Hz/cm

Compound A



<sup>1</sup>H spectrum of Trichloronitromethane  
(Solvent: CD<sub>2</sub>Cl<sub>2</sub>; Internal reference: CH<sub>2</sub>Cl<sub>2</sub> at 5.32 ppm)


## NMR METHOD AND ANALYSIS RESULTS

Laboratory code: Prep Lab    Sample code(s): Purity Check    Compound number: (solvent)

**Aliquot codes:** CW-CK-1-132-1

**Sample:** Solvent background CD<sub>2</sub>Cl<sub>2</sub>

**Compound structure:**



**NMR Method name:** <sup>1</sup>H NMR

### METHOD DESCRIPTION


|  |  |                                  |   |
|--|--|----------------------------------|---|
| <b>Instrument</b>                        | Bruker DRX 500   |                                  |   |
| <b>Manufacturer and Type:</b>            |  |                                  |   |
| <b>Frequency:</b>                        | 500.09 MHz   | <b>Temperature control unit:</b> | <input checked="" type="checkbox"/> Yes <input type="checkbox"/> No |
| <b>Probe head:</b>                       | TBI  | <b>Temperature:</b>              | 26.5 °C   |
| <b>Sample tube diameter:</b>             | 5 mm   | <b>Solvent:</b>                  | CD <sub>2</sub> Cl <sub>2</sub>                                     |
| <b>pH of:</b>                            | <b>Sample =</b>  | <b>Blank =</b>                   | <b>Reference =</b>  |
| <b>δ reference reagent in Sample:</b>    | <input checked="" type="checkbox"/> Internal = CD <sub>2</sub> Cl <sub>2</sub> <input type="checkbox"/> External = |                                  |   |
| <b>δ reference reagent in Reference:</b> | <input type="checkbox"/> Internal = <input type="checkbox"/> External =  |                                  |   |

### ANALYSIS

|  |   |                                 |                                |
|--|---|---------------------------------|--------------------------------|
| <input type="checkbox"/> Compared to reference chemical:     | Source : <input type="checkbox"/> Own Synthesis <input type="checkbox"/> Commercial   |                                 |                                |
| <input type="checkbox"/> Compared to library spectrum:       | Source : <input type="checkbox"/> OCAD (Code      ) <input type="checkbox"/> Other:   |                                 |                                |
| <input type="checkbox"/> Standard addition:                  | Source :  |                                 |                                |
| <input checked="" type="checkbox"/> Spectral interpretation: | E.Pretsch, P. Buhlmann, C. Affolter, <u>Structure Determination of Organic Compounds: Tables of Spectral Data</u> , (Berlin: Springer-Verlag) 2000.   |                                 |                                |
| <b>Experiments:</b>  | <input checked="" type="checkbox"/> <sup>1</sup> H <input type="checkbox"/> <sup>13</sup> C{ <sup>1</sup> H} <input type="checkbox"/> <sup>19</sup> F <input type="checkbox"/> <sup>31</sup> P{ <sup>1</sup> H} <input type="checkbox"/> <sup>31</sup> P <input type="checkbox"/> COSY<br><input type="checkbox"/> Other: |                                 |                                |
| <b>Assignment(s):</b>  | <b>Chemical shift(s)</b>  |                                 | <b>Coupling constant(s)</b>    |
|  | <b>Sample spectrum [ppm]</b>  | <b>Reference spectrum [ppm]</b> | <b>Sample spectrum [Hz]</b>    |
|  |   |                                 | <b>Reference spectrum [Hz]</b> |
|  |   |                                 |                                |
|  |   |                                 |                                |
|  |   |                                 |                                |
| <b>Interpretative comments</b>                               | Background <sup>1</sup> H NMR of solvent confirming solvent impurities shown in <sup>1</sup> H NMR of sample A.   |                                 |                                |

Current Data Parameters  
 NAME CWCK11321  
 EXPNO 1  
 PROCNO 1


Methylene chloride



F2 - Acquisition Parameters  
 Date\_ 20041123  
 Time 14.28  
 INSTRUM spect  
 PROBHD 5 mm TBI 1H/13  
 PULPROG zg  
 TD 131072  
 SOLVENT CD2Cl2  
 NS 16  
 DS 0  
 SWH 7507.507 Hz  
 FIDRES 0.057278 Hz  
 AQ 8.7295122 sec  
 RG 100  
 DW 66.600 usec  
 DE 6.00 usec  
 TE 297.6 K  
 D1 1.00000000 sec  
 MCREST 0.00000000 sec  
 MCWRK 0.01500000 sec

===== CHANNEL f1 =====  
 NUC1 1H  
 P1 4.00 usec  
 PL1 0.00 dB  
 SFO1 500.0894638 MHz

F2 - Processing parameters  
 SI 65536  
 SF 500.0870229 MHz  
 WDW no  
 SSB 0  
 LB 0.00 Hz  
 GB 0  
 PC 1.00



<sup>1</sup>H spectrum of methylene chloride - d<sub>2</sub>  
 (Solvent: CD<sub>2</sub>Cl<sub>2</sub>; Internal reference: CH<sub>2</sub>Cl<sub>2</sub> at 5.32 ppm)


# NMR METHOD AND ANALYSIS RESULTS

Laboratory code: Prep Lab    Sample code(s): Purity Check    Compound number: A

**Aliquot codes:** CW-CK-1-131-1

**Sample:** Trichloronitromethane in CD<sub>2</sub>Cl<sub>2</sub>

**Compound structure:**



**NMR Method name:** <sup>13</sup>C{<sup>1</sup>H} NMR

## METHOD DESCRIPTION

|  |  |                                  |   |
|--|--|----------------------------------|---|
| <b>Instrument</b>                        | Bruker DRX 500   |                                  |   |
| <b>Manufacturer and Type:</b>            |  |                                  |   |
| <b>Frequency:</b>                        | 125.75 MHz   | <b>Temperature control unit:</b> | <input checked="" type="checkbox"/> Yes <input type="checkbox"/> No |
| <b>Probe head:</b>                       | TBI  | <b>Temperature:</b>              | 26.5 °C   |
| <b>Sample tube diameter:</b>             | 5 mm   | <b>Solvent:</b>                  | CH <sub>2</sub> Cl <sub>2</sub> /CD <sub>2</sub> Cl <sub>2</sub>    |
| <b>pH of:</b>                            | <b>Sample =</b>  | <b>Blank =</b>                   | <b>Reference =</b>  |
| <b>δ reference reagent in Sample:</b>    | <input checked="" type="checkbox"/> Internal = CD <sub>2</sub> Cl <sub>2</sub> <input type="checkbox"/> External = |                                  |   |
| <b>δ reference reagent in Reference:</b> | <input type="checkbox"/> Internal = <input type="checkbox"/> External =  |                                  |   |

## ANALYSIS

|  |   |
|--|---|
| <input type="checkbox"/> Compared to reference chemical:     | Source : <input type="checkbox"/> Own Synthesis <input type="checkbox"/> Commercial   |
| <input type="checkbox"/> Compared to library spectrum:       | Source : <input type="checkbox"/> OCAD (Code      ) <input type="checkbox"/> Other:   |
| <input type="checkbox"/> Standard addition:                  | Source :  |
| <input checked="" type="checkbox"/> Spectral interpretation: | E.Pretsch, P. Buhlmann, C. Affolter, <u>Structure Determination of Organic Compounds: Tables of Spectral Data</u> , (Berlin: Springer-Verlag) 2000. |


|                     |   |
|---------------------|---|
| <b>Experiments:</b> | <input type="checkbox"/> <sup>1</sup> H <input checked="" type="checkbox"/> <sup>13</sup> C{ <sup>1</sup> H} <input type="checkbox"/> <sup>19</sup> F <input type="checkbox"/> <sup>31</sup> P{ <sup>1</sup> H} <input type="checkbox"/> <sup>31</sup> P <input type="checkbox"/> COSY<br><input type="checkbox"/> Other: |
|---------------------|---|

| Assignment(s):        | Chemical shift(s)   | Coupling constant(s)  |                          |   |            |  |  |  |  |  |  |   |                      |                         |  |  |  |  |  |  |  |  |
|-----------------------|---|-----------------------|--------------------------|---|------------|--|--|--|--|--|--|---|----------------------|-------------------------|--|--|--|--|--|--|--|--|
|                       | <table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <th style="width: 50%;">Sample spectrum [ppm]</th> <th style="width: 50%;">Reference spectrum [ppm]</th> </tr> <tr> <td>a</td> <td>114.92 ppm</td> </tr> <tr><td> </td><td> </td></tr> <tr><td> </td><td> </td></tr> <tr><td> </td><td> </td></tr> </table> | Sample spectrum [ppm] | Reference spectrum [ppm] | a | 114.92 ppm |  |  |  |  |  |  | <table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <th style="width: 50%;">Sample spectrum [Hz]</th> <th style="width: 50%;">Reference spectrum [Hz]</th> </tr> <tr><td> </td><td> </td></tr> <tr><td> </td><td> </td></tr> <tr><td> </td><td> </td></tr> <tr><td> </td><td> </td></tr> </table> | Sample spectrum [Hz] | Reference spectrum [Hz] |  |  |  |  |  |  |  |  |
| Sample spectrum [ppm] | Reference spectrum [ppm]  |                       |                          |   |            |  |  |  |  |  |  |   |                      |                         |  |  |  |  |  |  |  |  |
| a                     | 114.92 ppm  |                       |                          |   |            |  |  |  |  |  |  |   |                      |                         |  |  |  |  |  |  |  |  |
|                       |   |                       |                          |   |            |  |  |  |  |  |  |   |                      |                         |  |  |  |  |  |  |  |  |
|                       |   |                       |                          |   |            |  |  |  |  |  |  |   |                      |                         |  |  |  |  |  |  |  |  |
|                       |   |                       |                          |   |            |  |  |  |  |  |  |   |                      |                         |  |  |  |  |  |  |  |  |
| Sample spectrum [Hz]  | Reference spectrum [Hz]   |                       |                          |   |            |  |  |  |  |  |  |   |                      |                         |  |  |  |  |  |  |  |  |
|                       |   |                       |                          |   |            |  |  |  |  |  |  |   |                      |                         |  |  |  |  |  |  |  |  |
|                       |   |                       |                          |   |            |  |  |  |  |  |  |   |                      |                         |  |  |  |  |  |  |  |  |
|                       |   |                       |                          |   |            |  |  |  |  |  |  |   |                      |                         |  |  |  |  |  |  |  |  |
|                       |   |                       |                          |   |            |  |  |  |  |  |  |   |                      |                         |  |  |  |  |  |  |  |  |

|                                |               |
|--------------------------------|---------------|
| <b>Interpretative comments</b> | Purity = 100% |
|--------------------------------|---------------|




# NMR METHOD AND ANALYSIS RESULTS

Laboratory code: Prep Lab      Sample code(s): Purity Check      Compound number: B

**Aliquot codes:** CW-CK-1-89-2

**Sample:** 2-(N-Ethyl-N-propylamino)ethylchloride in CH<sub>2</sub>Cl<sub>2</sub>/CD<sub>2</sub>Cl<sub>2</sub>

**Compound structure:**



**NMR Method name:** <sup>1</sup>H NMR


## METHOD DESCRIPTION

|  |  |                                  |   |
|--|--|----------------------------------|---|
| <b>Instrument</b>                        | Bruker DRX 500   |                                  |   |
| <b>Manufacturer and Type:</b>            |  |                                  |   |
| <b>Frequency:</b>                        | 500.09 MHz   | <b>Temperature control unit:</b> | <input checked="" type="checkbox"/> Yes <input type="checkbox"/> No |
| <b>Probe head:</b>                       | TBI  | <b>Temperature:</b>              | 26.5 °C   |
| <b>Sample tube diameter:</b>             | 5 mm   | <b>Solvent:</b>                  | CH <sub>2</sub> Cl <sub>2</sub> /CD <sub>2</sub> Cl <sub>2</sub>    |
| <b>pH of:</b>                            | <b>Sample =</b>  | <b>Blank =</b>                   | <b>Reference =</b>  |
| <b>δ reference reagent in Sample:</b>    | <input checked="" type="checkbox"/> Internal = CD <sub>2</sub> Cl <sub>2</sub> <input type="checkbox"/> External = |                                  |   |
| <b>δ reference reagent in Reference:</b> | <input type="checkbox"/> Internal = <input type="checkbox"/> External =  |                                  |   |

## ANALYSIS

|  |   |   |  |
|--|---|---|--|
| <input type="checkbox"/> Compared to reference chemical:     |   | Source : <input type="checkbox"/> Own Synthesis <input type="checkbox"/> Commercial   |  |
| <input type="checkbox"/> Compared to library spectrum:       |   | Source : <input type="checkbox"/> OCAD (Code ) <input type="checkbox"/> Other:  |  |
| <input type="checkbox"/> Standard addition:                  |   | Source :  |  |
| <input checked="" type="checkbox"/> Spectral interpretation: |   | E.Pretsch, P. Buhlmann, C. Affolter, <u>Structure Determination of Organic Compounds: Tables of Spectral Data</u> , (Berlin: Springer-Verlag) 2000. |  |
| Experiments:   | <input checked="" type="checkbox"/> <sup>1</sup> H <input type="checkbox"/> <sup>13</sup> C{ <sup>1</sup> H} <input type="checkbox"/> <sup>19</sup> F <input type="checkbox"/> <sup>31</sup> P{ <sup>1</sup> H} <input type="checkbox"/> <sup>31</sup> P <input type="checkbox"/> COSY<br><input type="checkbox"/> Other: |   |  |
| Assignment(s):   | Chemical shift(s)   |   | Coupling constant(s)                         |
|  | Sample spectrum [ppm]   | Reference spectrum [ppm]  | Sample spectrum [Hz] Reference spectrum [Hz] |
| a  | 0.86  |   | 7.5 (b,a)                                    |
| e  | 0.99  |   | 7.1 (d,e)                                    |
| b  | 1.42  |   | 7.5 (a,b; c,b)                               |
| c  | 2.41  |   | 7.5 (b,c)                                    |
| d  | 2.53  |   | 7.1 (e,d)                                    |
| f  | 2.74  |   | 7.5 (g,f)                                    |
| g  | 3.49  |   | 7.5 (f,g)                                    |
| Interpretative comments                                      | Purity = 99.96%   |   |  |






# NMR METHOD AND ANALYSIS RESULTS

Laboratory code: Prep Lab      Sample code(s): Purity Check      Compound number: B

**Aliquot codes:** CW-CK-1-89-2

**Sample:** 2-(N-Ethyl-N-propylamino)ethylchloride in CH<sub>2</sub>Cl<sub>2</sub>/CD<sub>2</sub>Cl<sub>2</sub>

**Compound structure:**



**NMR Method name:** <sup>13</sup>C{<sup>1</sup>H} NMR

## METHOD DESCRIPTION


|  |  |                                  |   |
|--|--|----------------------------------|---|
| <b>Instrument</b>                        | Bruker DRX 500   |                                  |   |
| <b>Manufacturer and Type:</b>            |  |                                  |   |
| <b>Frequency:</b>                        | 125.75 MHz   | <b>Temperature control unit:</b> | <input checked="" type="checkbox"/> Yes <input type="checkbox"/> No |
| <b>Probe head:</b>                       | TBI  | <b>Temperature:</b>              | 26.5 °C   |
| <b>Sample tube diameter:</b>             | 5 mm   | <b>Solvent:</b>                  | CH <sub>2</sub> Cl <sub>2</sub> /CD <sub>2</sub> Cl <sub>2</sub>    |
| <b>pH of:</b>                            | <b>Sample =</b>  | <b>Blank =</b>                   | <b>Reference =</b>  |
| <b>δ reference reagent in Sample:</b>    | <input checked="" type="checkbox"/> Internal = CD <sub>2</sub> Cl <sub>2</sub> <input type="checkbox"/> External = |                                  |   |
| <b>δ reference reagent in Reference:</b> | <input type="checkbox"/> Internal = <input type="checkbox"/> External =  |                                  |   |

## ANALYSIS

|  |   |
|--|---|
| <input type="checkbox"/> Compared to reference chemical:     | Source : <input type="checkbox"/> Own Synthesis <input type="checkbox"/> Commercial   |
| <input type="checkbox"/> Compared to library spectrum:       | Source : <input type="checkbox"/> OCAD (Code      ) <input type="checkbox"/> Other:   |
| <input type="checkbox"/> Standard addition:                  | Source :  |
| <input checked="" type="checkbox"/> Spectral interpretation: | E.Pretsch, P. Buhlmann, C. Affolter, <u>Structure Determination of Organic Compounds: Tables of Spectral Data</u> , (Berlin: Springer-Verlag) 2000. |

|                     |   |
|---------------------|---|
| <b>Experiments:</b> | <input type="checkbox"/> <sup>1</sup> H <input checked="" type="checkbox"/> <sup>13</sup> C{ <sup>1</sup> H} <input type="checkbox"/> <sup>19</sup> F <input type="checkbox"/> <sup>31</sup> P{ <sup>1</sup> H} <input type="checkbox"/> <sup>31</sup> P <input type="checkbox"/> COSY<br><input type="checkbox"/> Other: |
|---------------------|---|

| Assignment(s):                 | Chemical shift(s)     | Coupling constant(s)     |
|--------------------------------|-----------------------|--------------------------|
|                                | Sample spectrum [ppm] | Reference spectrum [ppm] |
|                                | Sample spectrum [Hz]  | Reference spectrum [Hz]  |
| a                              | 11.51                 |                          |
| e                              | 11.94                 |                          |
| b                              | 20.63                 |                          |
| g                              | 42.39                 |                          |
| d                              | 48.08                 |                          |
| c                              | 55.63                 |                          |
| f                              | 55.99                 |                          |
| <b>Interpretative comments</b> |                       |                          |




# NMR METHOD AND ANALYSIS RESULTS

Laboratory code: Prep Lab      Sample code(s): Purity Check      Compound number: C

Aliquot codes: CW-CK-1-131-2

Sample: 2-(N,N-diisopropylamino)ethylchloride in CD<sub>2</sub>Cl<sub>2</sub>

Compound structure:




NMR Method name: <sup>1</sup>H NMR

## METHOD DESCRIPTION

|                                   |  |                           |   |
|-----------------------------------|--|---------------------------|---|
| Instrument                        | Bruker DRX 500   |                           |   |
| Manufacturer and Type:            |  |                           |   |
| Frequency:                        | 500.08 MHz   | Temperature control unit: | <input checked="" type="checkbox"/> Yes <input type="checkbox"/> No |
| Probe head:                       | TBI  | Temperature:              | 26.5 °C   |
| Sample tube diameter:             | 5 mm   | Solvent:                  | CD <sub>2</sub> Cl <sub>2</sub>                                     |
| pH of:                            | Sample =   | Blank =                   | Reference =   |
| δ reference reagent in Sample:    | <input checked="" type="checkbox"/> Internal = CD <sub>2</sub> Cl <sub>2</sub> <input type="checkbox"/> External = |                           |   |
| δ reference reagent in Reference: | <input type="checkbox"/> Internal = <input type="checkbox"/> External =  |                           |   |

## ANALYSIS

|  |   |   |   |
|--|---|---|---|
| <input type="checkbox"/> Compared to reference chemical:     |   | Source : <input type="checkbox"/> Own Synthesis <input type="checkbox"/> Commercial   |   |
| <input type="checkbox"/> Compared to library spectrum:       |   | Source : <input type="checkbox"/> OCAD (Code ) <input type="checkbox"/> Other:  |   |
| <input type="checkbox"/> Standard addition:                  |   | Source :  |   |
| <input checked="" type="checkbox"/> Spectral interpretation: |   | E.Pretsch, P. Buhlmann, C. Affolter, <u>Structure Determination of Organic Compounds: Tables of Spectral Data</u> , (Berlin: Springer-Verlag) 2000. |   |
| Experiments:   | <input checked="" type="checkbox"/> <sup>1</sup> H <input type="checkbox"/> <sup>13</sup> C{ <sup>1</sup> H} <input type="checkbox"/> <sup>19</sup> F <input type="checkbox"/> <sup>31</sup> P{ <sup>1</sup> H} <input type="checkbox"/> <sup>31</sup> P <input type="checkbox"/> COSY<br><input type="checkbox"/> Other: |   |   |
| Assignment(s):   | Chemical shift(s)   |   | Coupling constant(s)                                    |
|  | Sample spectrum<br>[ppm]  | Reference spectrum<br>[ppm]   | Sample spectrum<br>[Hz]      Reference spectrum<br>[Hz] |
| a,b,e,f  | 1.00, 0.99  |   | 6.5 (a,c; b,c; e,d; f;d)                                |
| g  | 2.72  |   | 8.1 (g,h)   |
| c,d  | 2.99  |   | 6.5 (c,a; c,b; d,e; d,f)                                |
| h  | 3.36  |   | 8.1 (h,g)   |
| Interpretative<br>comments                                   | Purity > 99.9%  |   |   |



<sup>1</sup>H spectrum of 2-(N,N-Diisopropylamino)ethylchloride  
(Solvent: CD<sub>2</sub>Cl<sub>2</sub>; Internal reference: CH<sub>2</sub>Cl<sub>2</sub> at 5.32 ppm)



Current Data Parameters  
NAME CWCK11312  
EXPNO 2  
PROCNO 1



F2 - Acquisition Parameters  
Date\_ 20041123  
Time 16.25  
INSTRUM spect  
PROBHD 5 mm TBI 1H/13  
PULPROG zgpg  
TD 131072  
SOLVENT CD2Cl2  
NS 9100  
DS 0  
SWH 50505.051 Hz  
FIDRES 0.385323 Hz  
AQ 1.2976727 sec  
RG 4096  
DE 9.900 usec  
TE 297.4 K  
D1 5.00000000 sec  
d11 0.03000000 sec  
MCREST 0.00000000 sec  
MCWRK 0.01500000 sec

===== CHANNEL f1 =====  
NUC1 13C  
P1 9.00 usec  
PL1 -2.00 dB  
SFO1 125.7533079 MHz

===== CHANNEL f2 =====  
CPDPRG2 waltz16  
NUC2 1H  
PCPD2 85.00 usec  
PL2 0.00 dB  
PL12 23.60 dB  
SFO2 500.0895004 MHz

F2 - Processing parameters  
SI 131072  
SF 125.7469265 MHz  
WDW EM  
SSB 0  
LB 0.30 Hz  
GB 0  
PC 1.00

Compound C



$^{13}\text{C}\{^1\text{H}\}$  spectrum of 2-(N,N-diisopropylamino)ethylchloride  
(Solvent:  $\text{CD}_2\text{Cl}_2$ ; Internal reference:  $\text{CH}_2\text{Cl}_2$  at 53.8 ppm)


# NMR METHOD AND ANALYSIS RESULTS

Laboratory code: Prep Lab    Sample code(s): Purity Check    Compound number: D

**Aliquot codes:** CW-CK-1-88-3

**Sample:** 2-(N,N-Diethylamino)ethylchloride in CH<sub>2</sub>Cl<sub>2</sub>/CD<sub>2</sub>Cl<sub>2</sub>

**Compound structure:**



**NMR Method name:** <sup>1</sup>H NMR

## METHOD DESCRIPTION

|  |  |                                  |   |
|--|--|----------------------------------|---|
| <b>Instrument</b>                        | Bruker DRX 500   |                                  |   |
| <b>Manufacturer and Type:</b>            |  |                                  |   |
| <b>Frequency:</b>                        | 500.09 MHz   | <b>Temperature control unit:</b> | <input checked="" type="checkbox"/> Yes <input type="checkbox"/> No |
| <b>Probe head:</b>                       | TBI  | <b>Temperature:</b>              | 26.5 °C   |
| <b>Sample tube diameter:</b>             | 5 mm   | <b>Solvent:</b>                  | CH <sub>2</sub> Cl <sub>2</sub> /CD <sub>2</sub> Cl <sub>2</sub>    |
| <b>pH of:</b>                            | <b>Sample =</b>  | <b>Blank =</b>                   | <b>Reference =</b>  |
| <b>δ reference reagent in Sample:</b>    | <input checked="" type="checkbox"/> Internal = CD <sub>2</sub> Cl <sub>2</sub> <input type="checkbox"/> External = |                                  |   |
| <b>δ reference reagent in Reference:</b> | <input type="checkbox"/> Internal = <input type="checkbox"/> External =  |                                  |   |


## ANALYSIS

| <input type="checkbox"/> Compared to reference chemical:     | Source : <input type="checkbox"/> Own Synthesis <input type="checkbox"/> Commercial   |                          |                      |                         |                      |  |                       |                          |                      |                         |     |      |  |                |  |     |      |  |                |  |   |      |  |           |  |   |      |  |           |  |
|--|---|--------------------------|----------------------|-------------------------|----------------------|--|-----------------------|--------------------------|----------------------|-------------------------|-----|------|--|----------------|--|-----|------|--|----------------|--|---|------|--|-----------|--|---|------|--|-----------|--|
| <input type="checkbox"/> Compared to library spectrum:       | Source : <input type="checkbox"/> OCAD (Code      ) <input type="checkbox"/> Other:   |                          |                      |                         |                      |  |                       |                          |                      |                         |     |      |  |                |  |     |      |  |                |  |   |      |  |           |  |   |      |  |           |  |
| <input type="checkbox"/> Standard addition:                  | Source :  |                          |                      |                         |                      |  |                       |                          |                      |                         |     |      |  |                |  |     |      |  |                |  |   |      |  |           |  |   |      |  |           |  |
| <input checked="" type="checkbox"/> Spectral interpretation: | E.Pretsch, P. Buhlmann, C. Affolter, <u>Structure Determination of Organic Compounds: Tables of Spectral Data</u> , (Berlin: Springer-Verlag) 2000.   |                          |                      |                         |                      |  |                       |                          |                      |                         |     |      |  |                |  |     |      |  |                |  |   |      |  |           |  |   |      |  |           |  |
| <b>Experiments:</b>  | <input checked="" type="checkbox"/> <sup>1</sup> H <input type="checkbox"/> <sup>13</sup> C{ <sup>1</sup> H} <input type="checkbox"/> <sup>19</sup> F <input type="checkbox"/> <sup>31</sup> P{ <sup>1</sup> H} <input type="checkbox"/> <sup>31</sup> P <input type="checkbox"/> COSY<br><input type="checkbox"/> Other:   |                          |                      |                         |                      |  |                       |                          |                      |                         |     |      |  |                |  |     |      |  |                |  |   |      |  |           |  |   |      |  |           |  |
| <b>Assignment(s):</b>  | <table border="1" style="width: 100%;"> <thead> <tr> <th rowspan="2"></th> <th colspan="2">Chemical shift(s)</th> <th colspan="2">Coupling constant(s)</th> </tr> <tr> <th>Sample spectrum [ppm]</th> <th>Reference spectrum [ppm]</th> <th>Sample spectrum [Hz]</th> <th>Reference spectrum [Hz]</th> </tr> </thead> <tbody> <tr> <td>a,d</td> <td>1.00</td> <td></td> <td>7.0 (b,a; c,d)</td> <td></td> </tr> <tr> <td>b,c</td> <td>2.54</td> <td></td> <td>7.0 (a,b; d,c)</td> <td></td> </tr> <tr> <td>e</td> <td>2.74</td> <td></td> <td>7.0 (f,e)</td> <td></td> </tr> <tr> <td>f</td> <td>3.50</td> <td></td> <td>7.0 (e,f)</td> <td></td> </tr> </tbody> </table> |                          | Chemical shift(s)    |                         | Coupling constant(s) |  | Sample spectrum [ppm] | Reference spectrum [ppm] | Sample spectrum [Hz] | Reference spectrum [Hz] | a,d | 1.00 |  | 7.0 (b,a; c,d) |  | b,c | 2.54 |  | 7.0 (a,b; d,c) |  | e | 2.74 |  | 7.0 (f,e) |  | f | 3.50 |  | 7.0 (e,f) |  |
|  | Chemical shift(s)   |                          | Coupling constant(s) |                         |                      |  |                       |                          |                      |                         |     |      |  |                |  |     |      |  |                |  |   |      |  |           |  |   |      |  |           |  |
|  | Sample spectrum [ppm]   | Reference spectrum [ppm] | Sample spectrum [Hz] | Reference spectrum [Hz] |                      |  |                       |                          |                      |                         |     |      |  |                |  |     |      |  |                |  |   |      |  |           |  |   |      |  |           |  |
| a,d  | 1.00  |                          | 7.0 (b,a; c,d)       |                         |                      |  |                       |                          |                      |                         |     |      |  |                |  |     |      |  |                |  |   |      |  |           |  |   |      |  |           |  |
| b,c  | 2.54  |                          | 7.0 (a,b; d,c)       |                         |                      |  |                       |                          |                      |                         |     |      |  |                |  |     |      |  |                |  |   |      |  |           |  |   |      |  |           |  |
| e  | 2.74  |                          | 7.0 (f,e)            |                         |                      |  |                       |                          |                      |                         |     |      |  |                |  |     |      |  |                |  |   |      |  |           |  |   |      |  |           |  |
| f  | 3.50  |                          | 7.0 (e,f)            |                         |                      |  |                       |                          |                      |                         |     |      |  |                |  |     |      |  |                |  |   |      |  |           |  |   |      |  |           |  |
| <b>Interpretative comments</b>                               | Purity = 99.9%  |                          |                      |                         |                      |  |                       |                          |                      |                         |     |      |  |                |  |     |      |  |                |  |   |      |  |           |  |   |      |  |           |  |



Compound D

Current Data Parameters  
 NAME CWCK1883  
 EXPNO 1  
 PROCNO 1  
 F2 - Acquisition Parameters  
 Date\_ 20041001  
 Time 19.13  
 INSTRUM spect  
 PROBHD 5 mm TBI 1H/13  
 PULPROG zg  
 TD 131072  
 SOLVENT CD2Cl2  
 NS 64  
 DS 0  
 SWH 7507.507 Hz  
 FIDRES 0.057278 Hz  
 AQ 8.7295122 sec  
 RG 3  
 DW 66.600 usec  
 DE 6.00 usec  
 TE 671.7 K  
 D1 1.00000000 sec  
 MCREST 0.00000000 sec  
 MCWRK 0.01500000 sec  
 ===== CHANNEL f1 =====  
 NUC1 1H  
 P1 5.40 usec  
 PL1 0.00 dB  
 SFO1 500.0894638 MHz  
 F2 - Processing parameters  
 SI 65536  
 SF 500.0870229 MHz  
 WDW no  
 SSB 0  
 LB 0.00 Hz  
 GB 0  
 PC 1.00



<sup>1</sup>H spectrum of 2-(N,N-diethylamino)ethylchloride  
 (Solvent: CH<sub>2</sub>Cl<sub>2</sub> + CD<sub>2</sub>Cl<sub>2</sub>; External reference: CD<sub>2</sub>Cl<sub>2</sub> at 5.32 ppm)


## NMR METHOD AND ANALYSIS RESULTS

Laboratory code: Prep Lab    Sample code(s): Purity Check    Compound number: D

**Aliquot codes:** CW-CK-1-88-3

**Sample:** 2-(N,N-Diethylamino)ethylchloride in CH<sub>2</sub>Cl<sub>2</sub>/CD<sub>2</sub>Cl<sub>2</sub>

**Compound structure:**



**NMR Method name:** <sup>13</sup>C{<sup>1</sup>H} NMR

### METHOD DESCRIPTION

|  |  |                                  |   |
|--|--|----------------------------------|---|
| <b>Instrument</b>                        | Bruker DRX 500   |                                  |   |
| <b>Manufacturer and Type:</b>            |  |                                  |   |
| <b>Frequency:</b>                        | 125.75 MHz   | <b>Temperature control unit:</b> | <input checked="" type="checkbox"/> Yes <input type="checkbox"/> No |
| <b>Probe head:</b>                       | TBI  | <b>Temperature:</b>              | 26.5 °C   |
| <b>Sample tube diameter:</b>             | 5 mm   | <b>Solvent:</b>                  | CH <sub>2</sub> Cl <sub>2</sub> /CD <sub>2</sub> Cl <sub>2</sub>    |
| <b>pH of:</b>                            | <b>Sample =</b>  | <b>Blank =</b>                   | <b>Reference =</b>  |
| <b>δ reference reagent in Sample:</b>    | <input checked="" type="checkbox"/> Internal = CD <sub>2</sub> Cl <sub>2</sub> <input type="checkbox"/> External = |                                  |   |
| <b>δ reference reagent in Reference:</b> | <input type="checkbox"/> Internal = <input type="checkbox"/> External =  |                                  |   |

### ANALYSIS

|  |   |  |  |
|--|---|--|--|
| <input type="checkbox"/> Compared to reference chemical:     | Source : <input type="checkbox"/> Own Synthesis <input type="checkbox"/> Commercial   |  |  |
| <input type="checkbox"/> Compared to library spectrum:       | Source : <input type="checkbox"/> OCAD (Code      ) <input type="checkbox"/> Other:   |  |  |
| <input type="checkbox"/> Standard addition:                  | Source :  |  |  |
| <input checked="" type="checkbox"/> Spectral interpretation: | E.Pretsch, P. Buhlmann, C. Affolter, <u>Structure Determination of Organic Compounds: Tables of Spectral Data</u> , (Berlin: Springer-Verlag) 2000. |  |  |

|                     |   |  |  |  |
|---------------------|---|--|--|--|
| <b>Experiments:</b> | <input type="checkbox"/> <sup>1</sup> H <input checked="" type="checkbox"/> <sup>13</sup> C{ <sup>1</sup> H} <input type="checkbox"/> <sup>19</sup> F <input type="checkbox"/> <sup>31</sup> P{ <sup>1</sup> H} <input type="checkbox"/> <sup>31</sup> P <input type="checkbox"/> COSY<br><input type="checkbox"/> Other: |  |  |  |
|---------------------|---|--|--|--|

| Assignment(s): | Chemical shift(s)     | Coupling constant(s)     |
|----------------|-----------------------|--------------------------|
|                | Sample spectrum [ppm] | Reference spectrum [ppm] |
|                | Sample spectrum [Hz]  | Reference spectrum [Hz]  |
| a,d            | 11.87                 |                          |
| f              | 42.36                 |                          |
| b,c            | 47.47                 |                          |
| e              | 54.99                 |                          |

|                                |  |
|--------------------------------|--|
| <b>Interpretative comments</b> |  |
|--------------------------------|--|

Compound D


Current Data Parameters  
NAME CWCK1883  
EXPNO 3  
PROCNO 1

F2 - Acquisition Parameters  
Date\_ 20041001  
Time 19.34  
INSTRUM spect  
PROBHD 5 mm TBI 1H/13  
PULPROG zgpg  
TD 131072  
SOLVENT CD2Cl2  
NS 15050  
DS 0  
SWH 50505.051 Hz  
FIDRES 0.385323 Hz  
AQ 1.2976727 sec  
RG 4096  
DE 9.900 usec  
TE 671.7 K  
D1 5.00000000 sec  
d11 0.03000000 sec  
MCREST 0.00000000 sec  
MCWRK 0.01500000 sec

===== CHANNEL f1 =====  
NUC1 13C  
P1 9.00 usec  
PL1 -2.00 dB  
SFO1 125.7533079 MHz

===== CHANNEL f2 =====  
CPDPRG2 waltz16  
NUC2 1H  
PCPD2 85.00 usec  
PL2 0.00 dB  
PL12 23.60 dB  
SFO2 500.0895004 MHz

F2 - Processing Parameters  
SI 131072  
SF 125.7469793 MHz  
WDW EM  
SSB 0  
LB 0.30 Hz  
GB 0  
PC 1.00



$^{13}\text{C}\{^1\text{H}\}$  spectrum of 2-(N,N-Diethylamino)ethylchloride  
(Solvent:  $\text{CH}_2\text{Cl}_2 + \text{CD}_2\text{Cl}_2$ ; Internal reference:  $\text{CH}_2\text{Cl}_2$  at 53.8 ppm)


# NMR METHOD AND ANALYSIS RESULTS

Laboratory code: Prep Lab    Sample code(s): Purity Check    Compound number: E, J

Aliquot codes: CW-CK-1-90-2

Sample: Triethanolamine in CD<sub>2</sub>Cl<sub>2</sub>

Compound structure:




NMR Method name: <sup>1</sup>H NMR

## METHOD DESCRIPTION

|                                   |  |                           |   |
|-----------------------------------|--|---------------------------|---|
| Instrument                        | Bruker Avance 400  |                           |   |
| Manufacturer and Type:            |  |                           |   |
| Frequency:                        | 400.10 MHz   | Temperature control unit: | <input checked="" type="checkbox"/> Yes <input type="checkbox"/> No |
| Probe head:                       | TBI  | Temperature:              | 26.5 °C   |
| Sample tube diameter:             | 5 mm   | Solvent:                  | CD <sub>2</sub> Cl <sub>2</sub>                                     |
| pH of:                            | Sample =   | Blank =                   | Reference =   |
| δ reference reagent in Sample:    | <input checked="" type="checkbox"/> Internal = CD <sub>2</sub> Cl <sub>2</sub> <input type="checkbox"/> External = |                           |   |
| δ reference reagent in Reference: | <input type="checkbox"/> Internal = <input type="checkbox"/> External =  |                           |   |

## ANALYSIS

|  |   |   |  |
|--|---|---|--|
| <input type="checkbox"/> Compared to reference chemical:     |   | Source : <input type="checkbox"/> Own Synthesis <input type="checkbox"/> Commercial   |  |
| <input type="checkbox"/> Compared to library spectrum:       |   | Source : <input type="checkbox"/> OCAD (Code ) <input type="checkbox"/> Other:  |  |
| <input type="checkbox"/> Standard addition:                  |   | Source :  |  |
| <input checked="" type="checkbox"/> Spectral interpretation: |   | E.Pretsch, P. Buhlmann, C. Affolter, <u>Structure Determination of Organic Compounds: Tables of Spectral Data</u> , (Berlin: Springer-Verlag) 2000. |  |
| Experiments:   | <input checked="" type="checkbox"/> <sup>1</sup> H <input type="checkbox"/> <sup>13</sup> C{ <sup>1</sup> H} <input type="checkbox"/> <sup>19</sup> F <input type="checkbox"/> <sup>31</sup> P{ <sup>1</sup> H} <input type="checkbox"/> <sup>31</sup> P <input type="checkbox"/> COSY<br><input type="checkbox"/> Other: |   |  |
| Assignment(s):   | Chemical shift(s)   |   | Coupling constant(s)                         |
|  | Sample spectrum [ppm]   | Reference spectrum [ppm]  | Sample spectrum [Hz] Reference spectrum [Hz] |
| c,d,g  | 2.56  |   | 4.8 (b,c; e,d; h,g)                          |
| b,e,h  | 3.575   |   | 4.8 (c,b; d,e; g,h)                          |
| a,f,i  | 4.87  |   |  |
|  |   |   |  |
| Interpretative comments                                      | Purity = 99.7%  |   |  |



<sup>1</sup>H spectrum of Triethanolamine  
(Solvent: CD<sub>2</sub>Cl<sub>2</sub>; Internal reference: CD<sub>2</sub>Cl<sub>2</sub> at 5.32 ppm)


# NMR METHOD AND ANALYSIS RESULTS

Laboratory code: Prep Lab      Sample code(s): Purity Check      Compound number: E, J

Aliquot codes: CW-CK-1-90-2

Sample: Triethanolamine in CD<sub>2</sub>Cl<sub>2</sub>

Compound structure:




NMR Method name: <sup>13</sup>C{<sup>1</sup>H} NMR

## METHOD DESCRIPTION

|                                   |  |                           |   |
|-----------------------------------|--|---------------------------|---|
| Instrument                        | Bruker DRX 500   |                           |   |
| Manufacturer and Type:            |  |                           |   |
| Frequency:                        | 125.75 MHz   | Temperature control unit: | <input checked="" type="checkbox"/> Yes <input type="checkbox"/> No |
| Probe head:                       | TBI  | Temperature:              | 26.5 °C   |
| Sample tube diameter:             | 5 mm   | Solvent:                  | CH <sub>2</sub> Cl <sub>2</sub> /CD <sub>2</sub> Cl <sub>2</sub>    |
| pH of:                            | Sample =                      Blank =                      Reference =   |                           |   |
| δ reference reagent in Sample:    | <input checked="" type="checkbox"/> Internal = CD <sub>2</sub> Cl <sub>2</sub> <input type="checkbox"/> External = |                           |   |
| δ reference reagent in Reference: | <input type="checkbox"/> Internal = <input type="checkbox"/> External =  |                           |   |

## ANALYSIS

|  |   |   |   |
|--|---|---|---|
| <input type="checkbox"/> Compared to reference chemical:     |   | Source : <input type="checkbox"/> Own Synthesis <input type="checkbox"/> Commercial   |   |
| <input type="checkbox"/> Compared to library spectrum:       |   | Source : <input type="checkbox"/> OCAD (Code        ) <input type="checkbox"/> Other:   |   |
| <input type="checkbox"/> Standard addition:                  |   | Source :  |   |
| <input checked="" type="checkbox"/> Spectral interpretation: |   | E.Pretsch, P. Buhlmann, C. Affolter, <u>Structure Determination of Organic Compounds: Tables of Spectral Data</u> , (Berlin: Springer-Verlag) 2000. |   |
| Experiments:   | <input type="checkbox"/> <sup>1</sup> H <input checked="" type="checkbox"/> <sup>13</sup> C{ <sup>1</sup> H} <input type="checkbox"/> <sup>19</sup> F <input type="checkbox"/> <sup>31</sup> P{ <sup>1</sup> H} <input type="checkbox"/> <sup>31</sup> P <input type="checkbox"/> COSY<br><input type="checkbox"/> Other: |   |   |
| Assignment(s):   | Chemical shift(s)   |   | Coupling constant(s)                              |
|  | Sample spectrum [ppm]   | Reference spectrum [ppm]  | Sample spectrum [Hz]      Reference spectrum [Hz] |
| c,d,g  | 57.39   |   |   |
| b,e,h  | 59.80   |   |   |
|  |   |   |   |
|  |   |   |   |
| Interpretative comments                                      |   |   |   |




# NMR METHOD AND ANALYSIS RESULTS

Laboratory code: Prep Lab      Sample code(s): Purity Check      Compound number: G

**Aliquot codes:** CW-CK-1-90-3

**Sample:** 2-(N,N-Diisopropylamino)ethanol in CD<sub>2</sub>Cl<sub>2</sub>

**Compound structure:**



**NMR Method name:** <sup>1</sup>H NMR


## METHOD DESCRIPTION

|  |  |                                  |   |
|--|--|----------------------------------|---|
| <b>Instrument</b>                        | Bruker Avance 400  |                                  |   |
| <b>Manufacturer and Type:</b>            |  |                                  |   |
| <b>Frequency:</b>                        | 400.10 MHz   | <b>Temperature control unit:</b> | <input checked="" type="checkbox"/> Yes <input type="checkbox"/> No |
| <b>Probe head:</b>                       | TBI  | <b>Temperature:</b>              | 26.5 °C   |
| <b>Sample tube diameter:</b>             | 5 mm   | <b>Solvent:</b>                  | CD <sub>2</sub> Cl <sub>2</sub>                                     |
| <b>pH of:</b>                            | <b>Sample =</b>  | <b>Blank =</b>                   | <b>Reference =</b>  |
| <b>δ reference reagent in Sample:</b>    | <input checked="" type="checkbox"/> Internal = CD <sub>2</sub> Cl <sub>2</sub> <input type="checkbox"/> External = |                                  |   |
| <b>δ reference reagent in Reference:</b> | <input type="checkbox"/> Internal = <input type="checkbox"/> External =  |                                  |   |

## ANALYSIS

|  |   |   |                          |                         |
|--|---|---|--------------------------|-------------------------|
| <input type="checkbox"/> Compared to reference chemical:     |   | Source : <input type="checkbox"/> Own Synthesis <input type="checkbox"/> Commercial   |                          |                         |
| <input type="checkbox"/> Compared to library spectrum:       |   | Source : <input type="checkbox"/> OCAD (Code ) <input type="checkbox"/> Other:  |                          |                         |
| <input type="checkbox"/> Standard addition:                  |   | Source :  |                          |                         |
| <input checked="" type="checkbox"/> Spectral interpretation: |   | E.Pretsch, P. Buhlmann, C. Affolter, <u>Structure Determination of Organic Compounds: Tables of Spectral Data</u> , (Berlin: Springer-Verlag) 2000. |                          |                         |
| Experiments:   | <input checked="" type="checkbox"/> <sup>1</sup> H <input type="checkbox"/> <sup>13</sup> C{ <sup>1</sup> H} <input type="checkbox"/> <sup>19</sup> F <input type="checkbox"/> <sup>31</sup> P{ <sup>1</sup> H} <input type="checkbox"/> <sup>31</sup> P <input type="checkbox"/> COSY<br><input type="checkbox"/> Other: |   |                          |                         |
| Assignment(s):   | Chemical shift(s)   |   | Coupling constant(s)     |                         |
|  | Sample spectrum [ppm]   | Reference spectrum [ppm]  | Sample spectrum [Hz]     | Reference spectrum [Hz] |
| a,b,d,e  | 1.01  |   | 6.8 (c,a; c,b; f,d; f,e) |                         |
| g  | 2.6   |   | 5.6 (h,g)                |                         |
| c,f  | 3.03  |   | 6.8 (a,c; b,c; d,f; e,f) |                         |
| h  | 3.39  |   | 5.6 (g,h)                |                         |
| Interpretative comments                                      | Purity = 99.8%  |   |                          |                         |





<sup>1</sup>H spectrum of 2-(N,N-diisopropylamino)ethanol  
(Solvent: CD<sub>2</sub>Cl<sub>2</sub>; Internal reference: CD<sub>2</sub>Cl<sub>2</sub> at 5.32 ppm)


# NMR METHOD AND ANALYSIS RESULTS

Laboratory code: Prep Lab      Sample code(s): Purity Check      Compound number: G

**Aliquot codes:** CW-CK-1-90-3

**Sample:** 2-(N,N-Diisopropylamino)ethanol in CD<sub>2</sub>Cl<sub>2</sub>

**Compound structure:**



**NMR Method name:** <sup>13</sup>C{<sup>1</sup>H} NMR


## METHOD DESCRIPTION

|  |  |                                  |   |
|--|--|----------------------------------|---|
| <b>Instrument</b>                        | Bruker DRX 500   |                                  |   |
| <b>Manufacturer and Type:</b>            |  |                                  |   |
| <b>Frequency:</b>                        | 125.75 MHz   | <b>Temperature control unit:</b> | <input checked="" type="checkbox"/> Yes <input type="checkbox"/> No |
| <b>Probe head:</b>                       | TBI  | <b>Temperature:</b>              | 26.5 °C   |
| <b>Sample tube diameter:</b>             | 5 mm   | <b>Solvent:</b>                  | CD <sub>2</sub> Cl <sub>2</sub>                                     |
| <b>pH of:</b>                            | <b>Sample =</b>  | <b>Blank =</b>                   | <b>Reference =</b>  |
| <b>δ reference reagent in Sample:</b>    | <input checked="" type="checkbox"/> Internal = CD <sub>2</sub> Cl <sub>2</sub> <input type="checkbox"/> External = |                                  |   |
| <b>δ reference reagent in Reference:</b> | <input type="checkbox"/> Internal = <input type="checkbox"/> External =  |                                  |   |

## ANALYSIS

|  |   |
|--|---|
| <input type="checkbox"/> Compared to reference chemical:     | Source : <input type="checkbox"/> Own Synthesis <input type="checkbox"/> Commercial   |
| <input type="checkbox"/> Compared to library spectrum:       | Source : <input type="checkbox"/> OCAD (Code      ) <input type="checkbox"/> Other:   |
| <input type="checkbox"/> Standard addition:                  | Source :  |
| <input checked="" type="checkbox"/> Spectral interpretation: | E.Pretsch, P. Buhlmann, C. Affolter, <u>Structure Determination of Organic Compounds: Tables of Spectral Data</u> , (Berlin: Springer-Verlag) 2000. |

|                                |   |                                 |                             |                                |
|--------------------------------|---|---------------------------------|-----------------------------|--------------------------------|
| <b>Experiments:</b>            | <input type="checkbox"/> <sup>1</sup> H <input checked="" type="checkbox"/> <sup>13</sup> C{ <sup>1</sup> H} <input type="checkbox"/> <sup>19</sup> F <input type="checkbox"/> <sup>31</sup> P{ <sup>1</sup> H} <input type="checkbox"/> <sup>31</sup> P <input type="checkbox"/> COSY<br><input type="checkbox"/> Other: |                                 |                             |                                |
| <b>Assignment(s):</b>          | <b>Chemical shift(s)</b>  |                                 | <b>Coupling constant(s)</b> |                                |
|                                | <b>Sample spectrum [ppm]</b>  | <b>Reference spectrum [ppm]</b> | <b>Sample spectrum [Hz]</b> | <b>Reference spectrum [Hz]</b> |
| a,b,d,e                        | 20.99   |                                 |                             |                                |
| g                              | 45.78   |                                 |                             |                                |
| c,f                            | 47.98   |                                 |                             |                                |
| h                              | 58.79   |                                 |                             |                                |
| <b>Interpretative comments</b> |   |                                 |                             |                                |




# NMR METHOD AND ANALYSIS RESULTS

Laboratory code: Prep Lab Sample code(s): Purity Check Compound number: H

Aliquot codes: CW-CK-1-91-2

Sample: Ethyldiethanolamine in CD<sub>2</sub>Cl<sub>2</sub>

Compound structure:




NMR Method name: <sup>1</sup>H NMR

## METHOD DESCRIPTION

|                                   |  |                           |   |
|-----------------------------------|--|---------------------------|---|
| Instrument                        | Bruker Avance 400  |                           |   |
| Manufacturer and Type:            |  |                           |   |
| Frequency:                        | 400.10 MHz   | Temperature control unit: | <input checked="" type="checkbox"/> Yes <input type="checkbox"/> No |
| Probe head:                       | TBI  | Temperature:              | 26.5 °C   |
| Sample tube diameter:             | 5 mm   | Solvent:                  | CD <sub>2</sub> Cl <sub>2</sub>                                     |
| pH of:                            | Sample =   | Blank =                   | Reference =   |
| δ reference reagent in Sample:    | <input checked="" type="checkbox"/> Internal = CD <sub>2</sub> Cl <sub>2</sub> <input type="checkbox"/> External = |                           |   |
| δ reference reagent in Reference: | <input type="checkbox"/> Internal = <input type="checkbox"/> External =  |                           |   |

## ANALYSIS

|  |   |   |  |
|--|---|---|--|
| <input type="checkbox"/> Compared to reference chemical:     |   | Source : <input type="checkbox"/> Own Synthesis <input type="checkbox"/> Commercial   |  |
| <input type="checkbox"/> Compared to library spectrum:       |   | Source : <input type="checkbox"/> OCAD (Code ) <input type="checkbox"/> Other:  |  |
| <input type="checkbox"/> Standard addition:                  |   | Source :  |  |
| <input checked="" type="checkbox"/> Spectral interpretation: |   | E.Pretsch, P. Buhlmann, C. Affolter, <u>Structure Determination of Organic Compounds: Tables of Spectral Data</u> , (Berlin: Springer-Verlag) 2000. |  |
| Experiments:   | <input checked="" type="checkbox"/> <sup>1</sup> H <input type="checkbox"/> <sup>13</sup> C{ <sup>1</sup> H} <input type="checkbox"/> <sup>19</sup> F <input type="checkbox"/> <sup>31</sup> P{ <sup>1</sup> H} <input type="checkbox"/> <sup>31</sup> P <input type="checkbox"/> COSY<br><input type="checkbox"/> Other: |   |  |
| Assignment(s):   | Chemical shift(s)   |   | Coupling constant(s)                         |
|  | Sample spectrum [ppm]   | Reference spectrum [ppm]  | Sample spectrum [Hz] Reference spectrum [Hz] |
| d  | 1.01  |   | 7.2 (c,d)                                    |
| b,e,c  | 2.59  |   |  |
| a,f  | 3.55  |   | 5.6 (b,a; e,f)                               |
|  |   |   |  |
| Interpretative comments                                      | Purity = 99.6%  |   |  |




# NMR METHOD AND ANALYSIS RESULTS

Laboratory code: Prep Lab    Sample code(s): Purity Check    Compound number: H

**Aliquot codes:** CW-CK-1-91-2

**Sample:** Ethyldiethanolamine in CD<sub>2</sub>Cl<sub>2</sub>

**Compound structure:**



**NMR Method name:** <sup>13</sup>C{<sup>1</sup>H} NMR

## METHOD DESCRIPTION


|  |  |                                  |   |
|--|--|----------------------------------|---|
| <b>Instrument</b>                        | Bruker DRX 500   |                                  |   |
| <b>Manufacturer and Type:</b>            |  |                                  |   |
| <b>Frequency:</b>                        | 125.75 MHz   | <b>Temperature control unit:</b> | <input checked="" type="checkbox"/> Yes <input type="checkbox"/> No |
| <b>Probe head:</b>                       | TBI  | <b>Temperature:</b>              | 26.5 °C   |
| <b>Sample tube diameter:</b>             | 5 mm   | <b>Solvent:</b>                  | CD <sub>2</sub> Cl <sub>2</sub>                                     |
| <b>pH of:</b>                            | <b>Sample =</b>  | <b>Blank =</b>                   | <b>Reference =</b>  |
| <b>δ reference reagent in Sample:</b>    | <input checked="" type="checkbox"/> Internal = CD <sub>2</sub> Cl <sub>2</sub> <input type="checkbox"/> External = |                                  |   |
| <b>δ reference reagent in Reference:</b> | <input type="checkbox"/> Internal = <input type="checkbox"/> External =  |                                  |   |

## ANALYSIS

|  |   |
|--|---|
| <input type="checkbox"/> Compared to reference chemical:     | Source : <input type="checkbox"/> Own Synthesis <input type="checkbox"/> Commercial   |
| <input type="checkbox"/> Compared to library spectrum:       | Source : <input type="checkbox"/> OCAD (Code      ) <input type="checkbox"/> Other:   |
| <input type="checkbox"/> Standard addition:                  | Source :  |
| <input checked="" type="checkbox"/> Spectral interpretation: | E.Pretsch, P. Buhlmann, C. Affolter, <u>Structure Determination of Organic Compounds: Tables of Spectral Data</u> , (Berlin: Springer-Verlag) 2000. |

|                     |   |
|---------------------|---|
| <b>Experiments:</b> | <input type="checkbox"/> <sup>1</sup> H <input checked="" type="checkbox"/> <sup>13</sup> C{ <sup>1</sup> H} <input type="checkbox"/> <sup>19</sup> F <input type="checkbox"/> <sup>31</sup> P{ <sup>1</sup> H} <input type="checkbox"/> <sup>31</sup> P <input type="checkbox"/> COSY<br><input type="checkbox"/> Other: |
|---------------------|---|

| Assignment(s):                 | Chemical shift(s)  | Coupling constant(s)  |                          |                      |                         |   |       |  |  |   |       |  |  |     |       |  |  |     |       |  |  |  |
|--------------------------------|--|-----------------------|--------------------------|----------------------|-------------------------|---|-------|--|--|---|-------|--|--|-----|-------|--|--|-----|-------|--|--|--|
|                                | <table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <th style="width: 30%;">Sample spectrum [ppm]</th> <th style="width: 30%;">Reference spectrum [ppm]</th> <th style="width: 30%;">Sample spectrum [Hz]</th> <th style="width: 30%;">Reference spectrum [Hz]</th> </tr> <tr> <td>d</td> <td>11.30</td> <td></td> <td></td> </tr> <tr> <td>c</td> <td>48.12</td> <td></td> <td></td> </tr> <tr> <td>b,e</td> <td>55.53</td> <td></td> <td></td> </tr> <tr> <td>a,f</td> <td>59.49</td> <td></td> <td></td> </tr> </table> | Sample spectrum [ppm] | Reference spectrum [ppm] | Sample spectrum [Hz] | Reference spectrum [Hz] | d | 11.30 |  |  | c | 48.12 |  |  | b,e | 55.53 |  |  | a,f | 59.49 |  |  |  |
| Sample spectrum [ppm]          | Reference spectrum [ppm]   | Sample spectrum [Hz]  | Reference spectrum [Hz]  |                      |                         |   |       |  |  |   |       |  |  |     |       |  |  |     |       |  |  |  |
| d                              | 11.30  |                       |                          |                      |                         |   |       |  |  |   |       |  |  |     |       |  |  |     |       |  |  |  |
| c                              | 48.12  |                       |                          |                      |                         |   |       |  |  |   |       |  |  |     |       |  |  |     |       |  |  |  |
| b,e                            | 55.53  |                       |                          |                      |                         |   |       |  |  |   |       |  |  |     |       |  |  |     |       |  |  |  |
| a,f                            | 59.49  |                       |                          |                      |                         |   |       |  |  |   |       |  |  |     |       |  |  |     |       |  |  |  |
| <b>Interpretative comments</b> |  |                       |                          |                      |                         |   |       |  |  |   |       |  |  |     |       |  |  |     |       |  |  |  |




# NMR METHOD AND ANALYSIS RESULTS

Laboratory code: Prep Lab      Sample code(s): Purity Check      Compound number: I

Aliquot codes: CW-CK-1-91-1

Sample: Methyldiethanolamine in CD<sub>2</sub>Cl<sub>2</sub>

Compound structure:



NMR Method name: <sup>1</sup>H NMR

## METHOD DESCRIPTION



|                                   |  |                           |   |
|-----------------------------------|--|---------------------------|---|
| Instrument                        | Bruker Avance 400  |                           |   |
| Manufacturer and Type:            |  |                           |   |
| Frequency:                        | 400.10 MHz   | Temperature control unit: | <input checked="" type="checkbox"/> Yes <input type="checkbox"/> No |
| Probe head:                       | TBI  | Temperature:              | 26.5 °C   |
| Sample tube diameter:             | 5 mm   | Solvent:                  | CD <sub>2</sub> Cl <sub>2</sub>                                     |
| pH of:                            | Sample =   | Blank =                   | Reference =   |
| δ reference reagent in Sample:    | <input checked="" type="checkbox"/> Internal = CD <sub>2</sub> Cl <sub>2</sub> <input type="checkbox"/> External = |                           |   |
| δ reference reagent in Reference: | <input type="checkbox"/> Internal = <input type="checkbox"/> External =  |                           |   |

## ANALYSIS

|  |   |   |  |
|--|---|---|--|
| <input type="checkbox"/> Compared to reference chemical:     |   | Source : <input type="checkbox"/> Own Synthesis <input type="checkbox"/> Commercial   |  |
| <input type="checkbox"/> Compared to library spectrum:       |   | Source : <input type="checkbox"/> OCAD (Code ) <input type="checkbox"/> Other:  |  |
| <input type="checkbox"/> Standard addition:                  |   | Source :  |  |
| <input checked="" type="checkbox"/> Spectral interpretation: |   | E.Pretsch, P. Buhlmann, C. Affolter, <u>Structure Determination of Organic Compounds: Tables of Spectral Data</u> , (Berlin: Springer-Verlag) 2000. |  |
| Experiments:   | <input checked="" type="checkbox"/> <sup>1</sup> H <input type="checkbox"/> <sup>13</sup> C{ <sup>1</sup> H} <input type="checkbox"/> <sup>19</sup> F <input type="checkbox"/> <sup>31</sup> P{ <sup>1</sup> H} <input type="checkbox"/> <sup>31</sup> P <input type="checkbox"/> COSY<br><input type="checkbox"/> Other: |   |  |
| Assignment(s):   | Chemical shift(s)   |   | Coupling constant(s)                         |
|  | Sample spectrum [ppm]   | Reference spectrum [ppm]  | Sample spectrum [Hz] Reference spectrum [Hz] |
| c  | 2.25  |   |  |
| b,d  | 2.51  |   | 5.6 (a,b; d,e)                               |
| a,e  | 3.57  |   | 5.6 (b,a; e,d)                               |
| f,g  | 4.13  |   |  |
| Interpretative comments                                      | Purity = 99.7%  |   |  |



Current Data Parameters  
 NAME CWCK1911  
 EXPNO 1  
 PROCNO 1  
 F2 - Acquisition Parameters  
 Date\_ 20041001  
 Time 14.46  
 INSTRUM spect  
 PROHD 5 mm TBI 1H-BB  
 PULPROG zg  
 TD 131072  
 SOLVENT CD2Cl2  
 NS 64  
 DS 0  
 SWH 6009.615 Hz  
 FIDRES 0.045850 Hz  
 AQ 10.9053240 sec  
 RG 18  
 DW 83.200 usec  
 DE 6.00 usec  
 TE 673.0 K  
 D1 1.00000000 sec  
 MCREST 0.00000000 sec  
 MCWRK 0.01500000 sec  
 ===== CHANNEL f1 =====  
 NUC1 1H  
 P1 2.27 usec  
 PL1 0.00 dB  
 SFO1 400.0980456 MHz  
 F2 - Processing parameters  
 SI 131072  
 SF 400.0960180 MHz  
 WDW no  
 SSB 0  
 LB 0.00 Hz  
 GB 0  
 PC 1.00



<sup>1</sup>H spectrum of Methylmethanolamine  
 (Solvent: CD<sub>2</sub>Cl<sub>2</sub>; Internal reference: CD<sub>2</sub>Cl<sub>2</sub> at 5.32 ppm)


## NMR METHOD AND ANALYSIS RESULTS

Laboratory code: Prep Lab      Sample code(s): Purity Check      Compound number: I

**Aliquot codes:** CW-CK-1-91-1

**Sample:** Methyldiethanolamine in CD<sub>2</sub>Cl<sub>2</sub>

**Compound structure:**



**NMR Method name:** <sup>13</sup>C{<sup>1</sup>H} NMR

### METHOD DESCRIPTION

|  |  |                                  |   |
|--|--|----------------------------------|---|
| <b>Instrument</b>                        | Bruker DRX 500   |                                  |   |
| <b>Manufacturer and Type:</b>            |  |                                  |   |
| <b>Frequency:</b>                        | 125.75 MHz   | <b>Temperature control unit:</b> | <input checked="" type="checkbox"/> Yes <input type="checkbox"/> No |
| <b>Probe head:</b>                       | TBI  | <b>Temperature:</b>              | 26.5 °C   |
| <b>Sample tube diameter:</b>             | 5 mm   | <b>Solvent:</b>                  | CD <sub>2</sub> Cl <sub>2</sub>                                     |
| <b>pH of:</b>                            | <b>Sample =</b>  | <b>Blank =</b>                   | <b>Reference =</b>  |
| <b>δ reference reagent in Sample:</b>    | <input checked="" type="checkbox"/> Internal = CD <sub>2</sub> Cl <sub>2</sub> <input type="checkbox"/> External = |                                  |   |
| <b>δ reference reagent in Reference:</b> | <input type="checkbox"/> Internal = <input type="checkbox"/> External =  |                                  |   |


### ANALYSIS

|  |   |
|--|---|
| <input type="checkbox"/> Compared to reference chemical:     | Source : <input type="checkbox"/> Own Synthesis <input type="checkbox"/> Commercial   |
| <input type="checkbox"/> Compared to library spectrum:       | Source : <input type="checkbox"/> OCAD (Code      ) <input type="checkbox"/> Other:   |
| <input type="checkbox"/> Standard addition:                  | Source :  |
| <input checked="" type="checkbox"/> Spectral interpretation: | E.Pretsch, P. Buhlmann, C. Affolter, <u>Structure Determination of Organic Compounds: Tables of Spectral Data</u> , (Berlin: Springer-Verlag) 2000. |

|                     |   |
|---------------------|---|
| <b>Experiments:</b> | <input type="checkbox"/> <sup>1</sup> H <input checked="" type="checkbox"/> <sup>13</sup> C{ <sup>1</sup> H} <input type="checkbox"/> <sup>19</sup> F <input type="checkbox"/> <sup>31</sup> P{ <sup>1</sup> H} <input type="checkbox"/> <sup>31</sup> P <input type="checkbox"/> COSY<br><input type="checkbox"/> Other: |
|---------------------|---|

| Assignment(s):        | Chemical shift(s)  | Coupling constant(s)  |                          |  |                      |                         |
|-----------------------|--|-----------------------|--------------------------|--|----------------------|-------------------------|
|                       | <table style="width: 100%; border-collapse: collapse;"> <tr> <th style="width: 50%;">Sample spectrum [ppm]</th> <th style="width: 50%;">Reference spectrum [ppm]</th> </tr> </table> | Sample spectrum [ppm] | Reference spectrum [ppm] | <table style="width: 100%; border-collapse: collapse;"> <tr> <th style="width: 50%;">Sample spectrum [Hz]</th> <th style="width: 50%;">Reference spectrum [Hz]</th> </tr> </table> | Sample spectrum [Hz] | Reference spectrum [Hz] |
| Sample spectrum [ppm] | Reference spectrum [ppm]   |                       |                          |  |                      |                         |
| Sample spectrum [Hz]  | Reference spectrum [Hz]  |                       |                          |  |                      |                         |
| c                     | 42.29  |                       |                          |  |                      |                         |
| a,e                   | 59.39  |                       |                          |  |                      |                         |
| b,d                   | 60.03  |                       |                          |  |                      |                         |
|                       |  |                       |                          |  |                      |                         |

|                                |  |
|--------------------------------|--|
| <b>Interpretative comments</b> |  |
|--------------------------------|--|





**ORGANISATION FOR THE PROHIBITION  
OF CHEMICAL WEAPONS**

**Report of the  
Sixteenth Official OPCW  
Proficiency Test**

***Part III: Qualitative Analysis***

Laboratory code: Sample Prep Lab

|                          |
|--------------------------|
| <b>TABLE OF CONTENTS</b> |
|--------------------------|

|  | <b>Page no.</b> |
|--|-----------------|
| <b>Cover page</b> .....                              | 57              |
| <b>Table of contents</b> .....                       | 58              |
| <b>Summary</b>                                       |                 |
| Participating Laboratory .....                       | 60              |
| Quality Assurance/Quality Control (QA/QC) .....      | 61              |
| Names and structures of all reported compounds ..... | 64              |
| Analytical techniques: <b>O/31</b> .....             | 65              |
| <b>L/31</b> .....                                    | 66              |
| <b>S/31</b> .....                                    | 67              |
| <b>Organic Sample</b>                                |                 |
| Sample preparation description/flowchart .....       | 68              |
| Compound number 1                                    |                 |
| Analytical technique 1: GC/MS-EI .....               | 70              |
| Analytical technique 2: GC/MS-CI .....               | 73              |
| Compound number 2                                    |                 |
| Analytical technique 1: GC/MS-EI .....               | 76              |
| Analytical technique 2: GC/MS-CI .....               | 79              |
| Compound number 3                                    |                 |
| Analytical technique 1: GC/MS-EI .....               | 82              |
| Analytical technique 2: GC/MS-CI .....               | 85              |
| Compound number 4                                    |                 |
| Analytical technique 1: GC/MS-EI .....               | 88              |
| Analytical technique 2: GC/MS-CI .....               | 91              |
| <b>Liquid Sample</b>                                 |                 |
| Sample preparation description/flowchart .....       | 94              |
| Compound number 5                                    |                 |
| Analytical technique 1: GC/MS-EI .....               | 96              |
| Analytical technique 2: GC/MS-CI .....               | 99              |
| Compound number 6                                    |                 |
| Analytical technique 1: GC/MS-EI .....               | 102             |
| Analytical technique 2: GC/MS-CI .....               | 105             |

**Soil Sample**

|  |     |
|--|-----|
| Sample preparation description/flowchart ..... | 108 |
| Compound number 7                              |     |
| Analytical technique 1: GC/MS-EI .....         | 110 |
| Analytical technique 2: GC/MS-CI .....         | 113 |
| Compound number 8                              |     |
| Analytical technique 1: GC/MS-EI .....         | 116 |
| Analytical technique 2: GC/MS-CI .....         | 119 |
| Compound number 9                              |     |
| Analytical technique 1: GC/MS-EI .....         | 122 |
| Analytical technique 2: GC/MS-CI .....         | 125 |
| <b>Comments</b> .....                          | 128 |

## SUMMARY: PARTICIPATING LABORATORY

### 1. Participating laboratory

|  |  |
|--|--|
| Laboratory code  | Sample Prep Lab, aliquot #31                                     |
| Name of the laboratory/institute participating in the test | Lawrence Livermore National Laboratory                           |
| Contact person   | Mr. Armando Alcaraz  |
| Address  | PO Box 808, M/S L-178<br>7000 East Avenue<br>Livermore, CA 94551 |
| Telephone number   | 925-423-6889   |
| Fax number   | 925-423-9014   |
| Email address  | Alcaraz1@llnl.gov  |

### 2. Analysts and authentication

|    | Name                       | Title                  | Pages*       | Date**        | Signature** |
|----|----------------------------|------------------------|--------------|---------------|-------------|
| 1  | Mr. Armando Alcaraz        | Principal Investigator | All          | Nov. 30, 2004 |             |
| 2  | Dr. Hugh Gregg             | Co-PI, Senior Chemist  | All          | Nov. 30, 2004 |             |
| 3  | Dr. Carolyn Koester        | Research Scientist     | Prep & GC/MS | Nov. 30, 2004 |             |
| 4  | Dr. Phil Pagoria           | Research Scientist     | Synthesis    | Nov. 30, 2004 |             |
| 5  | Dr. Andrew Vance           | Research Scientist     | Synthesis    | Nov. 30, 2004 |             |
| 6  | Dr. Robert Maxwell         | Research Scientist     | NMR analysis | Nov. 30, 2004 |             |
| 7  | Dr. Sarah Chinn            | Research Scientist     | NMR analysis | Nov. 30, 2004 |             |
| 8  | Mr. Rich Whipple           | Scientist              | Sample Prep  | Nov. 30, 2004 |             |
| 9  | Ms. Tuijauna Mitchell-Hall | QA Manager             | QA/QC        | Nov. 30, 2004 |             |
| 10 |                            |                        |              |               |             |
| 11 |                            |                        |              |               |             |
| 12 |                            |                        |              |               |             |

\* Page numbers defining the responsibility area of the analyst;

\*\* Date and signature of the responsible analyst;

## SUMMARY: QUALITY ASSURANCE / QUALITY CONTROL (QA/QC)

### 1. Status of the laboratory (tick where applicable)


- ☒ Accreditation accepted:  
Year 2001 Accreditation body: American Association for Laboratory Accreditation  
Scope of accreditation: Chemical
- ☐ Accreditation planned/pending:  
Target year \_\_\_\_\_ Accreditation body: \_\_\_\_\_  
Scope of accreditation: \_\_\_\_\_
- ☐ Not accredited.

### 2. Quality system (tick where applicable)

- ☒ Described in a Quality Assurance Manual/Handbook. Quality system in accordance with:  
☐ ISO 900\_\_\_\_\_, ☐ EN 4500\_\_\_\_\_, ☒ ISO Guide 17025, ☐ Other: \_\_\_\_\_
- ☐ No quality system. Please, fill in question number 3.

### 3. QA/QC Summary (Summary of the applied quality assurance and quality control (QA/QC) procedures concerning sample preparation, calibration, and analysis. Requested only from laboratories without a quality system).





THE AMERICAN  
ASSOCIATION  
FOR LABORATORY  
ACCREDITATION

## ACCREDITED LABORATORY

A2LA has accredited

**LAWRENCE LIVERMORE NATIONAL  
LABORATORY**  
Livermore, CA


for technical competence in the field of

### Chemical Testing

The accreditation covers the specific tests and types of tests listed on the agreed scope of accreditation. This laboratory meets the requirements of ISO/IEC 17025 - 1999 "General Requirements for the Competence of Testing and Calibration Laboratories" and any additional program requirements in the identified field of testing.

Presented this 12<sup>th</sup> day of April 2004.



  
President  
For the Accreditation Council  
Certificate Number 1914-01  
Valid to February 28, 2006

For tests or types of tests to which this accreditation applies,  
please refer to the laboratory's Chemical Scope of Accreditation.



**American Association for Laboratory Accreditation**SCOPE OF ACCREDITATION TO ISO/IEC 17025-1999

LAWRENCE LIVERMORE NATIONAL LABORATORY  
FORENSIC SCIENCE CENTER – OPCW PROJECT  
7000 East Avenue Mailstop L-178  
Livermore, CA 94550  
Armando Alcaraz Phone: 925 423 6889

**CHEMICAL**

Valid To: February 28, 2006

Certificate Number: 1914-01

In recognition of the successful completion of the A2LA evaluation process, accreditation is granted to this laboratory to perform the following types of Qualitative Tests for Chemicals related to Chemical Warfare Convention (CWC) in unknown samples:

Sample Preparation BB-SP5, BB-SP6, BB-SP8, BB-SP9,  
LL-SP1, LL-SP2, LL-SP3, LL-SP4

TestProceduresSpectroscopy

Nuclear magnetic resonance

BB-NMR1

Capillary zone electrophoresis / UV detection

LL-CE1

Gas chromatography / Fourier Transform  
Infrared Spectrometry

BB-IR1

Chromatography

Gas chromatography / Element Specific Detectors

BB-GC1

Gas chromatography / Mass Spectrometry

BB-MS1

Liquid chromatography / Atmospheric pressure  
Chemical Ionization / Mass Spectrometry

BB-MS4

Liquid Chromatography / Electrospray Ionization  
Mass Spectrometry

LL-MS1

Chain of Custody for Laboratory

LL-CC1

Work Instructions for the Preparation of Test Samples  
for OPCW Proficiency TestsQDOC/LAB/WI/ PT2  
(Minus Section 11 Confirmation)

(A2LA Cert. No. 1914-01) Revised 04/15/2004










Page 1 of 1

5301 Buckeystown Pike, Suite 350 • Frederick, MD 21704-8373 • Phone: 301-644 3248 • Fax: 301-662 2974



## SUMMARY: NAMES AND STRUCTURES OF ALL REPORTED COMPOUNDS

Laboratory code: 31

| Sam. code | Cmpd. no* | Compound name                          | Chemical Abstract number | Compound Structure  | Molecular formula                              | Schedule number | Comments ** |
|-----------|-----------|--|--------------------------|---|--|-----------------|-------------|
| O         | 1         | Trichloronitromethane                  | 76-06-2                  |    | CCl <sub>3</sub> NO <sub>2</sub>               | 3.A.04          |             |
| O         | 2         | 2-(N-Ethyl-N-propylamino)ethylchloride |                          |    | C <sub>7</sub> H <sub>16</sub> ClN             | 2.B.10          |             |
| O         | 3         | 2-(N,N-Diisopropylamino)ethylchloride  | 96-79-7                  |    | C <sub>8</sub> H <sub>18</sub> ClN             | 2.B.10          |             |
| O         | 4         | 2-(N,N-Diethylamino)ethylchloride      | 100-35-6                 |    | C <sub>6</sub> H <sub>14</sub> ClN             | 2.B.10          |             |
| L         | 5         | Triethanolamine                        | 102-71-6                 |    | C <sub>6</sub> H <sub>15</sub> NO <sub>3</sub> | 3.B.17          |             |
| L         | 6         | 2-(N,N-Diisopropylamino)ethanol        | 96-80-0                  |   | C <sub>8</sub> H <sub>19</sub> NO              | 2.B.11          |             |
| S         | 7         | Ethyldiethanolamine                    | 139-87-7                 |  | C <sub>6</sub> H <sub>15</sub> NO <sub>2</sub> | 3.B.15          |             |
| S         | 8         | Methyldiethanolamine                   | 105-59-9                 |  | C <sub>5</sub> H <sub>13</sub> NO <sub>2</sub> | 3.B.16          |             |
| S         | 9         | Triethanolamine                        | 102-71-6                 |  | C <sub>6</sub> H <sub>15</sub> NO <sub>3</sub> | 3.B.17          |             |

\* Compound number defined by the participating laboratory and used throughout the report for the reported compound.

\*\* Explanation for the reporting of non-scheduled compounds, details can be added in the comment section of the report.

Note: There must be an unbroken chain of evidence linking each reported chemical to the original sample from which an aliquot was prepared and analyzed for the identification of this chemical.

## SUMMARY: ANALYTICAL TECHNIQUES

Laboratory code: 31      Sample code(s): O/31

| Compound number* | Compound name                          | Compound analysed as  | Analytical technique | Method name | Method page no. | Aliquot name                 |
|------------------|--|---|----------------------|-------------|-----------------|------------------------------|
| 1                | Trichloronitromethane                  | <input checked="" type="checkbox"/> original compound<br><input type="checkbox"/> methylated<br><input type="checkbox"/> silylated<br><input type="checkbox"/> other: _____ | GC/MS-EI<br>GC/MS-CI | CW<br>CW-CI | 70<br>73        | CW-1-130-7-O<br>CW-1-130-7-O |
| 2                | 2-(N-Ethyl-N-propylamino)ethylchloride | <input checked="" type="checkbox"/> original compound<br><input type="checkbox"/> methylated<br><input type="checkbox"/> silylated<br><input type="checkbox"/> other: _____ | GC/MS-EI<br>GC/MS-CI | CW<br>CW-CI | 76<br>79        | CW-1-130-7-O<br>CW-1-130-7-O |
| 3                | 2-(N,N-Diisopropylamino)ethylchloride  | <input checked="" type="checkbox"/> original compound<br><input type="checkbox"/> methylated<br><input type="checkbox"/> silylated<br><input type="checkbox"/> other: _____ | GC/MS-EI<br>GC/MS-CI | CW<br>CW-CI | 82<br>85        | CW-1-130-7-O<br>CW-1-130-7-O |
| 4                | 2-(N,N-Diethylamino)ethylchloride      | <input checked="" type="checkbox"/> original compound<br><input type="checkbox"/> methylated<br><input type="checkbox"/> silylated<br><input type="checkbox"/> other: _____ | GC/MS-EI<br>GC/MS-CI | CW<br>CW-CI | 88<br>91        | CW-1-130-7-O<br>CW-1-130-7-O |

\* Compound number defined by the participating laboratory (see Summary: Names and Structures of All Reported Compounds);

## SUMMARY: ANALYTICAL TECHNIQUES

Laboratory code: 31      Sample code(s): L/31

| Compound number* | Compound name                   | Compound analysed as                          | Analytical technique | Method name | Method page no. | Aliquot name |
|------------------|---------------------------------|---|----------------------|-------------|-----------------|--------------|
| 5                | Triethanolamine                 | <input type="checkbox"/> original compound    | GC/MS-EI             | TMS_A       | 96              | CW-1-131-5-L |
|                  |                                 | <input type="checkbox"/> methylated           | GC/MS-CI             | CW-CI-TM    | 99              | CW-1-131-5-L |
|                  |                                 | <input checked="" type="checkbox"/> silylated |                      |             |                 |              |
|                  |                                 | <input type="checkbox"/> other: _____         |                      |             |                 |              |
| 6                | 2-(N,N-Diisopropylamino)ethanol | <input type="checkbox"/> original compound    | GC/MS-EI             | TMS_A       | 102             | CW-1-131-4-L |
|                  |                                 | <input type="checkbox"/> methylated           | GC/MS-CI             | CW-CI-TM    | 105             | CW-1-131-4-L |
|                  |                                 | <input checked="" type="checkbox"/> silylated |                      |             |                 |              |
|                  |                                 | <input type="checkbox"/> other: _____         |                      |             |                 |              |
|                  |                                 | <input type="checkbox"/> original compound    |                      |             |                 |              |
|                  |                                 | <input type="checkbox"/> methylated           |                      |             |                 |              |
|                  |                                 | <input type="checkbox"/> silylated            |                      |             |                 |              |
|                  |                                 | <input type="checkbox"/> other: _____         |                      |             |                 |              |

\* Compound number defined by the participating laboratory (see Summary: Names and Structures of All Reported Compounds);

## SUMMARY: ANALYTICAL TECHNIQUES

Laboratory code: 31    Sample code(s): S/31


| Compound number* | Compound name        | Compound analysed as  | Analytical technique | Method name       | Method page no. | Aliquot name                 |
|------------------|----------------------|---|----------------------|-------------------|-----------------|------------------------------|
| 7                | Ethyldiethanolamine  | <input type="checkbox"/> original compound<br><input type="checkbox"/> methylated<br><input checked="" type="checkbox"/> silylated<br><input type="checkbox"/> other: _____ | GC/MS-EI<br>GC/MS-CI | TMS_A<br>CW-CI-TM | 110<br>113      | CW-1-131-2-S<br>CW-1-131-2-S |
| 8                | Methyldiethanolamine | <input type="checkbox"/> original compound<br><input type="checkbox"/> methylated<br><input checked="" type="checkbox"/> silylated<br><input type="checkbox"/> other: _____ | GC/MS-EI<br>GC/MS-CI | TMS_A<br>CW-CI-TM | 116<br>119      | CW-1-131-2-S<br>CW-1-131-2-S |
| 9                | Triethanolamine      | <input type="checkbox"/> original compound<br><input type="checkbox"/> methylated<br><input checked="" type="checkbox"/> silylated<br><input type="checkbox"/> other: _____ | GC/MS-EI<br>GC/MS-CI | TMS_A<br>CW-CI-TM | 122<br>125      | CW-1-131-2-S<br>CW-1-131-2-S |
|                  |                      | <input type="checkbox"/> original compound<br><input type="checkbox"/> methylated<br><input type="checkbox"/> silylated<br><input type="checkbox"/> other: _____            |                      |                   |                 |                              |

\* Compound number defined by the participating laboratory (see Summary: Names and Structures of All Reported Compounds);

**SAMPLE PREPARATION DESCRIPTION**Laboratory code: 31 Sample code(s): O/31 Sample Blank code: OB/31**1. Sample preparation**

| <b>Sample/<br/>Aliquot Code</b> | <i>Specification of Sample/<br/>Type of Sample Preparation</i> | <b>Amount/<br/>Volume</b> | <i>Sample Preparation Procedures</i> | <b>End<br/>Volume</b> | <b>Resulting<br/>Aliquot Code</b> | <b>Analytical<br/>Technique(s)</b> |
|---------------------------------|--|---------------------------|--------------------------------------|-----------------------|-----------------------------------|------------------------------------|
| O/31<br>OB/31                   | 1 mL aliquot of each   | 1 mL                      | (none)                               | 1 mL                  | CW-1-130-7-O<br>CW-1-130-8-OB     | GC/MS-EI<br>GC/MS-CI               |
|                                 |  |                           |                                      |                       |                                   |                                    |
|                                 |  |                           |                                      |                       |                                   |                                    |
|                                 |  |                           |                                      |                       |                                   |                                    |
|                                 |  |                           |                                      |                       |                                   |                                    |
|                                 |  |                           |                                      |                       |                                   |                                    |
|                                 |  |                           |                                      |                       |                                   |                                    |
|                                 |  |                           |                                      |                       |                                   |                                    |
|                                 |  |                           |                                      |                       |                                   |                                    |

**2. Additional information**



Note: This flowchart is for visualization only; see the preceding sample preparation description page for sample aliquot numbers



# GC-EI-MS TECHNIQUE

## METHOD AND ANALYSIS DESCRIPTION

Laboratory code: 31 Sample code(s): O/31, OB/31 Compound number: 1

**Aliquot codes:**

**Sample:** CW-1-130-7-O

**Blank:** CW-1-130-8-OB

**GC-EI-MS Method name:** CW

### METHOD DESCRIPTION

|   |  |   |   |
|---|--|---|---|
| <b>Instrument Manufacturer and Type:</b>    | Agilent 6890/5973 GC/MSD   |   |   |
| <b>Carrier gas:</b>                         | <input checked="" type="checkbox"/> He   | <input type="checkbox"/> N <sub>2</sub>           | <input type="checkbox"/> H <sub>2</sub> <input type="checkbox"/> Other: |
| <b>Flow rate:</b>                           | <input type="checkbox"/> ml/min  | <input checked="" type="checkbox"/> 32 cm/s       |   |
| <b>Flow control:</b>                        | <input type="checkbox"/> Constant Pressure   | <input checked="" type="checkbox"/> Constant Flow |   |
| <b>Injection mode:</b>                      | <input type="checkbox"/> Split → Split ratio =<br><input checked="" type="checkbox"/> Splitless → Splitless time = 0.75 min. |   |   |
| <b>Injector temperature:</b>                | 250 °C   |   |   |
| <b>Column phase:</b>                        | 5% diphenyl 95% dimethyl polysiloxane  |   |   |
| <b>Column Length x ID x Film thickness:</b> | 30 m x 0.25 mm x 0.25 µm   |   |   |
| <b>GC temperature programme:</b>            | 40 °C (3 min), 8 °C/min, 300 °C (3 min)  |   |   |
| <b>Solvent delay time:</b>                  | 3 min  | <b>Scan range:</b>                                | 30-600 m/z  |
| <b>Electron energy:</b>                     | 70 eV  | <b>Scan time:</b>                                 | 0.7 s   |
| <b>Ionisation polarity:</b>                 | <input checked="" type="checkbox"/> Positive<br><input type="checkbox"/> Negative  | <b>Mass resolution:</b>                           | 0.7 u   |
| <b>Comments:</b>                            |  |   |   |

### ANALYSIS

|  |   |
|--|---|
| <b>Compound identified as:</b>   | <input checked="" type="checkbox"/> Original compound<br><input type="checkbox"/> Methyl ester derivative<br><input type="checkbox"/> TBDMS (t-Butyldimethylsilyl) derivative<br><input type="checkbox"/> TMS (Trimethylsilyl) derivative<br><input type="checkbox"/> Other derivative: |
| <b>Retention parameter used for (peak) identification:</b>                       | <input checked="" type="checkbox"/> Retention time (Rt)<br><input type="checkbox"/> Scan number   |
| <input checked="" type="checkbox"/> Compared to reference chemical:              | Source : <input checked="" type="checkbox"/> Own Synthesis <input type="checkbox"/> Commercial  |
| <input type="checkbox"/> Compared to library spectrum:                           | Source : <input type="checkbox"/> OCAD (code: ) <input type="checkbox"/> NIST<br><input type="checkbox"/> Wiley <input type="checkbox"/> Own <input type="checkbox"/> Other:  |
| <input type="checkbox"/> Not compared to reference chemical or library spectrum: | Intense ions in spectrum are explained; interpretation is supported by the spectral information derived from closely related chemical(s):   |
| <b>Comments:</b>   |   |



File : C:\DATA\16\J0003493.D

Acquired: 28 Oct 2004 9:55

Method : CW

Sample : CW-1-130-8-OB; 1 uL direct inject

Misc info: J&W DB-5MS; 30 m, 0.25 mm, 0.25 um





File : C:\DATA\16\J0003494.D

Acquired: 28 Oct 2004 10:43

Method : CW

Sample : CW-1-130-7-O; 1 uL direct inject

Misc info: J&W DB-5MS; 30 m, 0.25 mm, 0.25 um





File : C:\DATA\16\J0003495.D

Acquired: 28 Oct 2004 11:30

Method : CW

Sample : CW-CK-1-124-2; 1 uL direct inject A

Misc info: J&W DB-5MS; 30 m, 0.25 mm, 0.25 um




EI chromatograms supporting identification of compound 1; TIC on left; EIC (m/z 117) on right.

Top: Chromatograms of Organic blank, aliquot CW-1-130-8-OB from OB/31.


Center: Chromatograms of Organic sample, aliquot CW-1-130-7-O from O/31, retention time 5.49 min.

Bottom: Chromatograms of authentic reference standard of Trichloronitromethane, retention time 5.49 min.

File : C:\DATA\16\J0003494.D  
Acquired : 28 Oct 2004 10:43 using AcqMethod CW  
Sample Name: CW-1-130-7-O; 1 uL direct inject  
Misc Info : J&W DB-5MS; 30 m, 0.25 mm, 0.25 um



File : C:\DATA\16\J0003495.D  
Acquired : 28 Oct 2004 11:30 using AcqMethod CW  
Sample Name: CW-CK-1-124-2; 1 uL direct inject A  
Misc Info : J&W DB-5MS; 30 m, 0.25 mm, 0.25 um



El mass spectrum of:

Top: Compound **1** in Organic sample **O/31**, aliquot **CW-1-130-7-O**

Bottom: Authentic reference standard of **Trichloronitromethane** corresponding to compound **1** (MW: **163**)

# GC-CI-MS TECHNIQUE

## METHOD AND ANALYSIS DESCRIPTION

Laboratory code: 31 Sample code(s): O/31, OB/31 Compound number: 1

|                              |              |               |               |
|------------------------------|--------------|---------------|---------------|
| <b>Aliquot codes:</b>        |              |               |               |
| <b>Sample:</b>               | CW-1-130-7-O | <b>Blank:</b> | CW-1-130-8-OB |
| <b>GC-CI-MS Method name:</b> |              | CW-CI         |               |

### METHOD DESCRIPTION

|   |  |   |   |
|---|--|---|---|
| <b>Instrument Manufacturer and Type:</b>    | Agilent 6890/5973 GC/MSD   |   |   |
| <b>Carrier gas:</b>                         | <input checked="" type="checkbox"/> He   | <input type="checkbox"/> N <sub>2</sub>           | <input type="checkbox"/> H <sub>2</sub> <input type="checkbox"/> Other: |
| <b>Flow rate:</b>                           | <input type="checkbox"/> ml/min  | <input checked="" type="checkbox"/> 32 cm/s       |   |
| <b>Flow control:</b>                        | <input type="checkbox"/> Constant Pressure   | <input checked="" type="checkbox"/> Constant Flow |   |
| <b>Injection mode:</b>                      | <input type="checkbox"/> Split → Split ratio =<br><input checked="" type="checkbox"/> Splitless → Splitless time = 0.75 min. |   |   |
| <b>Injector temperature:</b>                | 250 °C   |   |   |
| <b>Column phase:</b>                        | 5% diphenyl 95% dimethyl polysiloxane  |   |   |
| <b>Column Length x ID x Film thickness:</b> | 30 m x 0.25 mm x 0.25 µm   |   |   |
| <b>GC temperature programme:</b>            | 40 °C (3 min), 8 °C/min, 300 °C (3 min)  |   |   |
| <b>Reaction gas:</b>                        | <input checked="" type="checkbox"/> Methane  | <input type="checkbox"/> Isobutane                | <input type="checkbox"/> Ammonia <input type="checkbox"/> Other:        |
| <b>Solvent delay time:</b>                  | 3 min  | <b>Scan range:</b>                                | 50-550 m/z  |
| <b>Electron energy:</b>                     | 235 eV   | <b>Scan time:</b>                                 | 0.35 s  |
| <b>Ionisation polarity:</b>                 | <input type="checkbox"/> Positive<br><input checked="" type="checkbox"/> Negative  | <b>Mass resolution:</b>                           | 0.7 u   |
| <b>Comments:</b>                            |  |   |   |

### ANALYSIS

|  |   |
|--|---|
| <b>Compound identified as:</b>   | <input checked="" type="checkbox"/> Original compound<br><input type="checkbox"/> Methyl ester derivative<br><input type="checkbox"/> TBDMS (t-Butyldimethylsilyl) derivative<br><input type="checkbox"/> TMS (Trimethylsilyl) derivative<br><input type="checkbox"/> Other derivative: |
| <b>Retention parameter used for (peak) identification:</b>                       | <input checked="" type="checkbox"/> Retention time (Rt)<br><input type="checkbox"/> Scan number   |
| <input checked="" type="checkbox"/> Compared to reference chemical:              | Source : <input checked="" type="checkbox"/> Own Synthesis <input type="checkbox"/> Commercial  |
| <input type="checkbox"/> Compared to library spectrum:                           | Source : <input type="checkbox"/> OCAD (code: ) <input type="checkbox"/> NIST<br><input type="checkbox"/> Wiley <input type="checkbox"/> Own <input type="checkbox"/> Other:  |
| <input type="checkbox"/> Not compared to reference chemical or library spectrum: | Intense ions in spectrum are explained; interpretation is supported by the spectral information derived from closely related chemical(s):   |
| <b>Comments:</b>   |   |



File : C:\DATA\16\HP00847.D

Acquired: 17 Nov 2004 17:49

Method : CW-CI

Sample : 1uL of CW-1-130-8-OB

Misc info: HP5-MS, 30m, 0.25 film, NCI-CH4





File : C:\DATA\16\HP00848.D

Acquired: 17 Nov 2004 18:36

Method : CW-CI

Sample : 1uL of CW-1-130-7-O

Misc info: HP5-MS, 30m, 0.25 film, NCI-CH4





File : C:\DATA\16\HP00849.D

Acquired: 17 Nov 2004 19:23

Method : CW-CI

Sample : 1uL of CW-CK-1-124-2-A (compound A)

Misc info: HP5-MS, 30m, 0.25 film, NCI-CH4




CI chromatograms supporting identification of compound **1**; TIC on left; EIC (m/z **117**) on right.

Top: Chromatograms of Organic blank, aliquot **CW-1-130-8-OB** from **OB/31**.


Center: Chromatograms of Organic sample, aliquot **CW-1-130-7-O** from **O/31**, retention time **5.10** min.

Bottom: Chromatograms of authentic reference standard of **Trichloronitromethane**, retention time **5.08** min.

File : C:\DATA\16\HP00848.D  
Acquired : 17 Nov 2004 18:36 using AcqMethod CW-CI  
Sample Name: 1uL of CW-1-130-7-O  
Misc Info : HP5-MS, 30m, 0.25 film, NCI-CH4



File : C:\DATA\16\HP00849.D  
Acquired : 17 Nov 2004 19:23 using AcqMethod CW-CI  
Sample Name: 1uL of CW-CK-1-124-2-A (compound A)  
Misc Info : HP5-MS, 30m, 0.25 film, NCI-CH4



CI mass spectrum of:

Top: Compound **1** in Organic sample **O/31**, aliquot **CW-1-130-7-O**

Bottom: Authentic reference standard of **Trichloronitromethane** corresponding to compound **1** (MW: **163**)

# GC-EI-MS TECHNIQUE

## METHOD AND ANALYSIS DESCRIPTION

Laboratory code: 31 Sample code(s): O/31, OB/31 Compound number: 2

**Aliquot codes:**

**Sample:** CW-1-130-7-O

**Blank:** CW-1-130-8-OB

**GC-EI-MS Method name:** CW

### METHOD DESCRIPTION

|   |  |   |   |
|---|--|---|---|
| <b>Instrument Manufacturer and Type:</b>    | Agilent 6890/5973 GC/MSD   |   |   |
| <b>Carrier gas:</b>                         | <input checked="" type="checkbox"/> He   | <input type="checkbox"/> N <sub>2</sub>           | <input type="checkbox"/> H <sub>2</sub> <input type="checkbox"/> Other: |
| <b>Flow rate:</b>                           | <input type="checkbox"/> ml/min  | <input checked="" type="checkbox"/> 32 cm/s       |   |
| <b>Flow control:</b>                        | <input type="checkbox"/> Constant Pressure   | <input checked="" type="checkbox"/> Constant Flow |   |
| <b>Injection mode:</b>                      | <input type="checkbox"/> Split → Split ratio =<br><input checked="" type="checkbox"/> Splitless → Splitless time = 0.75 min. |   |   |
| <b>Injector temperature:</b>                | 250 °C   |   |   |
| <b>Column phase:</b>                        | 5% diphenyl 95% dimethyl polysiloxane  |   |   |
| <b>Column Length x ID x Film thickness:</b> | 30 m x 0.25 mm x 0.25 µm   |   |   |
| <b>GC temperature programme:</b>            | 40 °C (3 min), 8 °C/min, 300 °C (3 min)  |   |   |
| <b>Solvent delay time:</b>                  | 3 min  | <b>Scan range:</b>                                | 30-600 m/z  |
| <b>Electron energy:</b>                     | 70 eV  | <b>Scan time:</b>                                 | 0.7 s   |
| <b>Ionisation polarity:</b>                 | <input checked="" type="checkbox"/> Positive<br><input type="checkbox"/> Negative  | <b>Mass resolution:</b>                           | 0.7 u   |
| <b>Comments:</b>                            |  |   |   |

### ANALYSIS

|  |   |
|--|---|
| <b>Compound identified as:</b>   | <input checked="" type="checkbox"/> Original compound<br><input type="checkbox"/> Methyl ester derivative<br><input type="checkbox"/> TBDMS (t-Butyldimethylsilyl) derivative<br><input type="checkbox"/> TMS (Trimethylsilyl) derivative<br><input type="checkbox"/> Other derivative: |
| <b>Retention parameter used for (peak) identification:</b>                       | <input checked="" type="checkbox"/> Retention time (Rt)<br><input type="checkbox"/> Scan number   |
| <input checked="" type="checkbox"/> Compared to reference chemical:              | Source : <input checked="" type="checkbox"/> Own Synthesis <input type="checkbox"/> Commercial  |
| <input type="checkbox"/> Compared to library spectrum:                           | Source : <input type="checkbox"/> OCAD (code: ) <input type="checkbox"/> NIST<br><input type="checkbox"/> Wiley <input type="checkbox"/> Own <input type="checkbox"/> Other:  |
| <input type="checkbox"/> Not compared to reference chemical or library spectrum: | Intense ions in spectrum are explained; interpretation is supported by the spectral information derived from closely related chemical(s):   |
| <b>Comments:</b>   |   |



File : C:\DATA\16\J0003493.D

Acquired: 28 Oct 2004 9:55

Method : CW

Sample : CW-1-130-8-OB; 1 uL direct inject

Misc info: J&amp;W DB-5MS; 30 m, 0.25 mm, 0.25 um





File : C:\DATA\16\J0003494.D

Acquired: 28 Oct 2004 10:43

Method : CW

Sample : CW-1-130-7-O; 1 uL direct inject

Misc info: J&amp;W DB-5MS; 30 m, 0.25 mm, 0.25 um





File : C:\DATA\16\J0003496.D

Acquired: 28 Oct 2004 12:17

Method : CW

Sample : CW-CK-1-124-3; 1 uL direct inject B

Misc info: J&amp;W DB-5MS; 30 m, 0.25 mm, 0.25 um



EI chromatograms supporting identification of compound **2**; TIC on left; EIC (m/z **120**) on right.


Top: Chromatograms of Organic blank, aliquot **CW-1-130-8-OB** from **OB/31**.

Center: Chromatograms of Organic sample, aliquot **CW-1-130-7-O** from **O/31**, retention time **10.29** min.


Bottom: Chromatograms of authentic reference standard of **2-(N-Ethyl-N-propylamino)ethylchloride**, retention time **10.29** min.



File : C:\DATA\16\J0003494.D  
Acquired : 28 Oct 2004 10:43 using AcqMethod CW  
Sample Name: CW-1-130-7-O; 1 uL direct inject  
Misc Info : J&W DB-5MS; 30 m, 0.25 mm, 0.25 um



File : C:\DATA\16\J0003496.D  
Acquired : 28 Oct 2004 12:17 using AcqMethod CW  
Sample Name: CW-CK-1-124-3; 1 uL direct inject B  
Misc Info : J&W DB-5MS; 30 m, 0.25 mm, 0.25 um



El mass spectrum of:

Top: Compound **2** in Organic sample **O/31**, aliquot **CW-1-130-7-O**

Bottom: Authentic reference standard of **2-(N-Ethyl-N-propylamino)ethylchloride** corresponding to compound **2** (MW: **149**)

# GC-CI-MS TECHNIQUE

## METHOD AND ANALYSIS DESCRIPTION

Laboratory code: 31 Sample code(s): O/31, OB/31 Compound number: 2

|                              |              |               |               |
|------------------------------|--------------|---------------|---------------|
| <b>Aliquot codes:</b>        |              |               |               |
| <b>Sample:</b>               | CW-1-130-7-O | <b>Blank:</b> | CW-1-130-8-OB |
| <b>GC-CI-MS Method name:</b> |              | CW-CI         |               |

### METHOD DESCRIPTION

|   |  |   |   |
|---|--|---|---|
| <b>Instrument Manufacturer and Type:</b>    | Agilent 6890/5973 GC/MSD   |   |   |
| <b>Carrier gas:</b>                         | <input checked="" type="checkbox"/> He   | <input type="checkbox"/> N <sub>2</sub>           | <input type="checkbox"/> H <sub>2</sub> <input type="checkbox"/> Other:     |
| <b>Flow rate:</b>                           | <input type="checkbox"/> ml/min  | <input checked="" type="checkbox"/> 32 cm/s       |   |
| <b>Flow control:</b>                        | <input type="checkbox"/> Constant Pressure   | <input checked="" type="checkbox"/> Constant Flow |   |
| <b>Injection mode:</b>                      | <input type="checkbox"/> Split → Split ratio =<br><input checked="" type="checkbox"/> Splitless → Splitless time = 0.75 min. |   |   |
| <b>Injector temperature:</b>                | 250 °C   |   |   |
| <b>Column phase:</b>                        | 5% diphenyl 95% dimethyl polysiloxane  |   |   |
| <b>Column Length x ID x Film thickness:</b> | 30 m x 0.25 mm x 0.25 µm   |   |   |
| <b>GC temperature programme:</b>            | 40 °C (3 min), 8 °C/min, 300 °C (3 min)  |   |   |
| <b>Reaction gas:</b>                        | <input type="checkbox"/> Methane   | <input type="checkbox"/> Isobutane                | <input checked="" type="checkbox"/> Ammonia <input type="checkbox"/> Other: |
| <b>Solvent delay time:</b>                  | 3 min  | <b>Scan range:</b>                                | 50-550 m/z  |
| <b>Electron energy:</b>                     | 235 eV   | <b>Scan time:</b>                                 | 0.35 s  |
| <b>Ionisation polarity:</b>                 | <input checked="" type="checkbox"/> Positive<br><input type="checkbox"/> Negative  | <b>Mass resolution:</b>                           | 0.7 u   |
| <b>Comments:</b>                            |  |   |   |

### ANALYSIS

|  |   |
|--|---|
| <b>Compound identified as:</b>   | <input checked="" type="checkbox"/> Original compound<br><input type="checkbox"/> Methyl ester derivative<br><input type="checkbox"/> TBDMS (t-Butyldimethylsilyl) derivative<br><input type="checkbox"/> TMS (Trimethylsilyl) derivative<br><input type="checkbox"/> Other derivative: |
| <b>Retention parameter used for (peak) identification:</b>                       | <input checked="" type="checkbox"/> Retention time (Rt)<br><input type="checkbox"/> Scan number   |
| <input checked="" type="checkbox"/> Compared to reference chemical:              | Source : <input checked="" type="checkbox"/> Own Synthesis <input type="checkbox"/> Commercial  |
| <input type="checkbox"/> Compared to library spectrum:                           | Source : <input type="checkbox"/> OCAD (code: ) <input type="checkbox"/> NIST<br><input type="checkbox"/> Wiley <input type="checkbox"/> Own <input type="checkbox"/> Other:  |
| <input type="checkbox"/> Not compared to reference chemical or library spectrum: | Intense ions in spectrum are explained; interpretation is supported by the spectral information derived from closely related chemical(s):   |
| <b>Comments:</b>   |   |



File : C:\DATA\16\HP00826.D

Acquired: 29 Oct 2004 14:59

Method : CW-CI

Sample : 1uL of CW-1-130-8-OB

Misc info: HP5-MS, 30m, 0.25 film, PCI-NH4





File : C:\DATA\16\HP00827.D

Acquired: 29 Oct 2004 15:47

Method : CW-CI

Sample : 1uL of CW-1-130-7-O

Misc info: HP5-MS, 30m, 0.25 film, PCI-NH4m





File : C:\DATA\16\HP00829.D

Acquired: 29 Oct 2004 17:21

Method : CW-CI

Sample : 1uL of CW-CK-1-124-3-B (compound B)

Misc info: HP5-MS, 30m, 0.25 film, PCI-NH4




CI chromatograms supporting identification of compound **2**; TIC on left; EIC (m/z **150**) on right.

Top: Chromatograms of Organic blank, aliquot **CW-1-130-8-OB** from **OB/31**.


Center: Chromatograms of Organic sample, aliquot **CW-1-130-7-O** from **O/31**, retention time **9.85** min.

Bottom: Chromatograms of authentic reference standard of **2-(N-Ethyl-N-propylamino)ethylchloride**, retention time **9.85** min.

File : C:\DATA\16\HP00827.D  
Acquired : 29 Oct 2004 15:47 using AcqMethod CW-CI  
Sample Name: 1uL of CW-1-130-7-O  
Misc Info : HP5-MS, 30m, 0.25 film, PCI-NH4m



File : C:\DATA\16\HP00829.D  
Acquired : 29 Oct 2004 17:21 using AcqMethod CW-CI  
Sample Name: 1uL of CW-CK-1-124-3-B (compound B)  
Misc Info : HP5-MS, 30m, 0.25 film, PCI-NH4



CI mass spectrum of:

Top: Compound **2** in Organic sample **O/31**, aliquot **CW-1-130-7-O**

Bottom: Authentic reference standard of **2-(N-Ethyl-N-propylamino)ethylchloride** corresponding to compound **2** (MW: **149**)

# GC-EI-MS TECHNIQUE

## METHOD AND ANALYSIS DESCRIPTION

Laboratory code: 31 Sample code(s): O/31, OB/31 Compound number: 3

**Aliquot codes:**

**Sample:** CW-1-130-7-O

**Blank:** CW-1-130-8-OB

**GC-EI-MS Method name:** CW

### METHOD DESCRIPTION

|   |  |   |            |
|---|--|---|------------|
| <b>Instrument Manufacturer and Type:</b>    | Agilent 6890/5973 GC/MSD   |   |            |
| <b>Carrier gas:</b>                         | <input checked="" type="checkbox"/> He <input type="checkbox"/> N <sub>2</sub> <input type="checkbox"/> H <sub>2</sub> <input type="checkbox"/> Other: |   |            |
| <b>Flow rate:</b>                           | <input type="checkbox"/> ml/min  | <input checked="" type="checkbox"/> 32 cm/s |            |
| <b>Flow control:</b>                        | <input type="checkbox"/> Constant Pressure <input checked="" type="checkbox"/> Constant Flow   |   |            |
| <b>Injection mode:</b>                      | <input type="checkbox"/> Split → Split ratio =<br><input checked="" type="checkbox"/> Splitless → Splitless time = 0.75 min.                           |   |            |
| <b>Injector temperature:</b>                | 250 °C   |   |            |
| <b>Column phase:</b>                        | 5% diphenyl 95% dimethyl polysiloxane  |   |            |
| <b>Column Length x ID x Film thickness:</b> | 30 m x 0.25 mm x 0.25 µm   |   |            |
| <b>GC temperature programme:</b>            | 40 °C (3 min), 8 °C/min, 300 °C (3 min)  |   |            |
| <b>Solvent delay time:</b>                  | 3 min  | <b>Scan range:</b>                          | 30-600 m/z |
| <b>Electron energy:</b>                     | 70 eV  | <b>Scan time:</b>                           | 0.7 s      |
| <b>Ionisation polarity:</b>                 | <input checked="" type="checkbox"/> Positive<br><input type="checkbox"/> Negative  | <b>Mass resolution:</b>                     | 0.7 u      |
| <b>Comments:</b>                            |  |   |            |

### ANALYSIS

|  |   |
|--|---|
| <b>Compound identified as:</b>   | <input checked="" type="checkbox"/> Original compound<br><input type="checkbox"/> Methyl ester derivative<br><input type="checkbox"/> TBDMS (t-Butyldimethylsilyl) derivative<br><input type="checkbox"/> TMS (Trimethylsilyl) derivative<br><input type="checkbox"/> Other derivative: |
| <b>Retention parameter used for (peak) identification:</b>                       | <input checked="" type="checkbox"/> Retention time (Rt)<br><input type="checkbox"/> Scan number   |
| <input checked="" type="checkbox"/> Compared to reference chemical:              | Source : <input type="checkbox"/> Own Synthesis <input checked="" type="checkbox"/> Commercial  |
| <input type="checkbox"/> Compared to library spectrum:                           | Source : <input type="checkbox"/> OCAD (code: ) <input type="checkbox"/> NIST<br><input type="checkbox"/> Wiley <input type="checkbox"/> Own <input type="checkbox"/> Other:  |
| <input type="checkbox"/> Not compared to reference chemical or library spectrum: | Intense ions in spectrum are explained; interpretation is supported by the spectral information derived from closely related chemical(s):   |
| <b>Comments:</b>   |   |



File : C:\DATA\16\J0003493.D

Acquired: 28 Oct 2004 9:55

Method : CW

Sample : CW-1-130-8-OB; 1 uL direct inject

Misc info: J&W DB-5MS; 30 m, 0.25 mm, 0.25 um





File : C:\DATA\16\J0003494.D

Acquired: 28 Oct 2004 10:43

Method : CW

Sample : CW-1-130-7-O; 1 uL direct inject

Misc info: J&W DB-5MS; 30 m, 0.25 mm, 0.25 um





File : C:\DATA\16\J0003499.D

Acquired: 28 Oct 2004 14:50

Method : CW

Sample : CW-CK-1-124-6; 1 uL direct inject C

Misc info: J&W DB-5MS; 30 m, 0.25 mm, 0.25 um




EI chromatograms supporting identification of compound **3**; TIC on left; EIC (m/z **106**) on right.

Top: Chromatograms of Organic blank, aliquot **CW-1-130-8-OB** from **OB/31**.


Center: Chromatograms of Organic sample, aliquot **CW-1-130-7-O** from **O/31**, retention time **11.45** min.

Bottom: Chromatograms of authentic reference standard of **2-(N,N-Diisopropylamino)ethylchloride**, retention time **11.45** min.

File : C:\DATA\16\J0003494.D  
Acquired : 28 Oct 2004 10:43 using AcqMethod CW  
Sample Name: CW-1-130-7-O; 1 uL direct inject  
Misc Info : J&W DB-5MS; 30 m, 0.25 mm, 0.25 um



File : C:\DATA\16\J0003499.D  
Acquired : 28 Oct 2004 14:50 using AcqMethod CW  
Sample Name: CW-CK-1-124-6; 1 uL direct inject C  
Misc Info : J&W DB-5MS; 30 m, 0.25 mm, 0.25 um



El mass spectrum of:

Top: Compound **3** in Organic sample **O/31**, aliquot **CW-1-130-7-O**

Bottom: Authentic reference standard of **2-(N,N-Diisopropylamino)ethylchloride** corresponding to compound **3** (MW: **163**)

# GC-CI-MS TECHNIQUE

## METHOD AND ANALYSIS DESCRIPTION

Laboratory code: 31 Sample code(s): O/31, OB/31 Compound number: 3

|                              |              |               |               |
|------------------------------|--------------|---------------|---------------|
| <b>Aliquot codes:</b>        |              |               |               |
| <b>Sample:</b>               | CW-1-130-7-O | <b>Blank:</b> | CW-1-130-8-OB |
| <b>GC-CI-MS Method name:</b> |              | CW-CI         |               |

### METHOD DESCRIPTION

|   |  |                         |            |
|---|--|-------------------------|------------|
| <b>Instrument Manufacturer and Type:</b>    | Agilent 6890/5973 GC/MSD   |                         |            |
| <b>Carrier gas:</b>                         | <input checked="" type="checkbox"/> He <input type="checkbox"/> N <sub>2</sub> <input type="checkbox"/> H <sub>2</sub> <input type="checkbox"/> Other: |                         |            |
| <b>Flow rate:</b>                           | <input type="checkbox"/> ml/min <input checked="" type="checkbox"/> 32 cm/s  |                         |            |
| <b>Flow control:</b>                        | <input type="checkbox"/> Constant Pressure <input checked="" type="checkbox"/> Constant Flow   |                         |            |
| <b>Injection mode:</b>                      | <input type="checkbox"/> Split → Split ratio =<br><input checked="" type="checkbox"/> Splitless → Splitless time = 0.75 min.                           |                         |            |
| <b>Injector temperature:</b>                | 250 °C   |                         |            |
| <b>Column phase:</b>                        | 5% diphenyl 95% dimethyl polysiloxane  |                         |            |
| <b>Column Length x ID x Film thickness:</b> | 30 m x 0.25 mm x 0.25 µm   |                         |            |
| <b>GC temperature programme:</b>            | 40 °C (3 min), 8 °C/min, 300 °C (3 min)  |                         |            |
| <b>Reaction gas:</b>                        | <input type="checkbox"/> Methane <input type="checkbox"/> Isobutane <input checked="" type="checkbox"/> Ammonia <input type="checkbox"/> Other:        |                         |            |
| <b>Solvent delay time:</b>                  | 3 min  | <b>Scan range:</b>      | 50-550 m/z |
| <b>Electron energy:</b>                     | 235 eV   | <b>Scan time:</b>       | 0.35 s     |
| <b>Ionisation polarity:</b>                 | <input checked="" type="checkbox"/> Positive<br><input type="checkbox"/> Negative  | <b>Mass resolution:</b> | 0.7 u      |
| <b>Comments:</b>                            |  |                         |            |

### ANALYSIS

|  |   |
|--|---|
| <b>Compound identified as:</b>   | <input checked="" type="checkbox"/> Original compound<br><input type="checkbox"/> Methyl ester derivative<br><input type="checkbox"/> TBDMS (t-Butyldimethylsilyl) derivative<br><input type="checkbox"/> TMS (Trimethylsilyl) derivative<br><input type="checkbox"/> Other derivative: |
| <b>Retention parameter used for (peak) identification:</b>                       | <input checked="" type="checkbox"/> Retention time (Rt)<br><input type="checkbox"/> Scan number   |
| <input checked="" type="checkbox"/> Compared to reference chemical:              | Source : <input type="checkbox"/> Own Synthesis <input checked="" type="checkbox"/> Commercial  |
| <input type="checkbox"/> Compared to library spectrum:                           | Source : <input type="checkbox"/> OCAD (code: ) <input type="checkbox"/> NIST<br><input type="checkbox"/> Wiley <input type="checkbox"/> Own <input type="checkbox"/> Other:  |
| <input type="checkbox"/> Not compared to reference chemical or library spectrum: | Intense ions in spectrum are explained; interpretation is supported by the spectral information derived from closely related chemical(s):   |
| <b>Comments:</b>   |   |





File : C:\DATA\16\HP00826.D

Acquired: 29 Oct 2004 14:59

Method : CW-CI

Sample : 1uL of CW-1-130-8-OB

Misc info: HP5-MS, 30m, 0.25 film, PCI-NH4





File : C:\DATA\16\HP00827.D

Acquired: 29 Oct 2004 15:47

Method : CW-CI

Sample : 1uL of CW-1-130-7-O

Misc info: HP5-MS, 30m, 0.25 film, PCI-NH4m





File : C:\DATA\16\HP00830.D

Acquired: 29 Oct 2004 18:09

Method : CW-CI

Sample : 1uL of CW-CK-1-124-6-C (compound C)

Misc info: HP5-MS, 30m, 0.25 film, PCI-NH4




CI chromatograms supporting identification of compound **3**; TIC on left; EIC (m/z **164**) on right.

Top: Chromatograms of Organic blank, aliquot **CW-1-130-8-OB** from **OB/31**.


Center: Chromatograms of Organic sample, aliquot **CW-1-130-7-O** from **O/31**, retention time **11.01** min.

Bottom: Chromatograms of authentic reference standard of **2-(N,N-Diisopropylamino)ethylchloride**, retention time **11.01** min.

File : C:\DATA\16\HP00827.D  
Acquired : 29 Oct 2004 15:47 using AcqMethod CW-CI  
Sample Name: 1uL of CW-1-130-7-O  
Misc Info : HP5-MS, 30m, 0.25 film, PCI-NH4m



File : C:\DATA\16\HP00830.D  
Acquired : 29 Oct 2004 18:09 using AcqMethod CW-CI  
Sample Name: 1uL of CW-CK-1-124-6-C (compound C)  
Misc Info : HP5-MS, 30m, 0.25 film, PCI-NH4



CI mass spectrum of:

Top: Compound **3** in Organic sample **O/31**, aliquot **CW-1-130-7-O**

Bottom: Authentic reference standard of **2-(N,N-Diisopropylamino)ethylchloride** corresponding to compound **3** (MW: **163**)

# GC-EI-MS TECHNIQUE

## METHOD AND ANALYSIS DESCRIPTION

Laboratory code: 31 Sample code(s): O/31, OB/31 Compound number: 4

**Aliquot codes:**

**Sample:** CW-1-130-7-O

**Blank:** CW-1-130-8-OB

**GC-EI-MS Method name:** CW

### METHOD DESCRIPTION

|   |  |   |   |
|---|--|---|---|
| <b>Instrument Manufacturer and Type:</b>    | Agilent 6890/5973 GC/MSD   |   |   |
| <b>Carrier gas:</b>                         | <input checked="" type="checkbox"/> He   | <input type="checkbox"/> N <sub>2</sub>           | <input type="checkbox"/> H <sub>2</sub> <input type="checkbox"/> Other: |
| <b>Flow rate:</b>                           | <input type="checkbox"/> ml/min  | <input checked="" type="checkbox"/> 32 cm/s       |   |
| <b>Flow control:</b>                        | <input type="checkbox"/> Constant Pressure   | <input checked="" type="checkbox"/> Constant Flow |   |
| <b>Injection mode:</b>                      | <input type="checkbox"/> Split → Split ratio =<br><input checked="" type="checkbox"/> Splitless → Splitless time = 0.75 min. |   |   |
| <b>Injector temperature:</b>                | 250 °C   |   |   |
| <b>Column phase:</b>                        | 5% diphenyl 95% dimethyl polysiloxane  |   |   |
| <b>Column Length x ID x Film thickness:</b> | 30 m x 0.25 mm x 0.25 µm   |   |   |
| <b>GC temperature programme:</b>            | 40 °C (3 min), 8 °C/min, 300 °C (3 min)  |   |   |
| <b>Solvent delay time:</b>                  | 3 min  | <b>Scan range:</b>                                | 30-600 m/z  |
| <b>Electron energy:</b>                     | 70 eV  | <b>Scan time:</b>                                 | 0.7 s   |
| <b>Ionisation polarity:</b>                 | <input checked="" type="checkbox"/> Positive<br><input type="checkbox"/> Negative  | <b>Mass resolution:</b>                           | 0.7 u   |
| <b>Comments:</b>                            |  |   |   |

### ANALYSIS

|  |   |
|--|---|
| <b>Compound identified as:</b>   | <input checked="" type="checkbox"/> Original compound<br><input type="checkbox"/> Methyl ester derivative<br><input type="checkbox"/> TBDMS (t-Butyldimethylsilyl) derivative<br><input type="checkbox"/> TMS (Trimethylsilyl) derivative<br><input type="checkbox"/> Other derivative: |
| <b>Retention parameter used for (peak) identification:</b>                       | <input checked="" type="checkbox"/> Retention time (Rt)<br><input type="checkbox"/> Scan number   |
| <input checked="" type="checkbox"/> Compared to reference chemical:              | Source : <input type="checkbox"/> Own Synthesis <input checked="" type="checkbox"/> Commercial  |
| <input type="checkbox"/> Compared to library spectrum:                           | Source : <input type="checkbox"/> OCAD (code: ) <input type="checkbox"/> NIST<br><input type="checkbox"/> Wiley <input type="checkbox"/> Own <input type="checkbox"/> Other:  |
| <input type="checkbox"/> Not compared to reference chemical or library spectrum: | Intense ions in spectrum are explained; interpretation is supported by the spectral information derived from closely related chemical(s):   |
| <b>Comments:</b>   |   |



File : C:\DATA\16\J0003493.D

Acquired: 28 Oct 2004 9:55

Method : CW

Sample : CW-1-130-8-OB; 1 uL direct inject

Misc info: J&W DB-5MS; 30 m, 0.25 mm, 0.25 um





File : C:\DATA\16\J0003494.D

Acquired: 28 Oct 2004 10:43

Method : CW

Sample : CW-1-130-7-O; 1 uL direct inject

Misc info: J&W DB-5MS; 30 m, 0.25 mm, 0.25 um





File : C:\DATA\16\J0003498.D

Acquired: 28 Oct 2004 13:52

Method : CW

Sample : CW-CK-1-124-5; 1 uL direct inject D

Misc info: J&W DB-5MS; 30 m, 0.25 mm, 0.25 um




EI chromatograms supporting identification of compound **4**; TIC on left; EIC (m/z **86**) on right.

Top: Chromatograms of Organic blank, aliquot **CW-1-130-8-OB** from **OB/31**.


Center: Chromatograms of Organic sample, aliquot **CW-1-130-7-O** from **O/31**, retention time **8.43** min.

Bottom: Chromatograms of authentic reference standard of **2-(N,N-Diethylamino)ethylchloride**, retention time **8.44** min.

File : C:\DATA\16\J0003494.D  
Acquired : 28 Oct 2004 10:43 using AcqMethod CW  
Sample Name: CW-1-130-7-O; 1 uL direct inject  
Misc Info : J&W DB-5MS; 30 m, 0.25 mm, 0.25 um



File : C:\DATA\16\J0003498.D  
Acquired : 28 Oct 2004 13:52 using AcqMethod CW  
Sample Name: CW-CK-1-124-5; 1 uL direct inject D  
Misc Info : J&W DB-5MS; 30 m, 0.25 mm, 0.25 um



El mass spectrum of:

Top: Compound **4** in Organic sample **O/31**, aliquot **CW-1-130-7-O**

Bottom: Authentic reference standard of **2-(N,N-Diethylamino)ethylchloride** corresponding to compound **4** (MW: **135**)

# GC-CI-MS TECHNIQUE

## METHOD AND ANALYSIS DESCRIPTION

Laboratory code: 31 Sample code(s): O/31, OB/31 Compound number: 4

|                              |              |               |               |
|------------------------------|--------------|---------------|---------------|
| <b>Aliquot codes:</b>        |              |               |               |
| <b>Sample:</b>               | CW-1-130-7-O | <b>Blank:</b> | CW-1-130-8-OB |
| <b>GC-CI-MS Method name:</b> |              | CW-CI         |               |

### METHOD DESCRIPTION

|   |  |   |   |
|---|--|---|---|
| <b>Instrument Manufacturer and Type:</b>    | Agilent 6890/5973 GC/MSD   |   |   |
| <b>Carrier gas:</b>                         | <input checked="" type="checkbox"/> He   | <input type="checkbox"/> N <sub>2</sub>           | <input type="checkbox"/> H <sub>2</sub> <input type="checkbox"/> Other:     |
| <b>Flow rate:</b>                           | <input type="checkbox"/> ml/min  | <input checked="" type="checkbox"/> 32 cm/s       |   |
| <b>Flow control:</b>                        | <input type="checkbox"/> Constant Pressure   | <input checked="" type="checkbox"/> Constant Flow |   |
| <b>Injection mode:</b>                      | <input type="checkbox"/> Split → Split ratio =<br><input checked="" type="checkbox"/> Splitless → Splitless time = 0.75 min. |   |   |
| <b>Injector temperature:</b>                | 250 °C   |   |   |
| <b>Column phase:</b>                        | 5% diphenyl 95% dimethyl polysiloxane  |   |   |
| <b>Column Length x ID x Film thickness:</b> | 30 m x 0.25 mm x 0.25 µm   |   |   |
| <b>GC temperature programme:</b>            | 40 °C (3 min), 8 °C/min, 300 °C (3 min)  |   |   |
| <b>Reaction gas:</b>                        | <input type="checkbox"/> Methane   | <input type="checkbox"/> Isobutane                | <input checked="" type="checkbox"/> Ammonia <input type="checkbox"/> Other: |
| <b>Solvent delay time:</b>                  | 3 min  | <b>Scan range:</b>                                | 50-550 m/z  |
| <b>Electron energy:</b>                     | 235 eV   | <b>Scan time:</b>                                 | 0.35 s  |
| <b>Ionisation polarity:</b>                 | <input checked="" type="checkbox"/> Positive<br><input type="checkbox"/> Negative  | <b>Mass resolution:</b>                           | 0.7 u   |
| <b>Comments:</b>                            |  |   |   |

### ANALYSIS

|  |   |
|--|---|
| <b>Compound identified as:</b>   | <input checked="" type="checkbox"/> Original compound<br><input type="checkbox"/> Methyl ester derivative<br><input type="checkbox"/> TBDMS (t-Butyldimethylsilyl) derivative<br><input type="checkbox"/> TMS (Trimethylsilyl) derivative<br><input type="checkbox"/> Other derivative: |
| <b>Retention parameter used for (peak) identification:</b>                       | <input checked="" type="checkbox"/> Retention time (Rt)<br><input type="checkbox"/> Scan number   |
| <input checked="" type="checkbox"/> Compared to reference chemical:              | Source : <input type="checkbox"/> Own Synthesis <input checked="" type="checkbox"/> Commercial  |
| <input type="checkbox"/> Compared to library spectrum:                           | Source : <input type="checkbox"/> OCAD (code: ) <input type="checkbox"/> NIST<br><input type="checkbox"/> Wiley <input type="checkbox"/> Own <input type="checkbox"/> Other:  |
| <input type="checkbox"/> Not compared to reference chemical or library spectrum: | Intense ions in spectrum are explained; interpretation is supported by the spectral information derived from closely related chemical(s):   |
| <b>Comments:</b>   |   |



File : C:\DATA\16\HP00826.D

Acquired: 29 Oct 2004 14:59

Method : CW-CI

Sample : 1uL of CW-1-130-8-OB

Misc info: HP5-MS, 30m, 0.25 film, PCI-NH4





File : C:\DATA\16\HP00827.D

Acquired: 29 Oct 2004 15:47

Method : CW-CI

Sample : 1uL of CW-1-130-7-O

Misc info: HP5-MS, 30m, 0.25 film, PCI-NH4m





File : C:\DATA\16\HP00831.D

Acquired: 29 Oct 2004 18:56

Method : CW-CI

Sample : 1uL of CW-CK-1-124-5-D (compound D)

Misc info: HP5-MS, 30m, 0.25 film, PCI-NH4




CI chromatograms supporting identification of compound 4; TIC on left; EIC (m/z 136) on right.

Top: Chromatograms of Organic blank, aliquot CW-1-130-8-OB from OB/31.


Center: Chromatograms of Organic sample, aliquot CW-1-130-7-O from O/31, retention time 8.02 min.

Bottom: Chromatograms of authentic reference standard of 2-(N,N-Diethylamino)ethylchloride, retention time 8.02 min.

File : C:\DATA\16\HP00827.D  
Acquired : 29 Oct 2004 15:47 using AcqMethod CW-CI  
Sample Name: 1uL of CW-1-130-7-O  
Misc Info : HP5-MS, 30m, 0.25 film, PCI-NH4m



File : C:\DATA\16\HP00831.D  
Acquired : 29 Oct 2004 18:56 using AcqMethod CW-CI  
Sample Name: 1uL of CW-CK-1-124-5-D (compound D)  
Misc Info : HP5-MS, 30m, 0.25 film, PCI-NH4



CI mass spectrum of:

Top: Compound **4** in Organic sample **O/31**, aliquot **CW-1-130-7-O**

Bottom: Authentic reference standard of **2-(N,N-Diethylamino)ethylchloride** corresponding to compound **4** (MW: **135**)



## SAMPLE PREPARATION DESCRIPTION


Laboratory code: 31 Sample code(s): L/31 Sample Blank code: LB/31

### 1. Sample preparation

| Sample/<br>Aliquot Code | Specification of Sample/<br>Type of Sample Preparation | Amount/<br>Volume | Sample Preparation Procedures  | End<br>Volume | Resulting<br>Aliquot Code     | Analytical<br>Technique(s) |
|-------------------------|--|-------------------|--|---------------|-------------------------------|----------------------------|
| L/31<br>LB/31           | TMS derivative of sample                               | 500 µL            | 500 µL evaporated to dryness (N <sub>2</sub> , no heat). Added 500 µL BSTFA. Reacted at 60°C for 30 minutes. | 500 µL        | CW-1-131-4-L<br>CW-1-131-3-LB | GC/MS-EI<br>GC/MS-CI       |
| CW-1-131-4-L            | Dilution   | 20 µL             | Diluted TMS derivitized samples with 180 µL BSTFA.   | 200 µL        | CW-1-131-5-L                  | GC/MS-EI<br>GC/MS-CI       |
|                         |  |                   |  |               |                               |                            |
|                         |  |                   |  |               |                               |                            |
|                         |  |                   |  |               |                               |                            |
|                         |  |                   |  |               |                               |                            |
|                         |  |                   |  |               |                               |                            |

### 2. Additional information

Note: sample CW-1-131-4-L was diluted with the same batch of BSTFA with which it had been derivitized to produce sample CW-1-131-5-L. Blank CW-1-131-3-LB was used directly as a blank (worst case – any scheduled compounds would be seen 10x as compared to the diluted sample).



Note: This flowchart is for visualization only; see the preceding sample preparation description page for sample aliquot numbers

# GC-EI-MS TECHNIQUE

## METHOD AND ANALYSIS DESCRIPTION

Laboratory code: 31 Sample code(s): L/31, LB/31 Compound number: 5

Aliquot codes:

Sample: CW-1-131-5-L Blank: CW-1-131-3-LB

GC-EI-MS Method name: TMS\_A



### METHOD DESCRIPTION

|   |  |                         |            |
|---|--|-------------------------|------------|
| <b>Instrument Manufacturer and Type:</b>    | Agilent 6890/5973 GC/MSD   |                         |            |
| <b>Carrier gas:</b>                         | <input checked="" type="checkbox"/> He <input type="checkbox"/> N <sub>2</sub> <input type="checkbox"/> H <sub>2</sub> <input type="checkbox"/> Other: |                         |            |
| <b>Flow rate:</b>                           | <input type="checkbox"/> ml/min <input checked="" type="checkbox"/> 38 cm/s  |                         |            |
| <b>Flow control:</b>                        | <input type="checkbox"/> Constant Pressure <input checked="" type="checkbox"/> Constant Flow   |                         |            |
| <b>Injection mode:</b>                      | <input type="checkbox"/> Split → Split ratio =<br><input checked="" type="checkbox"/> Splitless → Splitless time = 0.70 min.                           |                         |            |
| <b>Injector temperature:</b>                | 250 °C   |                         |            |
| <b>Column phase:</b>                        | 5% diphenyl 95% dimethyl polysiloxane  |                         |            |
| <b>Column Length x ID x Film thickness:</b> | 30 m x 0.25 mm x 0.25 µm   |                         |            |
| <b>GC temperature programme:</b>            | 70 °C (8 min), 8 °C/min, 300 °C (3 min)  |                         |            |
| <b>Solvent delay time:</b>                  | 8 min  | <b>Scan range:</b>      | 40-600 m/z |
| <b>Electron energy:</b>                     | 70 eV  | <b>Scan time:</b>       | 0.7 s      |
| <b>Ionisation polarity:</b>                 | <input checked="" type="checkbox"/> Positive<br><input type="checkbox"/> Negative  | <b>Mass resolution:</b> | 0.7 u      |
| <b>Comments:</b>                            |  |                         |            |



### ANALYSIS

|  |   |
|--|---|
| <b>Compound identified as:</b>   | <input type="checkbox"/> Original compound<br><input type="checkbox"/> Methyl ester derivative<br><input type="checkbox"/> TBDMS (t-Butyldimethylsilyl) derivative<br><input checked="" type="checkbox"/> TMS (Trimethylsilyl) derivative<br><input type="checkbox"/> Other derivative: |
| <b>Retention parameter used for (peak) identification:</b>                       | <input checked="" type="checkbox"/> Retention time (Rt)<br><input type="checkbox"/> Scan number   |
| <input checked="" type="checkbox"/> Compared to reference chemical:              | Source : <input type="checkbox"/> Own Synthesis <input checked="" type="checkbox"/> Commercial  |
| <input type="checkbox"/> Compared to library spectrum:                           | Source : <input type="checkbox"/> OCAD (code: ) <input type="checkbox"/> NIST<br><input type="checkbox"/> Wiley <input type="checkbox"/> Own <input type="checkbox"/> Other:  |
| <input type="checkbox"/> Not compared to reference chemical or library spectrum: | Intense ions in spectrum are explained; interpretation is supported by the spectral information derived from closely related chemical(s):   |
| <b>Comments:</b>   |   |



File : C:\DATA\16\J0003609.D  
Acquired: 29 Nov 2004 16:01 Sample : CW-1-131-3-LB  
Method : TMS\_A Misc info:



File : C:\DATA\16\J0003610.D  
Acquired: 29 Nov 2004 16:46 Sample : CW-1-131-5-L  
Method : TMS\_A Misc info:



File : C:\DATA\16\J0003613.D  
Acquired: 29 Nov 2004 19:02 Sample : CW-CK-1-127-6  
Method : TMS\_A Misc info:




EI chromatograms supporting identification of compound **5**; TIC on left; EIC (m/z **262**) on right.

Top: Chromatograms of Liquid blank, aliquot **CW-1-131-3-LB** from **LB/31**.


Center: Chromatograms of Liquid sample, aliquot **CW-1-131-5-L** from **L/31**, retention time **22.07** min.

Bottom: Chromatograms of TMS derivative of the authentic reference standard of **Triethanolamine** [tris(2-trimethylsiloxyethyl)amine], retention time **22.07** min.

File : C:\DATA\16\J0003610.D  
Acquired : 29 Nov 2004 16:46 using AcqMethod TMS\_A  
Sample Name: CW-1-131-5-L  
Misc Info :



File : C:\DATA\16\J0003613.D  
Acquired : 29 Nov 2004 19:02 using AcqMethod TMS\_A  
Sample Name: CW-CK-1-127-6  
Misc Info :



El mass spectrum of:

Top: Compound **5** in Liquid sample **L/31**, aliquot **CW-1-131-5-L**

Bottom: TMS derivative of the authentic reference standard of **Triethanolamine**  
[tris(2-trimethylsiloxyethyl)amine] corresponding to compound **5** (MW: **365**)

[tris(2-

# GC-CI-MS TECHNIQUE

## METHOD AND ANALYSIS DESCRIPTION

Laboratory code: 31 Sample code(s): L/31, LB/31 Compound number: 5

|                              |              |               |               |
|------------------------------|--------------|---------------|---------------|
| <b>Aliquot codes:</b>        |              |               |               |
| <b>Sample:</b>               | CW-1-131-5-L | <b>Blank:</b> | CW-1-131-3-LB |
| <b>GC-CI-MS Method name:</b> |              | CW-CI-TM      |               |

### METHOD DESCRIPTION

|   |  |   |   |
|---|--|---|---|
| <b>Instrument Manufacturer and Type:</b>    | Agilent 6890/5973 GC/MSD   |   |   |
| <b>Carrier gas:</b>                         | <input checked="" type="checkbox"/> He   | <input type="checkbox"/> N <sub>2</sub>           | <input type="checkbox"/> H <sub>2</sub> <input type="checkbox"/> Other:     |
| <b>Flow rate:</b>                           | <input type="checkbox"/> ml/min  | <input checked="" type="checkbox"/> 32 cm/s       |   |
| <b>Flow control:</b>                        | <input type="checkbox"/> Constant Pressure   | <input checked="" type="checkbox"/> Constant Flow |   |
| <b>Injection mode:</b>                      | <input type="checkbox"/> Split → Split ratio =<br><input checked="" type="checkbox"/> Splitless → Splitless time = 0.75 min. |   |   |
| <b>Injector temperature:</b>                | 250 °C   |   |   |
| <b>Column phase:</b>                        | 5% diphenyl 95% dimethyl polysiloxane  |   |   |
| <b>Column Length x ID x Film thickness:</b> | 30 m x 0.25 mm x 0.25 µm   |   |   |
| <b>GC temperature programme:</b>            | 70 °C (8 min), 8 °C/min, 300 °C (3 min)  |   |   |
| <b>Reaction gas:</b>                        | <input type="checkbox"/> Methane   | <input type="checkbox"/> Isobutane                | <input checked="" type="checkbox"/> Ammonia <input type="checkbox"/> Other: |
| <b>Solvent delay time:</b>                  | 8 min  | <b>Scan range:</b>                                | 50-550 m/z  |
| <b>Electron energy:</b>                     | 235 eV   | <b>Scan time:</b>                                 | 0.35 s  |
| <b>Ionisation polarity:</b>                 | <input checked="" type="checkbox"/> Positive<br><input type="checkbox"/> Negative  | <b>Mass resolution:</b>                           | 0.7 u   |
| <b>Comments:</b>                            |  |   |   |

### ANALYSIS

|  |   |
|--|---|
| <b>Compound identified as:</b>   | <input type="checkbox"/> Original compound<br><input type="checkbox"/> Methyl ester derivative<br><input type="checkbox"/> TBDMS (t-Butyldimethylsilyl) derivative<br><input checked="" type="checkbox"/> TMS (Trimethylsilyl) derivative<br><input type="checkbox"/> Other derivative: |
| <b>Retention parameter used for (peak) identification:</b>                       | <input checked="" type="checkbox"/> Retention time (Rt)<br><input type="checkbox"/> Scan number   |
| <input checked="" type="checkbox"/> Compared to reference chemical:              | Source : <input type="checkbox"/> Own Synthesis <input checked="" type="checkbox"/> Commercial  |
| <input type="checkbox"/> Compared to library spectrum:                           | Source : <input type="checkbox"/> OCAD (code: ) <input type="checkbox"/> NIST<br><input type="checkbox"/> Wiley <input type="checkbox"/> Own <input type="checkbox"/> Other:  |
| <input type="checkbox"/> Not compared to reference chemical or library spectrum: | Intense ions in spectrum are explained; interpretation is supported by the spectral information derived from closely related chemical(s):   |
| <b>Comments:</b>   |   |



File : C:\DATA\16\HP00857.D

Acquired: 23 Nov 2004 16:50

Method : CW-CI-TM

Sample : 1uL of CW-1-131-3-LB

Misc info: HP5-MS, 30m, 0.25 film, PCI-NH3





File : C:\DATA\16\HP00864.D

Acquired: 26 Nov 2004 11:48

Method : CW-CI-TM

Sample : 1uL of CW-1-131-5-L (dilute)

Misc info: HP5-MS, 30m, 0.25 film, PCI-NH3





File : C:\DATA\16\HP00861.D

Acquired: 24 Nov 2004 18:09

Method : CW-CI-TM

Sample : 1uL of CW-CK-1-127-6 (compound J)

Misc info: HP5-MS, 30m, 0.25 film, PCI-NH3




CI chromatograms supporting identification of compound **5**; TIC on left; EIC (m/z **366**) on right.

Top: Chromatograms of Liquid blank, aliquot **CW-1-131-3-LB** from **LB/31**.


Center: Chromatograms of Liquid sample, aliquot **CW-1-131-5-L** from **L/31**, retention time **22.30** min.

Bottom: Chromatograms of TMS derivative of the authentic reference standard of **Triethanolamine** [tris(2-trimethylsiloxyethyl)amine], retention time **22.29** min.

File : C:\DATA\16\HP00864.D  
Acquired : 26 Nov 2004 11:48 using AcqMethod CW-CI-TM  
Sample Name: 1uL of CW-1-131-5-L (dilute)  
Misc Info : HP5-MS, 30m, 0.25 film, PCI-NH3



File : C:\DATA\16\HP00861.D  
Acquired : 24 Nov 2004 18:09 using AcqMethod CW-CI-TM  
Sample Name: 1uL of CW-CK-1-127-6 (compound J)  
Misc Info : HP5-MS, 30m, 0.25 film, PCI-NH3



CI mass spectrum of:

Top: Compound **5** in Liquid sample **L/31**, aliquot **CW-1-131-5-L**

Bottom: TMS derivative of the authentic reference standard of **Triethanolamine**  
[tris(2-trimethylsiloxyethyl)amine] corresponding to compound **5** (MW: **365**)

[tris(2-



# GC-EI-MS TECHNIQUE

## METHOD AND ANALYSIS DESCRIPTION

Laboratory code: 31 Sample code(s): L/31, LB/31 Compound number: 6

**Aliquot codes:**

**Sample:** CW-1-131-4-L **Blank:** CW-1-131-3-LB

**GC-EI-MS Method name:** TMS\_A



### METHOD DESCRIPTION

|   |  |   |            |
|---|--|---|------------|
| <b>Instrument Manufacturer and Type:</b>    | Agilent 6890/5973 GC/MSD   |   |            |
| <b>Carrier gas:</b>                         | <input checked="" type="checkbox"/> He <input type="checkbox"/> N <sub>2</sub> <input type="checkbox"/> H <sub>2</sub> <input type="checkbox"/> Other: |   |            |
| <b>Flow rate:</b>                           | <input type="checkbox"/> ml/min  | <input checked="" type="checkbox"/> 38 cm/s |            |
| <b>Flow control:</b>                        | <input type="checkbox"/> Constant Pressure <input checked="" type="checkbox"/> Constant Flow   |   |            |
| <b>Injection mode:</b>                      | <input type="checkbox"/> Split → Split ratio =<br><input checked="" type="checkbox"/> Splitless → Splitless time = 0.70 min.                           |   |            |
| <b>Injector temperature:</b>                | 250 °C   |   |            |
| <b>Column phase:</b>                        | 5% diphenyl 95% dimethyl polysiloxane  |   |            |
| <b>Column Length x ID x Film thickness:</b> | 30 m x 0.25 mm x 0.25 µm   |   |            |
| <b>GC temperature programme:</b>            | 70 °C (8 min), 8 °C/min, 300 °C (3 min)  |   |            |
| <b>Solvent delay time:</b>                  | 8 min  | <b>Scan range:</b>                          | 40-600 m/z |
| <b>Electron energy:</b>                     | 70 eV  | <b>Scan time:</b>                           | 0.7 s      |
| <b>Ionisation polarity:</b>                 | <input checked="" type="checkbox"/> Positive<br><input type="checkbox"/> Negative  | <b>Mass resolution:</b>                     | 0.7 u      |
| <b>Comments:</b>                            |  |   |            |



### ANALYSIS

|  |   |
|--|---|
| <b>Compound identified as:</b>   | <input type="checkbox"/> Original compound<br><input type="checkbox"/> Methyl ester derivative<br><input type="checkbox"/> TBDMS (t-Butyldimethylsilyl) derivative<br><input checked="" type="checkbox"/> TMS (Trimethylsilyl) derivative<br><input type="checkbox"/> Other derivative: |
| <b>Retention parameter used for (peak) identification:</b>                       | <input checked="" type="checkbox"/> Retention time (Rt)<br><input type="checkbox"/> Scan number   |
| <input checked="" type="checkbox"/> Compared to reference chemical:              | Source : <input type="checkbox"/> Own Synthesis <input checked="" type="checkbox"/> Commercial  |
| <input type="checkbox"/> Compared to library spectrum:                           | Source : <input type="checkbox"/> OCAD (code: ) <input type="checkbox"/> NIST<br><input type="checkbox"/> Wiley <input type="checkbox"/> Own <input type="checkbox"/> Other:  |
| <input type="checkbox"/> Not compared to reference chemical or library spectrum: | Intense ions in spectrum are explained; interpretation is supported by the spectral information derived from closely related chemical(s):   |
| <b>Comments:</b>   |   |



File : C:\DATA\16\J0003609.D  
Acquired: 29 Nov 2004 16:01 Sample : CW-1-131-3-LB  
Method : TMS\_A Misc info:



File : C:\DATA\16\J0003612.D  
Acquired: 29 Nov 2004 18:17 Sample : CW-1-131-4-L  
Method : TMS\_A Misc info:



File : C:\DATA\16\J0003614.D  
Acquired: 29 Nov 2004 19:47 Sample : CW-CK-1-128-2  
Method : TMS\_A Misc info:




EI chromatograms supporting identification of compound **6**; TIC on left; EIC (m/z **114**) on right.

Top: Chromatograms of Liquid blank, aliquot **CW-1-131-3-LB** from **LB/31**.


Center: Chromatograms of Liquid sample, aliquot **CW-1-131-4-L** from **L/31**, retention time **13.74** min.

Bottom: Chromatograms of TMS derivative of the standard of **2-(N,N-Diisopropylamino)ethanol** [2-(N,N-Diisopropylamino)ethyl trimethylsilyl ether], retention time **13.75** min.

File : C:\DATA\16\J0003612.D  
Acquired : 29 Nov 2004 18:17 using AcqMethod TMS\_A  
Sample Name: CW-1-131-4-L  
Misc Info :



File : C:\DATA\16\J0003614.D  
Acquired : 29 Nov 2004 19:47 using AcqMethod TMS\_A  
Sample Name: CW-CK-1-128-2  
Misc Info :



El mass spectrum of:

Top: Compound **6** in Liquid sample **L/31**, aliquot **CW-1-131-4-L**

Bottom: TMS derivative of the authentic reference standard of **2-(N,N-Diisopropylamino)ethanol**  
[2-(N,N-Diisopropylamino)ethyl trimethylsilyl ether] corresponding to compound **6**  
(MW: **217**)

# GC-CI-MS TECHNIQUE

## METHOD AND ANALYSIS DESCRIPTION

Laboratory code: 31 Sample code(s): L/31, LB/31 Compound number: 6

|                              |              |               |               |
|------------------------------|--------------|---------------|---------------|
| <b>Aliquot codes:</b>        |              |               |               |
| <b>Sample:</b>               | CW-1-131-4-L | <b>Blank:</b> | CW-1-131-3-LB |
| <b>GC-CI-MS Method name:</b> |              | CW-CI-TM      |               |

### METHOD DESCRIPTION

|   |  |   |   |
|---|--|---|---|
| <b>Instrument Manufacturer and Type:</b>    | Agilent 6890/5973 GC/MSD   |   |   |
| <b>Carrier gas:</b>                         | <input checked="" type="checkbox"/> He   | <input type="checkbox"/> N <sub>2</sub>           | <input type="checkbox"/> H <sub>2</sub> <input type="checkbox"/> Other:     |
| <b>Flow rate:</b>                           | <input type="checkbox"/> ml/min  | <input checked="" type="checkbox"/> 32 cm/s       |   |
| <b>Flow control:</b>                        | <input type="checkbox"/> Constant Pressure   | <input checked="" type="checkbox"/> Constant Flow |   |
| <b>Injection mode:</b>                      | <input type="checkbox"/> Split → Split ratio =<br><input checked="" type="checkbox"/> Splitless → Splitless time = 0.75 min. |   |   |
| <b>Injector temperature:</b>                | 250 °C   |   |   |
| <b>Column phase:</b>                        | 5% diphenyl 95% dimethyl polysiloxane  |   |   |
| <b>Column Length x ID x Film thickness:</b> | 30 m x 0.25 mm x 0.25 µm   |   |   |
| <b>GC temperature programme:</b>            | 70 °C (8 min), 8 °C/min, 300 °C (3 min)  |   |   |
| <b>Reaction gas:</b>                        | <input type="checkbox"/> Methane   | <input type="checkbox"/> Isobutane                | <input checked="" type="checkbox"/> Ammonia <input type="checkbox"/> Other: |
| <b>Solvent delay time:</b>                  | 8 min  | <b>Scan range:</b>                                | 50-550 m/z  |
| <b>Electron energy:</b>                     | 235 eV   | <b>Scan time:</b>                                 | 0.35 s  |
| <b>Ionisation polarity:</b>                 | <input checked="" type="checkbox"/> Positive<br><input type="checkbox"/> Negative  | <b>Mass resolution:</b>                           | 0.7 u   |
| <b>Comments:</b>                            |  |   |   |

### ANALYSIS

|  |   |
|--|---|
| <b>Compound identified as:</b>   | <input type="checkbox"/> Original compound<br><input type="checkbox"/> Methyl ester derivative<br><input type="checkbox"/> TBDMS (t-Butyldimethylsilyl) derivative<br><input checked="" type="checkbox"/> TMS (Trimethylsilyl) derivative<br><input type="checkbox"/> Other derivative: |
| <b>Retention parameter used for (peak) identification:</b>                       | <input checked="" type="checkbox"/> Retention time (Rt)<br><input type="checkbox"/> Scan number   |
| <input checked="" type="checkbox"/> Compared to reference chemical:              | Source : <input type="checkbox"/> Own Synthesis <input checked="" type="checkbox"/> Commercial  |
| <input type="checkbox"/> Compared to library spectrum:                           | Source : <input type="checkbox"/> OCAD (code: ) <input type="checkbox"/> NIST<br><input type="checkbox"/> Wiley <input type="checkbox"/> Own <input type="checkbox"/> Other:  |
| <input type="checkbox"/> Not compared to reference chemical or library spectrum: | Intense ions in spectrum are explained; interpretation is supported by the spectral information derived from closely related chemical(s):   |
| <b>Comments:</b>   |   |



File : C:\DATA\16\HP00857.D

Acquired: 23 Nov 2004 16:50

Method : CW-CI-TM

Sample : 1uL of CW-1-131-3-LB

Misc info: HP5-MS, 30m, 0.25 film, PCI-NH3





File : C:\DATA\16\HP00858.D

Acquired: 23 Nov 2004 17:35

Method : CW-CI-TM

Sample : 1uL of CW-1-131-4-L

Misc info: HP5-MS, 30m, 0.25 film, PCI-NH3





File : C:\DATA\16\HP00859.D

Acquired: 23 Nov 2004 18:20

Method : CW-CI-TM

Sample : 1uL of CW-CK-1-128-2 (compound G)

Misc info: HP5-MS, 30m, 0.25 film, PCI-NH3




CI chromatograms supporting identification of compound 6; TIC on left; EIC (m/z 218) on right.

Top: Chromatograms of Liquid blank, aliquot **CW-1-131-3-LB** from **LB/31**.


Center: Chromatograms of Liquid sample, aliquot **CW-1-131-4-L** from **L/31**, retention time **14.01** min.

Bottom: Chromatograms of TMS derivative of the standard of **2-(N,N-Diisopropylamino)ethanol** [2-(N,N-Diisopropylamino)ethyl trimethylsilyl ether], retention time **14.02** min.

File : C:\DATA\16\HP00858.D  
Acquired : 23 Nov 2004 17:35 using AcqMethod CW-CI-TM  
Sample Name: 1uL of CW-1-131-4-L  
Misc Info : HP5-MS, 30m, 0.25 film, PCI-NH3



File : C:\DATA\16\HP00859.D  
Acquired : 23 Nov 2004 18:20 using AcqMethod CW-CI-TM  
Sample Name: 1uL of CW-CK-1-128-2 (compound G)  
Misc Info : HP5-MS, 30m, 0.25 film, PCI-NH3



CI mass spectrum of:


Top: Compound **6** in Liquid sample **L/31**, aliquot **CW-1-131-4-L**

Bottom: TMS derivative of the authentic reference standard of **2-(N,N-Diisopropylamino)ethanol** [2-(N,N-Diisopropylamino)ethyl trimethylsilyl ether] corresponding to compound **6** (MW: **217**)

**SAMPLE PREPARATION DESCRIPTION**Laboratory code: 31 Sample code(s): S/31 Sample Blank code: SB/31**1. Sample preparation**

| <b>Sample/<br/>Aliquot Code</b> | <i>Specification of Sample/<br/>Type of Sample Preparation</i> | <b>Amount/<br/>Volume</b> | <i>Sample Preparation Procedures</i>   | <b>End<br/>Volume</b> | <b>Resulting<br/>Aliquot Code</b> | <b>Analytical<br/>Technique(s)</b> |
|---------------------------------|--|---------------------------|--|-----------------------|-----------------------------------|------------------------------------|
| S/31<br>SB/31                   | TMS derivative of methanol<br>extract of sample                | 2 g                       | 2 mL Methanol was added to 2 g of soil and sonicated<br>for 15 minutes followed by 4 minutes of centrifuge.<br>500µL of the supernate was dried (N <sub>2</sub> , no heat). 200 µL<br>BSTFA was added and heated at 60°C for 30 minutes. | 200 µL                | CW-1-131-2-S<br>CW-1-131-1-SB     | GC/MS-EI<br>GC/MS-CI               |
|                                 |  |                           |  |                       |                                   |                                    |
|                                 |  |                           |  |                       |                                   |                                    |
|                                 |  |                           |  |                       |                                   |                                    |
|                                 |  |                           |  |                       |                                   |                                    |
|                                 |  |                           |  |                       |                                   |                                    |
|                                 |  |                           |  |                       |                                   |                                    |

**2. Additional information**



Note: This flowchart is for visualization only; see the preceding sample preparation description page for sample aliquot numbers



# GC-EI-MS TECHNIQUE

## METHOD AND ANALYSIS DESCRIPTION

Laboratory code: 31 Sample code(s): S/31, SB/31 Compound number: 7

**Aliquot codes:**

**Sample:** CW-1-131-2-S **Blank:** CW-1-131-1-SB

**GC-EI-MS Method name:** TMS\_A



### METHOD DESCRIPTION

|   |  |   |            |
|---|--|---|------------|
| <b>Instrument Manufacturer and Type:</b>    | Agilent 6890/5973 GC/MSD   |   |            |
| <b>Carrier gas:</b>                         | <input checked="" type="checkbox"/> He <input type="checkbox"/> N <sub>2</sub> <input type="checkbox"/> H <sub>2</sub> <input type="checkbox"/> Other: |   |            |
| <b>Flow rate:</b>                           | <input type="checkbox"/> ml/min  | <input checked="" type="checkbox"/> 38 cm/s |            |
| <b>Flow control:</b>                        | <input type="checkbox"/> Constant Pressure <input checked="" type="checkbox"/> Constant Flow   |   |            |
| <b>Injection mode:</b>                      | <input type="checkbox"/> Split → Split ratio =<br><input checked="" type="checkbox"/> Splitless → Splitless time = 0.70 min.                           |   |            |
| <b>Injector temperature:</b>                | 250 °C   |   |            |
| <b>Column phase:</b>                        | 5% diphenyl 95% dimethyl polysiloxane  |   |            |
| <b>Column Length x ID x Film thickness:</b> | 30 m x 0.25 mm x 0.25 µm   |   |            |
| <b>GC temperature programme:</b>            | 70 °C (8 min), 8 °C/min, 300 °C (3 min)  |   |            |
| <b>Solvent delay time:</b>                  | 8 min  | <b>Scan range:</b>                          | 40-600 m/z |
| <b>Electron energy:</b>                     | 70 eV  | <b>Scan time:</b>                           | 0.7 s      |
| <b>Ionisation polarity:</b>                 | <input checked="" type="checkbox"/> Positive<br><input type="checkbox"/> Negative  | <b>Mass resolution:</b>                     | 0.7 u      |
| <b>Comments:</b>                            |  |   |            |



### ANALYSIS

|  |   |
|--|---|
| <b>Compound identified as:</b>   | <input type="checkbox"/> Original compound<br><input type="checkbox"/> Methyl ester derivative<br><input type="checkbox"/> TBDMS (t-Butyldimethylsilyl) derivative<br><input checked="" type="checkbox"/> TMS (Trimethylsilyl) derivative<br><input type="checkbox"/> Other derivative: |
| <b>Retention parameter used for (peak) identification:</b>                       | <input checked="" type="checkbox"/> Retention time (Rt)<br><input type="checkbox"/> Scan number   |
| <input checked="" type="checkbox"/> Compared to reference chemical:              | Source : <input type="checkbox"/> Own Synthesis <input checked="" type="checkbox"/> Commercial  |
| <input type="checkbox"/> Compared to library spectrum:                           | Source : <input type="checkbox"/> OCAD (code: ) <input type="checkbox"/> NIST<br><input type="checkbox"/> Wiley <input type="checkbox"/> Own <input type="checkbox"/> Other:  |
| <input type="checkbox"/> Not compared to reference chemical or library spectrum: | Intense ions in spectrum are explained; interpretation is supported by the spectral information derived from closely related chemical(s):   |
| <b>Comments:</b>   |   |



File : C:\DATA\16\J0003603.D  
Acquired: 29 Nov 2004 11:29 Sample : CW-1-131-1-SB  
Method : TMS\_A Misc info:



File : C:\DATA\16\J0003604.D  
Acquired: 29 Nov 2004 12:14 Sample : CW-1-131-2-S  
Method : TMS\_A Misc info:



File : C:\DATA\16\J0003606.D  
Acquired: 29 Nov 2004 13:45 Sample : CW-CK-1-128-3  
Method : TMS\_A Misc info:




EI chromatograms supporting identification of compound **7**; TIC on left; EIC (m/z **174**) on right.

Top: Chromatograms of Soil blank, aliquot **CW-1-131-1-SB** from **SB/31**.


Center: Chromatograms of Soil sample, aliquot **CW-1-131-2-S** from **S/31**, retention time **17.62** min.

Bottom: Chromatograms of TMS derivative of the authentic reference standard of **Ethyl-diethanolamine** [Bis(2-trimethylsilyloxyethyl)ethylamine], retention time **17.63** min.

File : C:\DATA\16\J0003604.D  
Acquired : 29 Nov 2004 12:14 using AcqMethod TMS\_A  
Sample Name: CW-1-131-2-S  
Misc Info :



File : C:\DATA\16\J0003606.D  
Acquired : 29 Nov 2004 13:45 using AcqMethod TMS\_A  
Sample Name: CW-CK-1-128-3  
Misc Info :



El mass spectrum of:

Top: Compound **7** in Soil sample **S/31**, aliquot **CW-1-131-2-S**

Bottom: TMS derivative of the authentic reference standard of **Ethyl**diethanolamine  
[Bis(2-trimethylsilyloxyethyl)ethylamine] corresponding to compound **7** (MW: **277**)

# GC-CI-MS TECHNIQUE

## METHOD AND ANALYSIS DESCRIPTION

Laboratory code: 31 Sample code(s): S/31, SB/31 Compound number: 7

Aliquot codes:

Sample: CW-1-131-2-S Blank: CW-1-131-1-SB

GC-CI-MS Method name: CW-CI-TM

### METHOD DESCRIPTION

|   |  |   |   |
|---|--|---|---|
| <b>Instrument Manufacturer and Type:</b>    | Agilent 6890/5973 GC/MSD   |   |   |
| <b>Carrier gas:</b>                         | <input checked="" type="checkbox"/> He   | <input type="checkbox"/> N <sub>2</sub>           | <input type="checkbox"/> H <sub>2</sub> <input type="checkbox"/> Other:     |
| <b>Flow rate:</b>                           | <input type="checkbox"/> ml/min  | <input checked="" type="checkbox"/> 32 cm/s       |   |
| <b>Flow control:</b>                        | <input type="checkbox"/> Constant Pressure   | <input checked="" type="checkbox"/> Constant Flow |   |
| <b>Injection mode:</b>                      | <input type="checkbox"/> Split → Split ratio =<br><input checked="" type="checkbox"/> Splitless → Splitless time = 0.75 min. |   |   |
| <b>Injector temperature:</b>                | 250 °C   |   |   |
| <b>Column phase:</b>                        | 5% diphenyl 95% dimethyl polysiloxane  |   |   |
| <b>Column Length x ID x Film thickness:</b> | 30 m x 0.25 mm x 0.25 µm   |   |   |
| <b>GC temperature programme:</b>            | 70 °C (8 min), 8 °C/min, 300 °C (3 min)  |   |   |
| <b>Reaction gas:</b>                        | <input type="checkbox"/> Methane   | <input type="checkbox"/> Isobutane                | <input checked="" type="checkbox"/> Ammonia <input type="checkbox"/> Other: |
| <b>Solvent delay time:</b>                  | 8 min  | <b>Scan range:</b>                                | 50-550 m/z  |
| <b>Electron energy:</b>                     | 235 eV   | <b>Scan time:</b>                                 | 0.35 s  |
| <b>Ionisation polarity:</b>                 | <input checked="" type="checkbox"/> Positive<br><input type="checkbox"/> Negative  | <b>Mass resolution:</b>                           | 0.7 u   |
| <b>Comments:</b>                            |  |   |   |

### ANALYSIS

|  |   |
|--|---|
| <b>Compound identified as:</b>   | <input type="checkbox"/> Original compound<br><input type="checkbox"/> Methyl ester derivative<br><input type="checkbox"/> TBDMS (t-Butyldimethylsilyl) derivative<br><input checked="" type="checkbox"/> TMS (Trimethylsilyl) derivative<br><input type="checkbox"/> Other derivative: |
| <b>Retention parameter used for (peak) identification:</b>                       | <input checked="" type="checkbox"/> Retention time (Rt)<br><input type="checkbox"/> Scan number   |
| <input checked="" type="checkbox"/> Compared to reference chemical:              | Source : <input type="checkbox"/> Own Synthesis <input checked="" type="checkbox"/> Commercial  |
| <input type="checkbox"/> Compared to library spectrum:                           | Source : <input type="checkbox"/> OCAD (code: ) <input type="checkbox"/> NIST<br><input type="checkbox"/> Wiley <input type="checkbox"/> Own <input type="checkbox"/> Other:  |
| <input type="checkbox"/> Not compared to reference chemical or library spectrum: | Intense ions in spectrum are explained; interpretation is supported by the spectral information derived from closely related chemical(s):   |
| <b>Comments:</b>   |   |



File : C:\DATA\16\HP00853.D

Acquired: 23 Nov 2004 13:51

Method : CW-CI-TM

Sample : 1uL of CW-1-131-1-SB

Misc info: HP5-MS, 30m, 0.25 film, PCI-NH3





File : C:\DATA\16\HP00854.D

Acquired: 23 Nov 2004 14:36

Method : CW-CI-TM

Sample : 1uL of CW-1-131-2-S

Misc info: HP5-MS, 30m, 0.25 film, PCI-NH3





File : C:\DATA\16\HP00855.D

Acquired: 23 Nov 2004 15:21

Method : CW-CI-TM

Sample : 1uL of CW-CK-1-128-3 (compound H)

Misc info: HP5-MS, 30m, 0.25 film, PCI-NH3




CI chromatograms supporting identification of compound 7; TIC on left; EIC (m/z 278) on right.

Top: Chromatograms of Soil blank, aliquot **CW-1-131-1-SB** from **SB/31**.


Center: Chromatograms of Soil sample, aliquot **CW-1-131-2-S** from **S/31**, retention time **17.85** min.

Bottom: Chromatograms of TMS derivative of the authentic reference standard of **Ethyl-diethanolamine** [Bis(2-trimethylsilyloxyethyl)ethylamine], retention time **17.85** min.

File : C:\DATA\16\HP00854.D  
Acquired : 23 Nov 2004 14:36 using AcqMethod CW-CI-TM  
Sample Name: 1uL of CW-1-131-2-S  
Misc Info : HP5-MS, 30m, 0.25 film, PCI-NH3



File : C:\DATA\16\HP00855.D  
Acquired : 23 Nov 2004 15:21 using AcqMethod CW-CI-TM  
Sample Name: 1uL of CW-CK-1-128-3 (compound H)  
Misc Info : HP5-MS, 30m, 0.25 film, PCI-NH3



CI mass spectrum of:

Top: Compound **7** in Soil sample **S/31**, aliquot **CW-1-131-2-S**

Bottom: TMS derivative of the authentic reference standard of **Ethyldiethanolamine**  
[Bis(2-trimethylsilyloxyethyl)ethylamine]corresponding to compound **7** (MW: 277)

# GC-EI-MS TECHNIQUE

## METHOD AND ANALYSIS DESCRIPTION

Laboratory code: 31 Sample code(s): S/31, SB/31 Compound number: 8

Aliquot codes:

Sample: CW-1-131-2-S Blank: CW-1-131-1-SB

GC-EI-MS Method name: TMS\_A



### METHOD DESCRIPTION

|   |  |   |            |
|---|--|---|------------|
| <b>Instrument Manufacturer and Type:</b>    | Agilent 6890/5973 GC/MSD   |   |            |
| <b>Carrier gas:</b>                         | <input checked="" type="checkbox"/> He <input type="checkbox"/> N <sub>2</sub> <input type="checkbox"/> H <sub>2</sub> <input type="checkbox"/> Other: |   |            |
| <b>Flow rate:</b>                           | <input type="checkbox"/> ml/min  | <input checked="" type="checkbox"/> 38 cm/s |            |
| <b>Flow control:</b>                        | <input type="checkbox"/> Constant Pressure <input checked="" type="checkbox"/> Constant Flow   |   |            |
| <b>Injection mode:</b>                      | <input type="checkbox"/> Split → Split ratio =<br><input checked="" type="checkbox"/> Splitless → Splitless time = 0.70 min.                           |   |            |
| <b>Injector temperature:</b>                | 250 °C   |   |            |
| <b>Column phase:</b>                        | 5% diphenyl 95% dimethyl polysiloxane  |   |            |
| <b>Column Length x ID x Film thickness:</b> | 30 m x 0.25 mm x 0.25 µm   |   |            |
| <b>GC temperature programme:</b>            | 70 °C (8 min), 8 °C/min, 300 °C (3 min)  |   |            |
| <b>Solvent delay time:</b>                  | 8 min  | <b>Scan range:</b>                          | 40-600 m/z |
| <b>Electron energy:</b>                     | 70 eV  | <b>Scan time:</b>                           | 0.7 s      |
| <b>Ionisation polarity:</b>                 | <input checked="" type="checkbox"/> Positive<br><input type="checkbox"/> Negative  | <b>Mass resolution:</b>                     | 0.7 u      |
| <b>Comments:</b>                            |  |   |            |



### ANALYSIS

|  |   |
|--|---|
| <b>Compound identified as:</b>   | <input type="checkbox"/> Original compound<br><input type="checkbox"/> Methyl ester derivative<br><input type="checkbox"/> TBDMS (t-Butyldimethylsilyl) derivative<br><input checked="" type="checkbox"/> TMS (Trimethylsilyl) derivative<br><input type="checkbox"/> Other derivative: |
| <b>Retention parameter used for (peak) identification:</b>                       | <input checked="" type="checkbox"/> Retention time (Rt)<br><input type="checkbox"/> Scan number   |
| <input checked="" type="checkbox"/> Compared to reference chemical:              | Source : <input type="checkbox"/> Own Synthesis <input checked="" type="checkbox"/> Commercial  |
| <input type="checkbox"/> Compared to library spectrum:                           | Source : <input type="checkbox"/> OCAD (code: ) <input type="checkbox"/> NIST<br><input type="checkbox"/> Wiley <input type="checkbox"/> Own <input type="checkbox"/> Other:  |
| <input type="checkbox"/> Not compared to reference chemical or library spectrum: | Intense ions in spectrum are explained; interpretation is supported by the spectral information derived from closely related chemical(s):   |
| <b>Comments:</b>   |   |



File : C:\DATA\16\J0003603.D  
Acquired: 29 Nov 2004 11:29 Sample : CW-1-131-1-SB  
Method : TMS\_A Misc info:



File : C:\DATA\16\J0003604.D  
Acquired: 29 Nov 2004 12:14 Sample : CW-1-131-2-S  
Method : TMS\_A Misc info:



File : C:\DATA\16\J0003607.D  
Acquired: 29 Nov 2004 14:30 Sample : CW-CK-1-128-4  
Method : TMS\_A Misc info:



EI chromatograms supporting identification of compound **8**; TIC on left; EIC (m/z **160**) on right.


Top: Chromatograms of Soil blank, aliquot **CW-1-131-1-SB** from **SB/31**.

Center: Chromatograms of Soil sample, aliquot **CW-1-131-2-S** from **S/31**, retention time **16.69** min.


Bottom: Chromatograms of TMS derivative of the authentic reference standard of **Methyldiethanolamine** [Bis(2-trimethylsilyloxyethyl)methylamine], retention time **16.69** min.



File : C:\DATA\16\J0003604.D  
Acquired : 29 Nov 2004 12:14 using AcqMethod TMS\_A  
Sample Name: CW-1-131-2-S  
Misc Info :



File : C:\DATA\16\J0003607.D  
Acquired : 29 Nov 2004 14:30 using AcqMethod TMS\_A  
Sample Name: CW-CK-1-128-4  
Misc Info :



El mass spectrum of:

Top: Compound **8** in Soil sample **S/31**, aliquot **CW-1-131-2-S**

Bottom: TMS derivative of the authentic reference standard of **Methyldiethanolamine**  
[Bis(2-trimethylsilyloxyethyl)methylamine] corresponding to compound **8** (MW: **263**)

# GC-CI-MS TECHNIQUE

## METHOD AND ANALYSIS DESCRIPTION

Laboratory code: 31 Sample code(s): S/31, SB/31 Compound number: 8

|                              |              |               |               |
|------------------------------|--------------|---------------|---------------|
| <b>Aliquot codes:</b>        |              |               |               |
| <b>Sample:</b>               | CW-1-131-2-S | <b>Blank:</b> | CW-1-131-1-SB |
| <b>GC-CI-MS Method name:</b> |              | CW-CI-TM      |               |

### METHOD DESCRIPTION

|   |  |   |   |
|---|--|---|---|
| <b>Instrument Manufacturer and Type:</b>    | Agilent 6890/5973 GC/MSD   |   |   |
| <b>Carrier gas:</b>                         | <input checked="" type="checkbox"/> He   | <input type="checkbox"/> N <sub>2</sub>           | <input type="checkbox"/> H <sub>2</sub> <input type="checkbox"/> Other:     |
| <b>Flow rate:</b>                           | <input type="checkbox"/> ml/min  | <input checked="" type="checkbox"/> 32 cm/s       |   |
| <b>Flow control:</b>                        | <input type="checkbox"/> Constant Pressure   | <input checked="" type="checkbox"/> Constant Flow |   |
| <b>Injection mode:</b>                      | <input type="checkbox"/> Split → Split ratio =<br><input checked="" type="checkbox"/> Splitless → Splitless time = 0.75 min. |   |   |
| <b>Injector temperature:</b>                | 250 °C   |   |   |
| <b>Column phase:</b>                        | 5% diphenyl 95% dimethyl polysiloxane  |   |   |
| <b>Column Length x ID x Film thickness:</b> | 30 m x 0.25 mm x 0.25 µm   |   |   |
| <b>GC temperature programme:</b>            | 70 °C (8 min), 8 °C/min, 300 °C (3 min)  |   |   |
| <b>Reaction gas:</b>                        | <input type="checkbox"/> Methane   | <input type="checkbox"/> Isobutane                | <input checked="" type="checkbox"/> Ammonia <input type="checkbox"/> Other: |
| <b>Solvent delay time:</b>                  | 8 min  | <b>Scan range:</b>                                | 50-550 m/z  |
| <b>Electron energy:</b>                     | 235 eV   | <b>Scan time:</b>                                 | 0.35 s  |
| <b>Ionisation polarity:</b>                 | <input checked="" type="checkbox"/> Positive<br><input type="checkbox"/> Negative  | <b>Mass resolution:</b>                           | 0.7 u   |
| <b>Comments:</b>                            |  |   |   |

### ANALYSIS

|  |   |
|--|---|
| <b>Compound identified as:</b>   | <input type="checkbox"/> Original compound<br><input type="checkbox"/> Methyl ester derivative<br><input type="checkbox"/> TBDMS (t-Butyldimethylsilyl) derivative<br><input checked="" type="checkbox"/> TMS (Trimethylsilyl) derivative<br><input type="checkbox"/> Other derivative: |
| <b>Retention parameter used for (peak) identification:</b>                       | <input checked="" type="checkbox"/> Retention time (Rt)<br><input type="checkbox"/> Scan number   |
| <input checked="" type="checkbox"/> Compared to reference chemical:              | Source : <input type="checkbox"/> Own Synthesis <input checked="" type="checkbox"/> Commercial  |
| <input type="checkbox"/> Compared to library spectrum:                           | Source : <input type="checkbox"/> OCAD (code: ) <input type="checkbox"/> NIST<br><input type="checkbox"/> Wiley <input type="checkbox"/> Own <input type="checkbox"/> Other:  |
| <input type="checkbox"/> Not compared to reference chemical or library spectrum: | Intense ions in spectrum are explained; interpretation is supported by the spectral information derived from closely related chemical(s):   |
| <b>Comments:</b>   |   |



File : C:\DATA\16\HP00853.D

Acquired: 23 Nov 2004 13:51

Sample : 1uL of CW-1-131-1-SB

Method : CW-CI-TM

Misc info: HP5-MS, 30m, 0.25 film, PCI-NH3





File : C:\DATA\16\HP00854.D

Acquired: 23 Nov 2004 14:36

Sample : 1uL of CW-1-131-2-S

Method : CW-CI-TM

Misc info: HP5-MS, 30m, 0.25 film, PCI-NH3





File : C:\DATA\16\HP00856.D

Acquired: 23 Nov 2004 16:05

Sample : 1uL of CW-CK-1-128-4 (compound I)

Method : CW-CI-TM

Misc info: HP5-MS, 30m, 0.25 film, PCI-NH3




CI chromatograms supporting identification of compound **8**; TIC on left; EIC (m/z **264**) on right.

Top: Chromatograms of Soil blank, aliquot **CW-1-131-1-SB** from **SB/31**.


Center: Chromatograms of Soil sample, aliquot **CW-1-131-2-S** from **S/31**, retention time **16.93** min.

Bottom: Chromatograms of TMS derivative of the authentic reference standard of **Methyldiethanolamine** [Bis(2-trimethylsilyloxyethyl)methylamine], retention time **16.93** min.

File : C:\DATA\16\HP00854.D  
Acquired : 23 Nov 2004 14:36 using AcqMethod CW-CI-TM  
Sample Name: 1uL of CW-1-131-2-S  
Misc Info : HP5-MS, 30m, 0.25 film, PCI-NH3



File : C:\DATA\16\HP00856.D  
Acquired : 23 Nov 2004 16:05 using AcqMethod CW-CI-TM  
Sample Name: 1uL of CW-CK-1-128-4 (compound I)  
Misc Info : HP5-MS, 30m, 0.25 film, PCI-NH3



CI mass spectrum of:

Top: Compound **8** in Soil sample **S/31**, aliquot **CW-1-131-2-S**

Bottom: TMS derivative of the authentic reference standard of **Methyldiethanolamine**  
[Bis(2-trimethylsilyloxyethyl)methylamine] corresponding to compound **8** (MW: 263)

# GC-EI-MS TECHNIQUE

## METHOD AND ANALYSIS DESCRIPTION

Laboratory code: 31 Sample code(s): S/31, SB/31 Compound number: 9

**Aliquot codes:**

**Sample:** CW-1-131-2-S **Blank:** CW-1-131-1-SB

**GC-EI-MS Method name:** TMS\_A



### METHOD DESCRIPTION

|   |  |   |            |
|---|--|---|------------|
| <b>Instrument Manufacturer and Type:</b>    | Agilent 6890/5973 GC/MSD   |   |            |
| <b>Carrier gas:</b>                         | <input checked="" type="checkbox"/> He <input type="checkbox"/> N <sub>2</sub> <input type="checkbox"/> H <sub>2</sub> <input type="checkbox"/> Other: |   |            |
| <b>Flow rate:</b>                           | <input type="checkbox"/> ml/min  | <input checked="" type="checkbox"/> 38 cm/s |            |
| <b>Flow control:</b>                        | <input type="checkbox"/> Constant Pressure <input checked="" type="checkbox"/> Constant Flow   |   |            |
| <b>Injection mode:</b>                      | <input type="checkbox"/> Split → Split ratio =<br><input checked="" type="checkbox"/> Splitless → Splitless time = 0.70 min.                           |   |            |
| <b>Injector temperature:</b>                | 250 °C   |   |            |
| <b>Column phase:</b>                        | 5% diphenyl 95% dimethyl polysiloxane  |   |            |
| <b>Column Length x ID x Film thickness:</b> | 30 m x 0.25 mm x 0.25 µm   |   |            |
| <b>GC temperature programme:</b>            | 70 °C (8 min), 8 °C/min, 300 °C (3 min)  |   |            |
| <b>Solvent delay time:</b>                  | 8 min  | <b>Scan range:</b>                          | 40-600 m/z |
| <b>Electron energy:</b>                     | 70 eV  | <b>Scan time:</b>                           | 0.7 s      |
| <b>Ionisation polarity:</b>                 | <input checked="" type="checkbox"/> Positive<br><input type="checkbox"/> Negative  | <b>Mass resolution:</b>                     | 0.7 u      |
| <b>Comments:</b>                            |  |   |            |



### ANALYSIS

|  |   |
|--|---|
| <b>Compound identified as:</b>   | <input type="checkbox"/> Original compound<br><input type="checkbox"/> Methyl ester derivative<br><input type="checkbox"/> TBDMS (t-Butyldimethylsilyl) derivative<br><input checked="" type="checkbox"/> TMS (Trimethylsilyl) derivative<br><input type="checkbox"/> Other derivative: |
| <b>Retention parameter used for (peak) identification:</b>                       | <input checked="" type="checkbox"/> Retention time (Rt)<br><input type="checkbox"/> Scan number   |
| <input checked="" type="checkbox"/> Compared to reference chemical:              | Source : <input type="checkbox"/> Own Synthesis <input checked="" type="checkbox"/> Commercial  |
| <input type="checkbox"/> Compared to library spectrum:                           | Source : <input type="checkbox"/> OCAD (code: ) <input type="checkbox"/> NIST<br><input type="checkbox"/> Wiley <input type="checkbox"/> Own <input type="checkbox"/> Other:  |
| <input type="checkbox"/> Not compared to reference chemical or library spectrum: | Intense ions in spectrum are explained; interpretation is supported by the spectral information derived from closely related chemical(s):   |
| <b>Comments:</b>   |   |



File : C:\DATA\16\J0003603.D  
Acquired: 29 Nov 2004 11:29 Sample : CW-1-131-1-SB  
Method : TMS\_A Misc info:



File : C:\DATA\16\J0003604.D  
Acquired: 29 Nov 2004 12:14 Sample : CW-1-131-2-S  
Method : TMS\_A Misc info:



File : C:\DATA\16\J0003605.D  
Acquired: 29 Nov 2004 12:59 Sample : CW-CK-1-127-4  
Method : TMS\_A Misc info:




EI chromatograms supporting identification of compound **9**; TIC on left; EIC (m/z **262**) on right.

Top: Chromatograms of Soil blank, aliquot **CW-1-131-1-SB** from **SB/31**.


Center: Chromatograms of Soil sample, aliquot **CW-1-131-2-S** from **S/31**, retention time **22.07** min.

Bottom: Chromatograms of TMS derivative of the authentic reference standard of **Triethanolamine** [tris(2-trimethylsiloxyethyl)amine], retention time **22.07** min.

File : C:\DATA\16\J0003604.D  
Acquired : 29 Nov 2004 12:14 using AcqMethod TMS\_A  
Sample Name: CW-1-131-2-S  
Misc Info :



File : C:\DATA\16\J0003605.D  
Acquired : 29 Nov 2004 12:59 using AcqMethod TMS\_A  
Sample Name: CW-CK-1-127-4  
Misc Info :



El mass spectrum of:

Top: Compound **9** in Soil sample **S/31**, aliquot **CW-1-131-2-S**

Bottom: TMS derivative of the authentic reference standard of **Triethanolamine**  
[tris(2-trimethylsiloxyethyl)amine] corresponding to compound **9** (MW: **365**)

# GC-CI-MS TECHNIQUE

## METHOD AND ANALYSIS DESCRIPTION

Laboratory code: 31 Sample code(s): S/31, SB/31 Compound number: 9

|                              |              |               |               |
|------------------------------|--------------|---------------|---------------|
| <b>Aliquot codes:</b>        |              |               |               |
| <b>Sample:</b>               | CW-1-131-2-S | <b>Blank:</b> | CW-1-131-1-SB |
| <b>GC-CI-MS Method name:</b> |              | CW-CI-TM      |               |

### METHOD DESCRIPTION

|   |  |   |   |
|---|--|---|---|
| <b>Instrument Manufacturer and Type:</b>    | Agilent 6890/5973 GC/MSD   |   |   |
| <b>Carrier gas:</b>                         | <input checked="" type="checkbox"/> He   | <input type="checkbox"/> N <sub>2</sub>           | <input type="checkbox"/> H <sub>2</sub> <input type="checkbox"/> Other:     |
| <b>Flow rate:</b>                           | <input type="checkbox"/> ml/min  | <input checked="" type="checkbox"/> 32 cm/s       |   |
| <b>Flow control:</b>                        | <input type="checkbox"/> Constant Pressure   | <input checked="" type="checkbox"/> Constant Flow |   |
| <b>Injection mode:</b>                      | <input type="checkbox"/> Split → Split ratio =<br><input checked="" type="checkbox"/> Splitless → Splitless time = 0.75 min. |   |   |
| <b>Injector temperature:</b>                | 250 °C   |   |   |
| <b>Column phase:</b>                        | 5% diphenyl 95% dimethyl polysiloxane  |   |   |
| <b>Column Length x ID x Film thickness:</b> | 30 m x 0.25 mm x 0.25 µm   |   |   |
| <b>GC temperature programme:</b>            | 70 °C (8 min), 8 °C/min, 300 °C (3 min)  |   |   |
| <b>Reaction gas:</b>                        | <input type="checkbox"/> Methane   | <input type="checkbox"/> Isobutane                | <input checked="" type="checkbox"/> Ammonia <input type="checkbox"/> Other: |
| <b>Solvent delay time:</b>                  | 8 min  | <b>Scan range:</b>                                | 50-550 m/z  |
| <b>Electron energy:</b>                     | 235 eV   | <b>Scan time:</b>                                 | 0.35 s  |
| <b>Ionisation polarity:</b>                 | <input checked="" type="checkbox"/> Positive<br><input type="checkbox"/> Negative  | <b>Mass resolution:</b>                           | 0.7 u   |
| <b>Comments:</b>                            |  |   |   |

### ANALYSIS

|  |   |
|--|---|
| <b>Compound identified as:</b>   | <input type="checkbox"/> Original compound<br><input type="checkbox"/> Methyl ester derivative<br><input type="checkbox"/> TBDMS (t-Butyldimethylsilyl) derivative<br><input checked="" type="checkbox"/> TMS (Trimethylsilyl) derivative<br><input type="checkbox"/> Other derivative: |
| <b>Retention parameter used for (peak) identification:</b>                       | <input checked="" type="checkbox"/> Retention time (Rt)<br><input type="checkbox"/> Scan number   |
| <input checked="" type="checkbox"/> Compared to reference chemical:              | Source : <input type="checkbox"/> Own Synthesis <input checked="" type="checkbox"/> Commercial  |
| <input type="checkbox"/> Compared to library spectrum:                           | Source : <input type="checkbox"/> OCAD (code: ) <input type="checkbox"/> NIST<br><input type="checkbox"/> Wiley <input type="checkbox"/> Own <input type="checkbox"/> Other:  |
| <input type="checkbox"/> Not compared to reference chemical or library spectrum: | Intense ions in spectrum are explained; interpretation is supported by the spectral information derived from closely related chemical(s):   |
| <b>Comments:</b>   |   |





File : C:\DATA\16\HP00853.D

Acquired: 23 Nov 2004 13:51

Sample : 1uL of CW-1-131-1-SB

Method : CW-CI-TM

Misc info: HP5-MS, 30m, 0.25 film, PCI-NH3





File : C:\DATA\16\HP00854.D

Acquired: 23 Nov 2004 14:36

Sample : 1uL of CW-1-131-2-S

Method : CW-CI-TM

Misc info: HP5-MS, 30m, 0.25 film, PCI-NH3





File : C:\DATA\16\HP00862.D

Acquired: 24 Nov 2004 18:54

Sample : 1uL of CW-CK-1-127-4 (compound J)

Method : CW-CI-TM

Misc info: HP5-MS, 30m, 0.25 film, PCI-NH3




CI chromatograms supporting identification of compound **9**; TIC on left; EIC (m/z **366**) on right.

Top: Chromatograms of Soil blank, aliquot **CW-1-131-1-SB** from **SB/31**.


Center: Chromatograms of Soil sample, aliquot **CW-1-131-2-S** from **S/31**, retention time **22.29** min.

Bottom: Chromatograms of TMS derivative of the authentic reference standard of **Triethanolamine** [tris(2-trimethylsiloxyethyl)amine], retention time **22.29** min.

File : C:\DATA\16\HP00854.D  
Acquired : 23 Nov 2004 14:36 using AcqMethod CW-CI-TM  
Sample Name: 1uL of CW-1-131-2-S  
Misc Info : HP5-MS, 30m, 0.25 film, PCI-NH3



File : C:\DATA\16\HP00862.D  
Acquired : 24 Nov 2004 18:54 using AcqMethod CW-CI-TM  
Sample Name: 1uL of CW-CK-1-127-4 (compound J)  
Misc Info : HP5-MS, 30m, 0.25 film, PCI-NH3



CI mass spectrum of:

Top: Compound **9** in Soil sample **S/31**, aliquot **CW-1-131-2-S**

Bottom: TMS derivative of the authentic reference standard of **Triethanolamine**  
[tris(2-trimethylsiloxyethyl)amine] corresponding to compound **9** (MW: **365**)

[tris(2-

**COMMENTS****1. General****2. Sample preparation****3. Analysis****4. Report**

The field length in the Chemstation software often did not allow for the full chemical name and may be truncated or abbreviated.