



LAWRENCE
LIVERMORE
NATIONAL
LABORATORY

Chemical Imaging of Lipid Domains by High-Resolution Secondary Ion Mass Spectrometry

M. L. Kraft, P. K. Weber, M. L. Longo, I. D.
Hutcheon, S. G. Boxer

October 3, 2005

Biophysical Society 50th Annual Meeting
Salt Lake, UT, United States
February 18, 2005 through February 22, 2005

This document was prepared as an account of work sponsored by an agency of the United States Government. Neither the United States Government nor the University of California nor any of their employees, makes any warranty, express or implied, or assumes any legal liability or responsibility for the accuracy, completeness, or usefulness of any information, apparatus, product, or process disclosed, or represents that its use would not infringe privately owned rights. Reference herein to any specific commercial product, process, or service by trade name, trademark, manufacturer, or otherwise, does not necessarily constitute or imply its endorsement, recommendation, or favoring by the United States Government or the University of California. The views and opinions of authors expressed herein do not necessarily state or reflect those of the United States Government or the University of California, and shall not be used for advertising or product endorsement purposes.

Chemical Imaging of Lipid Domains by High-Resolution Secondary Ion Mass Spectrometry

Mary L. Kraft^{*}, Peter K. Weber[†], Marjorie L. Longo[‡], Ian D. Hutcheon[†], Steven G. Boxer^{*}

^{*}Department of Chemistry, Stanford University, Stanford, CA 94305

[†]GTSI & BSNL, Lawrence Livermore National Laboratory, Livermore, CA 94551

[‡]Department of Chemical Engineering and Material Science, University of California, Davis, CA 95616

Lipid microdomains within supported lipid bilayers composed of binary phosphocholine mixtures were chemically imaged by high-resolution secondary ion mass spectrometry performed with the NanoSIMS 50 (Cameca Instruments). This instrument images the sample components based on the elemental or isotopic composition of their atomic and small molecular secondary ions. Up to five different secondary ions can be simultaneously detected, and a lateral resolution of 50 nm can be achieved with high sensitivity at high mass resolution. In our experiments, the NanoSIMS 50 extensively fragmented the supported membrane, therefore an isotopic labeling strategy was used to encode the identities of the lipid components. Supported lipid membranes that contained distinct lipid microdomains were freeze-dried to preserve their lateral organization and analyzed with the NanoSIMS 50. Lipid microdomains as small as 100 nm in diameter were successfully imaged, and this was validated by comparison to AFM images taken at the same region prior to chemical imaging. Quantitative information on the lipid distribution within the domain was also determined by calibrating against supported membranes of known composition. We believe this will be a valuable approach for analyzing the composition of complex membrane domains with high spatial resolution.

Portions of this work were performed under the auspices of the U.S. Department of Energy by the University of California, Lawrence Livermore National Laboratory under Contract W-7405-Eng-48.