



LAWRENCE
LIVERMORE
NATIONAL
LABORATORY

X-Ray Digital Radiography and Computed Tomography of High Energy Density Physics (HEDP) Material, Subassemblies and Targets

William D Brown, Harry E Martz Jr.

July 19, 2006

ASNT Digital Imaging IX Topical Conference
Mashantucket, CT, United States
July 24, 2006 through July 26, 2006

Disclaimer

This document was prepared as an account of work sponsored by an agency of the United States Government. Neither the United States Government nor the University of California nor any of their employees, makes any warranty, express or implied, or assumes any legal liability or responsibility for the accuracy, completeness, or usefulness of any information, apparatus, product, or process disclosed, or represents that its use would not infringe privately owned rights. Reference herein to any specific commercial product, process, or service by trade name, trademark, manufacturer, or otherwise, does not necessarily constitute or imply its endorsement, recommendation, or favoring by the United States Government or the University of California. The views and opinions of authors expressed herein do not necessarily state or reflect those of the United States Government or the University of California, and shall not be used for advertising or product endorsement purposes.

X-ray Digital Radiography and Computed Tomography of High Energy Density Physics (HEDP) Material, Subassemblies and Targets

Presented by:


William D. Brown

Lawrence Livermore National Laboratory


Collaborator:

Harry Martz

UCRL-PRES-



DR image (I/I_0)



CT image (x-y plane)

**ASNT Digital Imaging IX Topical Conference
Mashantucket, CT
July 24, 2006**




UCRL-PRES-#

Presentation outline




- **LLNL DR/CT mesoscale systems**
- **X-ray DR system performance**
 - MTF
 - SNR
- **Mesoscale phase effects**
- **Programmatic applications**
 - Low-Temperature Raleigh-Taylor Targets
 - Double-Shell Target
- **Summary**

NIF will use 192 focused beams to perform energy research and HEDP experiments



We are benchmarking several different x-ray imaging systems for nondestructive characterization of HEDP targets



System Comparison

	LLNL KCAT	ALS Synchrotron	Xradia μ XCT
X-ray Source	Keve μ focus	Synchrotron	Hamamatsu μ -focus
Energy	20 - 160 kV	6 - 30 keV	20 - 150 kV
Detector type	CCD	CCD	CCD
Scintillator	Tb doped glass	Cd ₂ WO ₄	CsI
Camera pitch (μ m)	9	9	13.5
Pixel size at detector (μ m)	3	9.0 - 1.7	2.7/1.35/0.68

We are using MTF and SNR parameters of an edge as a measure of x-ray DR system performance





- **Modulation Transfer Function (MTF)**

- Frequency domain description of spatial resolution and contrast
- System MTF is the product of components


- $MTF_{\text{system}} = MTF_{\text{scintillator}} * MTF_{\text{optics}} * MTF_{\text{detector}}$

- **Signal to Noise Ratio (SNR)**


- Contrast measurement of images
- SNR factors
 - Source noise
 - Detector noise
 - Scattered radiation



Edge



Derivative



Fourier transform

$$SNR = \frac{S_1 - S_2}{\sqrt{\sigma_1^2 + \sigma_2^2}}$$

S_1 = Mean of air

S_2 = Mean of Cu


σ_1 = Standard deviation of air

σ_2 = Standard deviation of Cu



Comparison of DR MTF's showed the Xradia high resolution mode has the best MTF



Au(Cu) edge



DR of edge



We are using MTF and SNR parameters of tubes as a measure of x-ray CT system performance





• Modulation Transfer Function (MTF)


- Frequency domain description of spatial resolution
- System MTF is the product of components
 - $MTF_{system} = MTF_{scintillator} * MTF_{optics} * MTF_{detector} * MTF_{staging}$

• Signal to Noise Ratio (SNR)


- Contrast measurement of images
- SNR factors
 - Source noise
 - Detector noise
 - Scattered radiation



Edge



Derivative



Fourier transform

$$SNR = \frac{S_1 - S_2}{\sqrt{\sigma_1^2 + \sigma_2^2}}$$



S_1 = Mean of air

σ_1 = Standard deviation of air


S_2 = Mean of Au

σ_2 = Standard deviation of Au


Three tubes are being used to measure x-ray CT system MTF's




LDPE




Copper
DR Transmission Images




Gold



LDPE




Copper
KCAT CT Images



Gold

Tube	Goodfellow specified outside diameter (mm)	Micrometer measured* outside diameter (mm)	CT measured outside diameter (mm)	Goodfellow specified wall thick (mm)	Mean CT measured wall thickness (mm)	Goodfellow Density (g/cm ³)
LDPE _Δ	1.100	1.075±0.007	1.08±0.01	0.350	0.389±0.009	0.92
Cu	0.750	0.625±0.007	0.625±0.002	0.220	0.197±0.008	8.96
Au	0.300	0.270±0.007 ⁺	0.28±0.02	0.050	0.042±0.002	19.32



Comparison of CT OD and ID showed the OD has a greater MTF than the ID due to scattered radiation




Copper

Gold

KCAT System



We need to quantitatively account for x-ray phase effects for accurate image analysis results




Phase effects can generate dimensional errors

Phase effects change with

- Object materials & geometry
- Source-object-detector geometry
- Source energy
- Spatial resolution



Phase effects can generate

- Dimensional errors
- Fictitious gaps
- Wrong material identification



Phase effects impact both radiographic and tomographic x-ray imaging

A transmission (I/I₀) radiograph of a Be strip with Lineout Illustrates x-ray phase effects




DR parameters: 60kV, 0.066 mA
100 sec/frame, 1 frame average
SDD=68.2mm, SOD=55.7mm
20X objective FOV: 1.2 mm X 1.2 mm
1X1 Binning; 0.582 μ m pixel pitch
WD Brown-7/18/2006- 11 UCRL-PRES-#

Image binned 4X4 for VG

UNCLASSIFIED


Schematic of gradient density reservoir





Side View



Top View




Photolithography is used to make graded density SU8





SEM image of SU8


DR image and lineout thru sample showing graded density with peak and valley anomalies



Xradia image
Image scaled to fit VG




Computed tomography image reveals a density gradient and artifacts from the phase effects



Artifacts make it difficult to quantify the density gradient




HADES simulations are being used to interpret the empirical data




The results reveal more work is needed

Both model and Xradia data normalized


A double shell is a complicated target with many critical dimensions 2004 design included 2 aerogel hemis



Attenuation radiograph of double shell targets revealed small voids between aerogel hemis in limited angles



Attenuation radiograph 0 degrees




Click on image to begin movie


Review of the manufacturing process reveals why the flaw may have occurred



μ XCT Attenuation images ($\ln I_0/I$) DS #7




Photograph of assembly DS #7




μ XCT data acquired on Dec 13, 2004
pixel size ~600 nm (1 x 1 binning)
Both I and I_0 images were 1 DR avg. 60 sec. ea.
sod 170 mm odd 10 mm

We use tomography to obtain full volumetric characterization






XY plane tomography image




Click image to start movie

We were able to accurately measure wall thickness



TARGET	GLOBAL ¹ [185°, 355°]				POLE ² [269°, 271°]	
	MEAN [μm]	STD DEV [μm]	MAX [μm]	MIN [μm]	MEAN [μm]	STD DEV ³ [μm]
1	51.14	0.47	52.20	49.52	51.22	0.00
2	49.59	0.35	50.44	48.46	49.73	0.25
3	51.56	0.61	52.52	49.68	51.52	0.28
6	50.37	0.46	51.63	48.95	49.94	0.18
7	54.39	0.46	55.86	53.24	53.86	0.00
8	49.92	0.53	51.22	48.38	50.07	0.06


As a result of our DR/CT, a new fabrication technique was developed to remove the wedge flaw by casting the inner sphere in aerogel



SiO₂ of
density
50 mg/cc

Capsule was centered using a vision system on the machine tool.



SiO₂ was machined concentric to the capsule.



The vision system also functions as a measuring microscope.



Volumetric rendering of double shell with new design confirmed the elimination of the wedge type flaws



Summary





- We have used x-ray DR and CT measurements to develop fabrication processes that meet target specifications
- We have also used x-ray images for metrology and quality control
- We are using DR and CT MTF's and SNR parameters to benchmark many different x-ray point projection systems
- Phase effects need to be considered in mesoscale imaging. We are developing object recovery algorithms to obtain an accurate object recovery

Backup Slides





Transmission image (I/I_0) of Double Shell showing inner SiO_2 and outer aerogel sphere with line-out



120 kVp, 0.033 mA
30 Sec Integration
20X FOV Setting (1.4 mm X 1.4 mm)
1 X 1 Binning, 1 Frame Avg
68.0 mm SDD, 55.7 mm SOD, pixel pitch 0.554 μm
WD Brown-7/18/2006- 26 UCRL-PRES-#

UNCLASSIFIED

Transmission* (I/I_0) radiograph of Be strip with laser drilled holes – Images cropped for VG





DR parameters: 60kV, 0.066 mA
100 sec/frame, 1 frame average
SDD=68.2mm, SOD=55.7mm
20X objective FOV: 1.2 mm X 1.2 mm
1X1 Binning; 0.582 μ m pixel pitch

WD Brown-7/18/2006- 27 UCRL-PRES-#



$$^*T=I/I_0=\exp-(\mu l)$$

UNCLASSIFIED

CT images (x-y plane) of hole #1 after volume rotation - images cropped for VG



LoTRT drawing and DR image (I/I₀) 90 degrees




90 kV, 0.067 mA
 120 Sec Integration
 5X FOV Setting (5.4 mm X 5.4 mm)
 1 X 1 Binning, 1 Frame Avg
 202 mm SDD, 197 mm SOD, pixel pitch 2.7µm
 WD Brown-7/18/2006- 29 UCRL-PRES-#


Transmission image (I/I₀)
 Image cropped and 4x4 binned

UNCLASSIFIED

We use tomography to obtain full volumetric characterization



YZ plane tomography image



Click image to start movie

Concentricity values for 6 Double Shell targets



Target Number	Avg Pixel Size (Microns)	Avg X Coord Outer Sphere	Avg Y Coord Outer Sphere	Avg Z Coord Outer Sphere	Avg X Coord Inner Sphere	Avg Y Coord Inner Sphere	Avg Z Coord Inner Sphere	RMS vector (Pixels)	RMS Vector (μmeters)
1	0.612	512.90	432.23	513.75	511.78	432.54	511.88	2.21	1.35
2	0.591	514.25	396.90	511.56	512.32	395.85	513.13	2.69	1.59
3	0.612	514.32	385.29	512.67	511.66	385.96	512.48	2.75	1.68
6	0.612	512.45	380.18	511.26	511.79	379.17	513.32	2.39	1.46
7	0.612	546.68	412.36	549.67	549.30	412.77	550.16	2.70	1.65
8	0.612	512.70	393.25	512.53	511.84	393.14	512.41	0.87	0.53

LoTRT DR Transmission Images are used to validate critical dimensions






Image cropped and 1X1 binned

90 kV, 0.067 ma
120 Sec Integration
5X FOV Setting
1 X 1 Binning, 1 Frame Avg
202 mm SDD, 197 mm SOD
WD Brown-7/18/2006- 32 UCRL-PRES-#

Line-out values




Derivative values



UNCLASSIFIED

We showed as-built target met specifications or were slightly out of spec




L1

L2

L3

L4

Derivative values



	L1 (mm)	L2 (mm)	L3 (mm)	L4 (mm)
Drawing Values	0.5 ± 0.003	0.3 ± 0.005	$0.1 - 0.15^* \pm 0.005$	0.05 ± 0.005
Measured Values	0.496	0.281	0.139	0.053

- Layered values are taken from Deg AAA05-500076AB
- L1 is not indicated on the drawing
- * L4 values vary in size 0.1, 0.125 and 0.15 mm

BrCH/SU8 gradient modeling

