

AFM and NanoSIMS analyses of Vaccinia virions

P. K. Weber, S. D. Gates, A. J. Malkin, R. C. Condit, N. Moussatche

March 4, 2013

Microscopy and Microanalysis Indianapolis, IN, United States August 4, 2013 through August 8, 2013

Disclaimer

This document was prepared as an account of work sponsored by an agency of the United States government. Neither the United States government nor Lawrence Livermore National Security, LLC, nor any of their employees makes any warranty, expressed or implied, or assumes any legal liability or responsibility for the accuracy, completeness, or usefulness of any information, apparatus, product, or process disclosed, or represents that its use would not infringe privately owned rights. Reference herein to any specific commercial product, process, or service by trade name, trademark, manufacturer, or otherwise does not necessarily constitute or imply its endorsement, recommendation, or favoring by the United States government or Lawrence Livermore National Security, LLC. The views and opinions of authors expressed herein do not necessarily state or reflect those of the United States government or Lawrence Livermore National Security, LLC, and shall not be used for advertising or product endorsement purposes.

AFM and NanoSIMS analyses of Vaccinia virions

Peter K. Weber¹, Sean D. Gates^{1,2}, Richard C. Condit³, Nissin Moussatche³ and Alexander J. Malkin¹

^{2.} Department of Mechanical Engineering, Stanford University, Stanford, California, USA.

Characterization of virus virion structure remains a topic of interest in virology. While electron microscopy and X-ray diffraction have provided a wealth of structural data on virions, questions remain regarding the structure and assembly of Pox viruses [1]. Here we explore the potential for high-resolution secondary ion mass spectrometry (SIMS) to be applied to the analysis of the virion structure. This work is technically of interest in the context of a mass spectrometry imaging session because the size of virions pushes the limits of resolution and sensitivity of imaging mass spectrometry. In our work, we are studying Vaccinia, a member of the Pox virus group that is used as the live vaccine for the small pox virus. Like other Pox viruses, the Vaccinia virion is relatively large for a virus (~350 x 300 x 260 nm hydrated, and ~330 x 260 x 125 nm dehydrated), making it an easier virus test subject for high-resolution SIMS. Our intention is to use dynamic SIMS with a monatomic primary ion beam, and therefore we are not seeking to detect molecules directly but rather intend to use isotopic (e.g., ¹⁵N) tags to map the distribution of labeled molecules within the virion.

Vaccinia virions in water were dispersed on single-crystal silicon wafers and air-dried for analysis. Samples were first imaged with atomic force microscopy (AFM) to provide virion dimensions and quality control, particularly to control for debris on the virions. For high-resolution SIMS, we used a Cameca NanoSIMS 50 with nominal lateral resolution of the Cs⁺ primary ion beam of 50 nm. The Cs⁺ beam was rastered over the sample to generate secondary ions, which were collected simultaneously on discrete electron multipliers to and used to generate ion-specific digital images of the sample. For this work, up to five species were collected simultaneously. Because the lateral resolution of the primary ion beam was at the scale of the virions, we used an imaging depth profile approach, which has the potential to provide higher spatial resolution in depth than could be achieved laterally [2]. This approach also avoided the need for sectioning. The challenge with this approach for the Vaccinia virions was that a substantial portion of the virion was consumed by sputtering before equilibrium could be achieved, and ion count rates constantly evolved during the course of analysis.

Co-registration between AFM and SIMS was achieved by imaging. The NanoSIMS secondary electron (e⁻) and ion imaging were able to resolve and image individual virions for positive identification of areas imaged by AFM (Fig. 1). Depth profiles for individual virions were extracted from serial NanoSIMS ion images for individual virions (Fig. 2). Our tests indicated that the CN⁻ species had the highest useful yield of the C-bearing species for iridium-coated virions, and the C₂⁻ species had the highest useful yield of the C-bearing species for uncoated virions. We also noted that yield of CN⁻ at the virion surface was enhanced by iridium coating.

^{1.} Physical and Life Sciences Directorate, Lawrence Livermore National Laboratory, Livermore, California, USA.

^{3.} Department of Molecular Genetics & Microbiology, University of Florida, Gainesville, Florida, USA.

Total ions collected for these species were 100,000 to 1,500,000. Based on these results, we anticipated that sufficient ion counts can be collected from isotopic labeling with ¹³C and ¹⁵N can provide information on the location of the labeled.

To test our ability to determine the location of specific molecules within the virion, we labeled Vaccinia virions with ¹⁵N-thymidine. Initial SIMS results showed that individual virions were labeled to the expected level of ¹⁵N enrichment. However the distribution of ¹⁵N within the virion was skewed relative to expectations. Follow up studies suggest that the rate of sputtering is changing rapidly while sputtering equilibrium is approached. We are currently developing a model of sputtering rate for use with our isotopic depth profiles. The latest results will be presented.

References:

- [1] RC Condit, N Moussatche, and P Traktman, Advances in Virus Research **66** (2006), pp. 31-124.
- [2] S Ghosal, SJ Fallon, TJ Leighton, KD Wheeler, MJ Kristo, ID Hutcheon, and PK Weber, Analytical Chemistry **80** (2008) pp. 5986-5992.
- [3] This work was funded by Lawrence Livermore National Laboratory (LLNL) Laboratory Directed Research and Development and performed under the auspices of the U.S. Department of Energy by LLNL under Contract DE-AC52-07NA27344. LLNL-CONF-624061




Figure 1. Correlated AFM and depth averaged NanoSIMS secondary electron (e⁻) and ¹²C¹⁴N⁻ images of vaccinia virions. The virions are consumed in a series of analysis scans.




Figure 2. Depth profile of single virions— uncoated (left) and iridium coated (right)— showing counts per cycle for four carbonbearing species.