



LAWRENCE
LIVERMORE
NATIONAL
LABORATORY

Report of the Third Biomedical Confidence Building Exercise

A. Alcaraz, A. Williams

April 30, 2013

Disclaimer

This document was prepared as an account of work sponsored by an agency of the United States government. Neither the United States government nor Lawrence Livermore National Security, LLC, nor any of their employees makes any warranty, expressed or implied, or assumes any legal liability or responsibility for the accuracy, completeness, or usefulness of any information, apparatus, product, or process disclosed, or represents that its use would not infringe privately owned rights. Reference herein to any specific commercial product, process, or service by trade name, trademark, manufacturer, or otherwise does not necessarily constitute or imply its endorsement, recommendation, or favoring by the United States government or Lawrence Livermore National Security, LLC. The views and opinions of authors expressed herein do not necessarily state or reflect those of the United States government or Lawrence Livermore National Security, LLC, and shall not be used for advertising or product endorsement purposes.

This work performed under the auspices of the U.S. Department of Energy by Lawrence Livermore National Laboratory under Contract DE-AC52-07NA27344.

Report of the Third Biomedical Confidence Building Exercise

Laboratory code: 07

Total number of pages: 77

TABLE OF CONTENTS













	Page no.
Cover page	1
Table of contents.....	2
Summary: Participating Laboratory	3
Summary: Quality System	4
Sample A	5
Sample B	9
Sample C	13
Sample D	29
Sample E.....	39
Sample F.....	49
Sample G	59

SUMMARY: PARTICIPATING LABORATORY

1. Participating laboratory

Name of the laboratory/institute:	Lawrence Livermore National Laboratory
Contact person:	Mr. Armando Alcaraz
Address:	7000 East Ave, L-091 Livermore, CA 94550 USA
Telephone number:	1-925-423-6889
Email address:	alcaraz1@llnl.gov
Date of sample receipt:	February 5 th , 2013
Date of report:	April 30 th , 2013

2. Analysts and authentication

	Name	Title	Signature
1	Armando Alcaraz	Principal Investigator	
2	Deon Anex	Research Scientist	
3	Bradley Hart	FSC Director	
4	Saphon Hok	Research Scientist	
5	Carolyn Koester	Research Scientist	
6	Roald Leif	Research Scientist	
7	Brian Mayer	Research Scientist	
8	Tuijuana Mitchell-Hall	QA Manager	
9	Heather Mulcahy	Research Scientist	
10	Carlos Valdez	Research Scientist	
11	Alexander Vu	Research Scientist	
12	Audrey Williams	Research Scientist	

SUMMARY: QUALITY SYSTEM

Quality system:

☒ Described in a Quality Assurance Manual/Handbook. Quality system in accordance with:

☒ ISO/IEC 17025

☐ Other/none: (describe below)

☒ Accreditation accepted; Year of accreditation: 2001

☐ Accreditation planned/pending; Year accreditation expected: _____

Accreditation body: A2LA

Scope of accreditation: CW analysis

Summary for laboratories with either “other” accreditation or no planned accreditation:

SAMPLE SUMMARY: A

Sample Code: P-301/07	Laboratory Assigned Code: CW-4-147-1
Description and condition of sample: Approximately 5 mL of plasma	

No chemicals relevant to the purpose of analysis were found.


SAMPLE PREPARATION DESCRIPTION: A

1. Sample preparation

Initial Aliquot Code	Type of Sample Preparation	Amount/ Volume	Sample/Blank Preparation Procedures	End Volume	Resulting Aliquot Code
CW-4-147-1	Fluoride Reactivation	1 mL	Added 3 mL of acetate buffer and 220 µL of 4.51 M KF. Incubated 30 minutes at room temperature. Neutralized with 0.8M Sodium Bicarbonate. ABS ELUT-Nexus cartridge cleanup. Eluted using 2 mL ethyl acetate and reduced sample volume to approximately 100 µL using nitrogen gas.	100 µL	CW-4-149-1-P-301
CW-4-149-1-P-301	Sample Split	50 µL	Split of sample CW-4-149-1-P-301.	50 µL	CW-4-149-3-P-301
CW-4-147-1	Fluoride Reactivation	2 mL	Added 3 mL of acetate buffer and 277 µL of 4.51 M KF. Incubated 30 minutes at room temperature. Neutralized with 0.8 M Sodium Bicarbonate. ABS ELUT-Nexus cartridge cleanup. Eluted using 2 mL ethyl acetate.	2 mL	CW-4-154-1-P-301
CW-4-147-1	Fluoride Reactivation	2 mL	Added 3mL of acetate buffer and 277 µL of 4.51 M KF. Incubated 30 minutes at room temperature. Neutralized with 0.8 M Sodium Bicarbonate. ABS ELUT-Nexus cartridge cleanup. Eluted using 2 mL ethyl acetate.	2 mL	CW-4-154-2-P-301
CW-4-154-1-P-301 CW-4-154-2-P-301	Concentrate combined samples	2 mLx2	Combined elution solvent from CW-4-154-1-P-301 and CW-4-154-2-P-301. Reduced sample volume to approximately 100 µL using nitrogen gas.	100 µL	CW-4-154-3-P-301

2. Additional information

--



Description of sample preparation and analysis methods

Sample preparation and analysis methods were developed using an in-house standard made by incubating the nerve agent, Sarin, with commercially procured human blood plasma. Cholinesterase activities of the plasma sample were measured using method outlined in Ellman (1961) before and after incubation to ensure agent-adduction occurred. The resultant agent-adducted plasma was used for subsequent method development.

Adducted plasma samples were reactivated using the process outlined in Degenhardt *et al* (2004). 1 mL of plasma was added to 3 mL of acetate buffer solution [0.189 M acetic acid; 10.8 mM sodium acetate]. Appropriate amount of Potassium Fluoride [4.51M] was added to achieve a final concentration of 0.25 M KF in the unknown sample work-up. Sample was allowed to incubate at room temperature for 30 minutes. 0.5 mL of 0.8 M sodium bicarbonate was added to neutralize acid.

Several methods were attempted to isolate the reactivated sarin from the sample:

1. Direct dichloromethane extraction
2. SPE preparation #1
Agilent ABS Elut-NEXUS (200 mg/ 6 mL)
Conditioning: 2 mL N-hexane; 4 mL Ethyl Acetate; 4 mL H₂O
Load 4.2 mL of reactivated sample.
Elution: 2 mL Ethyl Acetate
3. SPE preparation #2
Phenomenex Strata-X 33µm (30mg/ 3mL)
Conditioning: 2 mL N-hexane; 4 mL Ethyl Acetate; 4 mL H₂O
Load 4.2 mL of reactivated sample.
Elution: 2 mL Ethyl Acetate

The following procedure describes, in detail, the method that was determined to be most successful and used for the sample (SPE preparation #1):

An Agilent ABS Elut-Nexus SPE cartridge was conditioned, by gravity, with 2 mL of *N*-hexane, 2x2 mL ethyl acetate, followed by 2x2 mL of H₂O. 1 mL reactivated plasma sample was then loaded into cartridge and allowed to pass through. Sample was eluted with 2 mL of ethyl acetate into a clean vial. Sample was reduced in volume to approximately 100 µL for analysis.

When the target compound was not detected in CW-4-147-1, a deviation was made to this method. CW-4-147-1 was reactivated by using two replicates of 2 mL instead of the 1 mL described above and the resulting eluted samples were combined for their respective samples to increase the analytical concentration of the target analyte, should it be present. The combined samples were then reduced in volume to approximately 100 µL for analysis.

References

- G.L. Ellman. A new and rapid colorimetric determination of acetylcholinesterase activity. *Biochem. Pharm* 7: 88-95(1961).
- C. E. A. M. Degenhardt, K. Pleijsier, M. J. van der Schans, J. P. Langenberg, K. E. Preston, M. I. Solano, V. L. Maggio, J. R. Barr. Improvements of the Fluoride Reactivation Method for the verification of nerve agent exposure. *Journal of Analytical Toxicology*. 28(5): 364-371 (2004).

SAMPLE SUMMARY: B

Sample Code: P-302/07	Laboratory Assigned Code: CW-4-147-2
Description and condition of sample: Approximately 5 mL of plasma	

No chemicals relevant to the purpose of analysis were found.


SAMPLE PREPARATION DESCRIPTION: B

1. Sample preparation

Initial Aliquot Code	Type of Sample Preparation	Amount/ Volume	Sample/Blank Preparation Procedures	End Volume	Resulting Aliquot Code
CW-4-147-2	Fluoride Reactivation	1 mL	Added 3 mL of acetate buffer/0.25 M KF. Incubated 30 minutes at room temperature. Neutralized with 0.8M Sodium Bicarbonate. ABS ELUT-Nexus cartridge cleanup. Eluted using 2 mL ethyl acetate and reduced sample volume approximately 100 µL with nitrogen gas.	100 µL	CW-4-149-2-P-302
CW-4-149-2-P-302	Sample Split	100 µL	Split of sample CW-4-149-2-P-302.	50 µL	CW-4-149-4-P-302
CW-4-147-2	Fluoride Reactivation	2 mL	Added 3 mL of acetate buffer and 277 µL of 4.51 M KF. Incubated 30 minutes at room temperature. Neutralized with 0.8 M Sodium Bicarbonate. ABS ELUT-Nexus cartridge cleanup. Eluted using 2 mL ethyl acetate.	2 mL	CW-4-155-1-P-302
CW-4-147-2	Fluoride Reactivation	2 mL	Added 3 mL of acetate buffer and 277 µL of 4.51 M KF. Incubated 30 minutes at room temperature. Neutralized with 0.8 M Sodium Bicarbonate. ABS ELUT-Nexus cartridge cleanup. Eluted using 2 mL ethyl acetate.	2 mL	CW-4-155-2-P-302
CW-4-155-1-P-302 CW-4-155-2-P-302	Concentrate combined samples	2 mLx2	Combined elution solvent from CW-4-155-1-P-302 and CW-4-155-2-P-302. Reduced sample volume to approximately 100 µL using nitrogen gas.	100 µL	CW-4-155-3-P-302

2. Additional information

--



Description of sample preparation and analysis methods

Sample preparation and analysis methods were developed using an in-house standard made by incubating the nerve agent, Sarin, with commercially procured human blood plasma. Cholinesterase activities of the plasma sample were measured using method outlined in Ellman (1961) before and after incubation to ensure agent-adduction occurred. The resultant agent-adducted plasma was used for subsequent method development.

Adducted plasma samples were reactivated using the process outlined in Degenhardt *et al* (2004). 1 mL of plasma was added to 3 mL of acetate buffer solution [0.189 M acetic acid; 10.8 mM sodium acetate]. Appropriate amount of Potassium Fluoride [4.51M] was added to achieve a final concentration of 0.25 M KF in the unknown sample work-up. Sample was allowed to incubate at room temperature for 30 minutes. 0.5 mL of 0.8 M sodium bicarbonate was added to neutralize acid.

Several methods were attempted to isolate the reactivated sarin from the sample:

1. Direct dichloromethane extraction
2. SPE preparation #1
Agilent ABS Elut-NEXUS (200 mg/ 6 mL)
Conditioning: 2 mL N-hexane; 4 mL Ethyl Acetate; 4 mL H₂O
Load 4.2 mL of reactivated sample.
Elution: 2 mL Ethyl Acetate
3. SPE preparation #2
Phenomenex Strata-X 33µm (30mg/ 3mL)
Conditioning: 2 mL N-hexane; 4 mL Ethyl Acetate; 4 mL H₂O
Load 4.2 mL of reactivated sample.
Elution: 2 mL Ethyl Acetate

The following procedure describes, in detail, the method that was determined to be most successful and used for the sample (SPE preparation #1):

An Agilent ABS Elut-Nexus SPE cartridge was conditioned, by gravity, with 2 mL of *N*-hexane, 2x2 mL ethyl acetate, followed by 2x2 mL of H₂O. 1 mL reactivated plasma sample was then loaded into cartridge and allowed to pass through. Sample was eluted with 2 mL of ethyl acetate into a clean vial. Sample was reduced in volume to approximately 100 µL for analysis.

When the target compound was not detected in CW-4-147-2, a deviation was made to this method. CW-4-147-2 was reactivated by using two replicates of 2 mL instead of the 1 mL described above and the resulting eluted samples were combined for their respective samples to increase the analytical concentration of the target analyte, should it be present. The combined samples were then reduced in volume to approximately 100 µL for analysis.


References

- G.L. Ellman. A new and rapid colorimetric determination of acetylcholinesterase activity. *Biochem. Pharm* 7: 88-95(1961).
- C. E. A. M. Degenhardt, K. Pleijsier, M. J. van der Schans, J. P. Langenberg, K. E. Preston, M. I. Solano, V. L. Maggio, J. R. Barr. Improvements of the Fluoride Reactivation Method for the verification of nerve agent exposure. *Journal of Analytical Toxicology*. 28(5): 364-371 (2004).

SAMPLE SUMMARY: C

Sample Code: P-303/07	Laboratory Assigned Code: CW-4-147-3
Description and condition of sample: Approximately 5 mL of plasma	

Chemical: C-1


Chemical name & Structure		CAS #	Schedule
O-Isopropyl methylphosphonofluoridate		107-44-8	1.A.01
			
Aliquot(s)	Original/derivative	Analysis technique	
CW-4-150-4-P-303	Original	GC-MS/MS (CI)	
CW-4-150-4-P-303	Original	GC-MS (CI)	
CW-4-150-1-P-303	Original	GCxGC-MS (EI)	
Comments:			

SAMPLE PREPARATION DESCRIPTION: C**1. Sample preparation**

Initial Aliquot Code	Type of Sample Preparation	Amount/ Volume	Sample/Blank Preparation Procedures	End Volume	Resulting Aliquot Code
CW-4-147-3	Fluoride Reactivation	1 mL	Added 3 mL of acetate buffer and 220 µL of 4.51 M KF. Incubated 30 minutes at room temperature. Neutralized with 0.8 M Sodium Bicarbonate. ABS ELUT-Nexus cartridge cleanup. Eluted using 2 mL ethyl acetate and reduced sample volume to approximately 100µL with nitrogen gas.	100 µL	CW-4-150-1-P-303
CW-4-150-1-P-303	Sample Split	50 µL	Split of sample CW-4-150-1-P-303.	50 µL	CW-4-150-4-P-303

2. Additional information

--



Description of sample preparation and analysis methods

Sample preparation and analysis methods were developed using an in-house standard made by incubating the nerve agent, Sarin, with commercially procured human blood plasma. Cholinesterase activities of the plasma sample were measured using method outlined in Ellman (1961) before and after incubation to ensure agent-adduction occurred. The resultant agent-adducted plasma was used for subsequent method development.

Adducted plasma samples were reactivated using the process outlined in Degenhardt *et al* (2004). 1 mL of plasma was added to 3 mL of acetate buffer solution [0.189 M acetic acid; 10.8 mM sodium acetate]. Appropriate amount of Potassium Fluoride [4.51M] was added to achieve a final concentration of 0.25 M KF in the unknown sample work-up. Sample was allowed to incubate at room temperature for 30 minutes. 0.5 mL of 0.8 M sodium bicarbonate was added to neutralize acid.

Several methods were attempted to isolate the reactivated sarin from the sample:

1. Direct dichloromethane extraction
2. SPE preparation #1
Agilent ABS Elut-NEXUS (200 mg/ 6 mL)
Conditioning: 2 mL N-hexane; 4 mL Ethyl Acetate; 4 mL H₂O
Load 4.2mL of reactivated sample.
Elution: 2 mL Ethyl Acetate
3. SPE preparation #2
Phenomenex Strata-X 33µm (30mg/ 3mL)
Conditioning: 2 mL N-hexane; 4 mL Ethyl Acetate; 4 mL H₂O
Load 4.2mL of reactivated sample.
Elution: Ethyl Acetate

The following procedure describes, in detail, the method that was determined to be most successful and used for the sample (SPE preparation #1):

An Agilent ABS Elut-Nexus SPE cartridge was conditioned, by gravity, with 2 mL of N-hexane, 2x2 mL ethyl acetate, followed by 2x2 mL of H₂O. 1 mL reactivated plasma sample was then loaded into cartridge and allowed to pass through. Sample was eluted with 2 mL of ethyl acetate into a clean vial. Sample was reduced in volume to approximately 100 µL for analysis.

References

- G.L. Ellman. A new and rapid colorimetric determination of acetylcholinesterase activity. *Biochem. Pharm* 7: 88-95(1961).
- C. E. A. M. Degenhardt, K. Pleijsier, M. J. van der Schans, J. P. Langenberg, K. E. Preston, M. I. Solano, V. L. Maggio, J. R. Barr. Improvements of the Fluoride Reactivation Method for the verification of nerve agent exposure. *Journal of Analytical Toxicology*. 28(5): 364-371 (2004).

GC-MS/MS (CI) TECHNIQUE METHOD AND ANALYSIS DESCRIPTION

Identification

Chemical: C-1

Aliquot code: CW-4-150-4-P-303

Datafile name: TSQ0328

Compound identified as: Original Compound

Compound reference: Reference Chemical (OPCW)

Match algorithm and match factor: NIST, 980/999

Analysis Method

GC Instrument manufacturer and type: Thermo Fisher Trace GC Ultra

Carrier gas: Helium

Flow control/rate: Constant flow, 1 mL/min

Injection mode: Splitless, 0.60 min

Injector temperature: 230 °C

Column brand/phase: Agilent HP-5MS: (5%-Phenyl)-methylpolysiloxane

Column Length x ID x Film thickness: 30 m x 0.25 mm x 0.25 µm

GC temperature programme: 40 °C (3 min), 8 °C/min to 90 °C, 20 °C/min to 300 °C (1 min)

MS Instrument manufacturer and type: Thermo Fisher TSQ Quantum

Solvent delay time: 3 min

Electron energy: 100 eV

Reaction gas: Ammonia

Ionisation polarity: Positive

Scan range/time: 60-200 m/z in 1 second

Mass resolution: 0.7

Type of MS/MS scan: Product ion scan

Precursor ion(s): m/z 158

Collision gas: Argon

Collision Energy: 10

Ions/transitions	Sample		Standard		Criteria	Result
	Area*	Ion Ratio	Area*	Ion Ratio	Tolerance [#]	
99	450163	1.0	532327	1.0	N/A	N/A
141	62037	0.138	57494	0.108	±30%	27.6%

*Peak area of the ion, % intensity compared to the most abundant ion

[#] Compare with the relative ion intensity of the standard: >50%: ±20%; >20% to 50%: ±25%; >10% to 20%: ±30%; ≤10%: ±50%


Remarks

--

C:\OPCWBioMed2013\TSQ0327
CW-4-148-3-P-BLK

4/22/2013 7:35:34 PM
C:\Xcalibur\methods\BioMed 3rd\CW-PCI-NH3_MS-MS_GB.meth


RT: 4.16 - 8.16 SM: 9G



C:\OPCWBioMed2013\TSQ0328
CW-4-150-3-P-303

4/22/2013 8:10:02 PM
C:\Xcalibur\methods\BioMed 3rd\CW-PCI-NH3_MS-MS_GB.meth


RT: 4.18 - 8.18 SM: 9G



C:\OPCWBioMed2013\TSQ0331
CW-4-148-2-STD

4/22/2013 9:06:40 PM
C:\Xcalibur\methods\BioMed 3rd\CW-PCI-NH3_MS-MS_GB.meth

RT: 4.18 - 8.18 SM: 9G



GC-MS/MS (CI) chromatograms supporting identification of **Chemical C-1**; EIC


Top: Chromatogram of the blank.

Center: Chromatogram of the sample (aliquot code listed in header).

Bottom: Chromatogram of reference chemical of **O-Isopropyl methylphosphonofluoridate**.


C:\OPCWBioMed2013\TSQ0328 4/22/2013 8:10:02 PM
CW-4-150-4-P-303 C:\xcalibur\methods\BioMed 3rd\CW-PCI-NH3_MS-MS_GB.meth

TSQ0328 #185-186 RT: 6.15-6.16 AV: 2 SB: 6 5.98-6.06 NL: 1.63E5
T: + c EI Full ms2 158.000 [60.000-200.000]



C:\OPCWBioMed2013\TSQ0331 4/22/2013 9:06:40 PM
CW-4-148-2-STD C:\xcalibur\methods\BioMed 3rd\CW-PCI-NH3_MS-MS_GB.meth

TSQ0331 #185-186 RT: 6.14-6.16 AV: 2 SB: 6 5.98-6.06 NL: 7.52E4
T: + c EI Full ms2 158.000 [60.000-200.000]



CI product mass spectra of:

Top: **Chemical C-1** (aliquot code listed in header).

Bottom: Reference chemical of **O-Isopropyl methylphosphonofluoridate**.

GC-MS(CI) TECHNIQUE METHOD AND ANALYSIS DESCRIPTION

Identification

Chemical: C-1

Aliquot code: CW000626.D

Datafile name: CW-4-150-4-P-303

Compound identified as: Original Compound

Compound reference: Reference Chemical (OPCW)

Analysis Method

GC Instrument manufacturer and type: Agilent 6890

Carrier gas: Helium

Flow control/rate: Constant Flow, 32 cm/sec

Injection mode: Pulsed Splitless, 0.75 min

Injector temperature: 250 °C

Column brand/phase: Agilent HP-5MS: (5%-Phenyl)-methylpolysiloxane

Column Length x ID x Film thickness: 30 m x 0.25 mm x 0.25 µm

GC temperature programme: 40 °C (3 min), 8 °C/min to 100 °C, 25 °C/min to 300 °C (3 min)

MS Instrument manufacturer and type: Agilent 5975 MSD

Solvent delay time: 5.4 min

Electron energy: 82 eV

Scan range/time: SIM, m/z 141, 158, 100 µs dwell time

Source temperature: 230 °C

Reaction gas: Ammonia and **flow:** 21%

Ionisation polarity: Positive

Ions/transitions	Sample		Standard		Criteria	Result
	Area*	Ion Ratio	Area*	Ion Ratio	Tolerance [#]	
141	21504	0.038	28172	0.029	±50%	27.6%
158	571553	1.0	955561	1.0	N/A	N/A

*Peak area of the ion, % intensity compared to the most abundant ion



[#] Compare with the relative ion intensity of the standard: >50%: ±20%; >20% to 50%: ±25%; >10% to 20%: ±30%; ≤10%: ±50%

Remarks

Due to the solvent delay of 5.4 min, the chromatogram could not be presented with a window of the retention time of the compound minus 2 min.



File : C:\OPCWBioMed2013\CW000625.D
Sample : CW-4-148-3-P-BLK
Misc info: +NH3, CI, SIM, plasma method blank

Acquired: 22 Apr 2013 19:34
Method : CWNH3POS_SIM.M





File : C:\OPCWBioMed2013\CW000626.D
Sample : CW-4-150-4-P-303
Misc info: +NH3, CI, SIM, plasma

Acquired: 22 Apr 2013 20:04
Method : CWNH3POS_SIM.M



File : C:\OPCWBioMed2013\CW000609.D
Sample : CW-AV-1-128-1-A; with 50 ppb GB
Misc info: +NH3, SIM

Acquired: 19 Apr 2013 18:48
Method : CWNH3POS_SIM.M



EI chromatograms supporting identification of **Chemical C-1**;

TIC on left, EIC (m/z 158) on right.



Top: Chromatogram of the blank.

Center: Chromatogram of the sample (aliquot code listed in header).

Bottom: Chromatogram of reference chemical of **O-Isopropyl methylphosphonofluoridate**.



File : C:\OPCWBioMed2013\CW000625.D
Sample : CW-4-148-3-P-BLK
Misc info: +NH₃, CI, SIM, plasma method blank

Acquired: 22 Apr 2013 19:34
Method : CWNH3POS_SIM.M





File : C:\OPCWBioMed2013\CW000626.D
Sample : CW-4-150-4-P-303
Misc info: +NH₃, CI, SIM, plasma

Acquired: 22 Apr 2013 20:04
Method : CWNH3POS_SIM.M



File : C:\OPCWBioMed2013\CW000609.D
Sample : CW-AV-1-128-1-A; with 50 ppb GB
Misc info: +NH₃, SIM

Acquired: 19 Apr 2013 18:48
Method : CWNH3POS_SIM.M



EI chromatograms supporting identification of **Chemical C-1**;

TIC on left, EIC (m/z 141) on right.

Top: Chromatogram of the blank.

Center: Chromatogram of the sample (aliquot code listed in header).

Bottom: Chromatogram of reference chemical of **O-Isopropyl methylphosphonofluoridate**.

GCXGC-MS(EI) TECHNIQUE METHOD AND ANALYSIS DESCRIPTION

Identification

Chemical: C-1

Aliquot code: CW-4-150-1-P-303

Datafile name: TOF1335

Chemical identified as: Original Compound

Chemical reference: Reference Chemical

Match algorithm and match factor: 914/999

Analysis Method

GC Instrument manufacturer and type: Agilent 6890 GC

Carrier gas: Helium

Flow rate: 1.10 mL/min

☒ Corrected constant flow via pressure ramps

Injection mode: Splitless, 0.17 min

Injector temperature: 250° C

Primary Column:

Brand/phase: Agilent HP-5 MS UI: (5% diphenyl 95% dimethyl polysiloxane)

Column Length x ID x Film thickness: 30 m x 0.25 mm x 0.25 µm

GC temperature programme: 40°C (3 min), 8°C/min, 300°C (3 min)

Secondary Column:

Brand/phase: Restek Rxi-17: (50%-Phenyl)-methylpolysiloxane

Column Length x ID x Film thickness: 1.5 m x 0.18 mm x 0.18 µm

GC temperature programme: 50°C (3 min), 8°C/min, 310°C (3 min)

☒ Modulator enabled

Modulator temperature offset: 20 °C, relative to the GC secondary oven temperature

Modulation period: 3.5 sec

Hot pulse time: 0.6 sec

Cool time between stages: 1.15 sec

MS Instrument manufacturer and type: LECO Pegasus 4D

Solvent delay time: 3 minutes

Detector voltage: 1750 V

Electron energy: 70 eV

Mass resolution: 0.6 u

Scan range: 29-600 m/z

Scan rate: 200 spectra/sec

Source temperature: 250 °C

Ions/transitions	Sample		Standard		Criteria	Result
	Area*	Ion Ratio	Area*	Ion Ratio	Tolerance [#]	
39	3989.4	0.162	27388	0.152	±30%	6.6%
81	3972.7	0.162	27547	0.153	±30%	5.5%
99	24559	1.00	179664	1.00	N/A	N/A
125	4239.4	0.173	36495	0.203	±25%	15.0%

*Peak area of the ion, % intensity compared to the most abundant ion

[#] Compare with the relative ion intensity of the standard: >50%: ±20%; >20% to 50%: ±25%; >10% to 20%: ±30%; ≤10%: ±50%

Remarks


A small peak is present in the blank within 0.1 min of the target chemical. This peak has been tentatively identified as hexamethylcyclotrisiloxane through library matching (NIST, 960/999). The mass spectrum of this compound and the library spectrum are included.

File : C:\LLNL\TOF1334

Method: 2D OPCW020613

Acquired: 24 Apr 2013 1:10:02 PM

Sample: CW-4-148-1-P-BLK




File : C:\LLNL\TOF1335

Method: 2D OPCW020613

Acquired: 24 Apr 2013 1:52:49 PM

Sample: CW-4-150-1-P-303




File : C:\LLNL\TOF1336

Method: 2D OPCW020613

Acquired: 24 Apr 2013 4:34:02 PM

Sample: CW-4-150-3-SARINSTD

EI chromatograms supporting identification of **Chemical C-1**; EIC

Top: Chromatogram of the blank.

Center: Chromatogram of the sample (aliquot code listed in header).


Bottom: Chromatogram of reference chemical of **O-Isopropyl methylphosphonofluoridate**.

File : C:\LLNL\TOF1334

Method: 2D OPCW020613

Acquired: 24 Apr 2013 1:10:02 PM

Sample: CW-4-148-1-P-BLK




File : C:\LLNL\TOF1335

Method: 2D OPCW020613

Acquired: 24 Apr 2013 1:52:49 PM

Sample: CW-4-150-1-P-303




File : C:\LLNL\TOF1336

Method: 2D OPCW020613

Acquired: 24 Apr 2013 4:34:02 PM

Sample: CW-4-150-3-SARINSTD


GCxGC EI contour plot supporting identification of **Chemical C-1**; EIC

Top: Contour plot of the blank.


Center: Contour plot of the sample (aliquot code listed in header).

Bottom: Contour plot of reference chemical of **O-Isopropyl methylphosphonofluoridate**.

Peak True - sample "TOF:1335", peak 115, at 427 , 1.905 sec , sec



Peak True - sample "TOF:1336", peak 89, at 423.5 , 1.920 sec , sec




EI mass spectra of:


Top: **Chemical C-1** (file name listed in header).

Bottom: Reference chemical of **O-Isopropyl methylphosphonofluoridate**.

Peak True - sample "TOF:1334", peak 102, at 423.5 , 1.315 sec , sec



Library Hit - similarity 960, "Cyclotrisiloxane, hexamethyl"




EI mass spectra of:

Top: Background contaminant found in blank within ± 0.1 min of reported chemical (file name listed in header).

Bottom: Library mass spectrum of hexamethylcyclotrisiloxane.

SAMPLE SUMMARY: D**Sample Code:** P-304/07**Laboratory Assigned Code:** CW-4-147-4**Description and condition of sample:** Approximately 5 mL of plasma**Chemical:** D-1


Chemical name & Structure		CAS #	Schedule
O-Isopropyl methylphosphonofluoridate		107-44-8	1.A.01
			
Aliquot(s)	Original/derivative	Analysis technique	
CW-4-150-5-P-304	Original	GC-MS/MS (CI)	
CW-4-150-5-P-304	Original	GC-MS (CI)	
Comments:			

SAMPLE PREPARATION DESCRIPTION: D**1. Sample preparation**

Initial Aliquot Code	Type of Sample Preparation	Amount/ Volume	Sample/Blank Preparation Procedures	End Volume	Resulting Aliquot Code
CW-4-147-4	Fluoride Reactivation	1 mL	Added 3 mL of acetate buffer and 220 µL of 4.51 M KF. Incubated 30 minutes at room temperature. Neutralized with 0.8 M Sodium Bicarbonate. ABS ELUT-Nexus cartridge cleanup. Eluted using 2 mL ethyl acetate and reduced sample volume approximately 100 µL.	100 µL	CW-4-150-2-P-304
CW-4-150-2-P-304	Sample Split	50 µL	Split of sample CW-4-150-2-P-304.	50 µL	CW-4-150-5-P-304

2. Additional information

--



Description of sample preparation and analysis methods

Sample preparation and analysis methods were developed using an in-house standard made by incubating the nerve agent, Sarin, with commercially procured human blood plasma. Cholinesterase activities of the plasma sample were measured using method outline in Ellman (1961) before and after incubation to ensure agent-adduction occurred. The resultant agent-adducted plasma was used for subsequent method development.

Adducted plasma samples were reactivated using the process outlined in Degenhardt *et al* (2004). 1 mL of plasma was added to 3 mL of acetate buffer solution [0.189 M acetic acid; 10.8 mM sodium acetate]. Appropriate amount of Potassium Fluoride [4.51M] was added to achieve a final concentration of 0.25 M KF in the unknown sample work-up. Sample was allowed to incubate at room temperature for 30 minutes. 0.5 mL of 0.8 M sodium bicarbonate was added to neutralize acid.

Several methods were attempted to isolate the reactivated sarin from the sample:

4. Direct dichloromethane extraction
5. SPE preparation #1
Agilent ABS Elut-NEXUS (200 mg/ 6 mL)
Conditioning: 2 mL N-hexane; 4 mL Ethyl Acetate; 4 mL H₂O
Load 4.2 mL of reactivated sample.
Elution: 2 mL Ethyl Acetate
6. SPE preparation #2
Phenomenex Strata-X 33 μ m (30mg/ 3mL)
Conditioning: 2 mL N-hexane; 4 mL Ethyl Acetate; 4 mL H₂O
Load 4.2 mL of reactivated sample.
Elution: 2 mL Ethyl Acetate

The following procedure describes, in detail, the method that was determined to be most successful and used for the sample (SPE preparation #1):

An Agilent ABS Elut-Nexus SPE cartridge was conditioned, by gravity, with 2 mL of *N*-hexane, 2x2 mL ethyl acetate, followed by 2x2 mL of H₂O. 1 mL reactivated plasma sample was then loaded into cartridge and allowed to pass through. Sample was eluted with 2 mL of ethyl acetate into a clean vial. Sample was reduced in volume to approximately 100 μ L for analysis.

References

- G.L. Ellman. A new and rapid colorimetric determination of acetylcholinesterase activity. *Biochem. Pharm* 7: 88-95(1961).
- C. E. A. M. Degenhardt, K. Pleijsier, M. J. van der Schans, J. P. Langenberg, K. E. Preston, M. I. Solano, V. L. Maggio, J. R. Barr. Improvements of the Fluoride Reactivation Method for the verification of nerve agent exposure. *Journal of Analytical Toxicology*. 28(5): 364-371 (2004).

GC-MS/MS (CI) TECHNIQUE METHOD AND ANALYSIS DESCRIPTION

Identification

Chemical: D-1

Aliquot code: CW-4-150-5-P-304

Datafile name: TSQ0334

Compound identified as: Original Compound

Compound reference: Reference Chemical (OPCW)

Match algorithm and match factor: NIST 948/999

Analysis Method

GC Instrument manufacturer and type: Thermo Fisher Trace GC Ultra

Carrier gas: Helium

Flow control/rate: Constant flow, 1 mL/min

Injection mode: Splitless, 0.60 min

Injector temperature: 230 °C

Column brand/phase: Agilent HP-5MS: (5%-Phenyl)-methylpolysiloxane

Column Length x ID x Film thickness: 30 m x 0.25 mm x 0.25 µm

GC temperature programme: 40 °C (3 min), 8 °C/min to 90 °C, 20 °C/min to 300 °C (1 min)

MS Instrument manufacturer and type: Thermo Fisher TSQ Quantum

Solvent delay time: 3 min

Electron energy: 100 eV

Reaction gas: Ammonia

Ionisation polarity: Positive

Scan range/time: 60-200 m/z in 1 second

Mass resolution: 0.7

Type of MS/MS scan: Product ion scan

Precursor ion(s): m/z 158

Collision gas: Argon




Collision Energy: 10

Ions/transitions	Sample		Standard		Criteria	Result
	Area*	Ion Ratio	Area*	Ion Ratio	Tolerance [#]	
99	23240	1.0	532327	1.0	N/A	N/A
141	3087	0.133	57494	0.108	±30%	23.0%

*Peak area of the ion, % intensity compared to the most abundant ion

[#] Compare with the relative ion intensity of the standard: >50%: ±20%; >20% to 50%: ±25%; >10% to 20%: ±30%; ≤10%: ±50%

Remarks




GC-MS/MS (CI) chromatograms supporting identification of **Chemical D-1**; EIC

Top: Chromatogram of the blank.

Center: Chromatogram of the sample (aliquot code listed in header).

Bottom: Chromatogram of reference chemical of **O-Isopropyl methylphosphonofluoridate**.



CI product mass spectra of:

Top: **Chemical D-1** (aliquot code listed in header).

Bottom: Reference chemical of **O-Isopropyl methylphosphonofluoridate**.

GC-MS(CI) TECHNIQUE METHOD AND ANALYSIS DESCRIPTION

Identification

Chemical: D-1

Aliquot code: CW-4-150-5-P-304

Datafile name: CW000630.D

Compound identified as: Original Compound

Compound reference: Reference Chemical (OPCW)

Analysis Method

GC Instrument manufacturer and type: Agilent 6890

Carrier gas: Helium

Flow control/rate: Constant Flow, 32 cm/sec

Injection mode: Pulsed Splitless, 0.75 min

Injector temperature: 250 °C

Column brand/phase: Agilent HP-5MS: (5%-Phenyl)-methylpolysiloxane

Column Length x ID x Film thickness: 30 m x 0.25 mm x 0.25 µm

GC temperature programme: 40 °C (3 min), 8 °C/min to 100 °C, 25 °C/min to 300 °C (3 min)

MS Instrument manufacturer and type: Agilent 5975 MSD

Solvent delay time: 5.4 min

Electron energy: 82 eV

Scan range/time: SIM, m/z 141, 158, 100 µs dwell time

Source temperature: 230 °C

Reaction gas: Ammonia and flow: 21%

Ionisation polarity: Positive

Ions/transitions	Sample		Standard		Criteria	Result
	Area*	Ion Ratio	Area*	Ion Ratio	Tolerance [#]	
141	6306	0.044	28724	0.047	±50%	4.9%
158	142119	1.000	615664	1.000	N/A	N/A

*Peak area of the ion, % intensity compared to the most abundant ion



[#] Compare with the relative ion intensity of the standard: >50%: ±20%; >20% to 50%: ±25%; >10% to 20%: ±30%; ≤10%: ±50%

Remarks

Due to the solvent delay of 5.4 min, the chromatogram could not be presented with a window of the retention time of the compound minus 2 min.



File : C:\OPCWBioMed2013\CW000629.D
Sample : CW-4-148-3-P-BLK
Misc info: +NH3, CI, SIM, plasma method blank

Acquired: 22 Apr 2013 21:24
Method : CWNH3POS_SIM_MAN.M





File : C:\OPCWBioMed2013\CW000630.D
Sample : CW-4-150-5-P-304
Misc info: +NH3, CI, SIM, plasma (concentrated)

Acquired: 22 Apr 2013 21:55
Method : CWNH3POS_SIM_MAN.M



File : C:\OPCWBioMed2013\CW000633.D
Sample : CW-4-148-2-STD
Misc info: +NH3, CI, SIM, GB spiked in clean plasma

Acquired: 22 Apr 2013 23:13
Method : CWNH3POS_SIM_MAN.M



EI chromatograms supporting identification of **Chemical D-1**;

TIC on left, EIC (m/z 158) on right.



Top: Chromatogram of the blank.

Center: Chromatogram of the sample (aliquot code listed in header).

Bottom: Chromatogram of reference chemical of **O-Isopropyl methylphosphonofluoridate**.



File : C:\OPCWBioMed2013\CW000629.D
Sample : CW-4-148-3-P-BLK
Misc info: +NH3, CI, SIM, plasma method blank

Acquired: 22 Apr 2013 21:24
Method : CWNH3POS_SIM_MAN.M





File : C:\OPCWBioMed2013\CW000630.D
Sample : CW-4-150-5-P-304
Misc info: +NH3, CI, SIM, plasma (concentrated)

Acquired: 22 Apr 2013 21:55
Method : CWNH3POS_SIM_MAN.M



File : C:\OPCWBioMed2013\CW000633.D
Sample : CW-4-148-2-STD
Misc info: +NH3, CI, SIM, GB spiked in clean plasma

Acquired: 22 Apr 2013 23:13
Method : CWNH3POS_SIM_MAN.M



EI chromatograms supporting identification of **Chemical D-1**;


TIC on left, EIC (m/z 141) on right.

Top: Chromatogram of the blank.

Center: Chromatogram of the sample (aliquot code listed in header).

Bottom: Chromatogram of reference chemical of **O-Isopropyl methylphosphonofluoridate**.

SAMPLE SUMMARY: E**Sample Code:** U-305/07**Laboratory Assigned Code:** CW-4-147-5**Description and condition of sample:** Approximately 10 mL of urine**Chemical:** E-1

Chemical name & Structure		CAS #	Schedule
Isopropyl methylphosphonate		1832-54-8	2.B.04
			
Aliquot(s)	Original/derivative	Analysis technique	
CW-4-151-2-U-305	Pentafluorobenzyl derivative	GC-MS/MS(CI)	
CW-4-151-2-U-305	Pentafluorobenzyl derivative	GC-MS/MS(CI)	
Comments:			


SAMPLE PREPARATION DESCRIPTION: E

1. Sample preparation

Initial Aliquot Code	Type of Sample Preparation	Amount/ Volume	Sample/Blank Preparation Procedures	End Volume	Resulting Aliquot Code
CW-4-147-5	Urine Clean-Up and Pentafluorobenzyl Derivatization	3 mL	Silica SPE cartridge cleanup. Eluted using 3 mL 25% H ₂ O in acetonitrile. Reduced sample volume to complete dryness with nitrogen gas. Added 300 µL of acetone, 5 µL pentafluorobenzyl bromide and 5 µL diisopropylamine. Sample heated at 70°C for one hour.	310 µL	CW-4-151-1-U-305
CW-4-151-1-U-305	Sample Split	100 µL	Split of sample CW-4-151-1-U-305.	100 µL	CW-4-151-2-U-305
CW-4-151-1-U-305	Sample Split	100 µL	Split of sample CW-4-151-1-U-305.	100 µL	CW-4-151-3-U-305
CW-4-151-3-U-305	Preparation for LC-MS Analysis	100 µL	Reduced volume to ~25 µL with nitrogen gas. Added 75 µL of H ₂ O.	100 µL	CW-4-153-4-U-305

2. Additional information

SPE=solid phase extraction



Description of sample preparation and analysis methods

Sample preparation and analysis methods were developed using an in-house standard made by spiking a mixture of methylphosphonic acid, ethyl methylphosphonic acid, isopropyl methylphosphonic acid, and pinacolyl methylphosphonic acid into commercially procured human urine. The resultant acid-containing urine was used for subsequent method development.

Several SPE methods were attempted to isolate the phosphonic acids from the urine sample:

1. Alltech Silica (500 mg/4 mL)
Conditioning: 25% H₂O in acetonitrile; 3 mL of acetonitrile
Load 3x1 mL of urine sample.
Wash: 2mL acetonitrile; 2 mL 10% H₂O in acetonitrile.
Elution: 25% H₂O (in acetonitrile)
2. Agilent ABS Elut-NEXUS (200 mg/ 6 mL)
Conditioning: N-hexane; Ethyl Acetate; H₂O
Load 3x1 mL of urine sample
Elution: Ethyl Acetate
3. Phenomenex Strata-X 33µm (30 mg/ 3 mL)
Sample pretreatment: Dilute sample 1:1 with Acetate Buffer (pH ~3.5)
Conditioning: Methanol; Acetate Buffer (pH~3.5)
Load: 1 mL of urine sample
Wash: Acetate Buffer (pH ~3.5); methanol
Elution: 5% ammonium hydroxide in methanol

The following procedure describes, in detail, the method that was determined to be most successful and used for the sample (SPE preparation #1):

Spiked urine samples were processed following cleanup conditions found in Mawhinney (2007). An Alltech silica SPE was conditioned with 4mL 25% water in acetonitrile, followed by 3 mL of acetonitrile. Sample, 3x1 mL, was loaded onto cartridge and washed with 2 mL acetonitrile and 2 mL of 10% water in acetonitrile. Samples were eluted with using 25% water in acetonitrile.

Mawhinney took urine sample to dryness and reconstituted before introduction to SPE cartridge; however, in a side by side comparison, we found loading urine directly onto a conditioned cartridge resulted in the best recovery. Additionally, success has been reported with using polymeric SPE for cleanup and although we found success with one polymeric SPE, our best recovery was using silica SPE.

The urine sample eluted from the SPE cartridge was derivatized before analysis. The eluted sample was taken to dryness. Following the process outlined in Palit (2004), the dried samples were derivatized by adding 5 µL pentafluorobenzyl bromide, 5 µL diisopropylamine and 300 µL acetone and heating to 70°C for 1 hour.

References

- Mawhinney, D.B., Hameli, E.I., Fraser, R., Silva, S.S., Pavlopoulos, A.J. Kobelski, R.J. J.; The determination of organophosphate nerve agent metabolites in human urine by hydrophilic interaction liquid chromatography tandem mass spectrometry. *J. Chromatogr. B.* 852 (2007) 235-243.
- Palit, M., Gupta, A.K., Jain, R. and Raza, S.K.; Determination of pentafluorobenzyl derivatives of phosphonic and phosphonothionic acids by gas chromatography-mass spectrometry. *J Chromatogr. A*, 1043 (2004) 275-284.

GC-MS/MS (CI) TECHNIQUE METHOD AND ANALYSIS DESCRIPTION

Identification

Chemical: E-1
Aliquot code: CW-4-151-2-U-305
Datafile name: TSQ0358
Compound identified as: Pentafluorobenzyl derivative
Compound reference: Reference Chemical (own synthesis)
Match algorithm and match factor: NIST, 981/999

Analysis Method

GC Instrument manufacturer and type: Thermo Fisher Trace GC Ultra
Carrier gas: Helium
Flow control/rate: Constant flow, 1 mL/min
Injection mode: Splitless, 0.60 min
Injector temperature: 230 °C
Column brand/phase: Agilent HP-5MS: (5%-Phenyl)-methylpolysiloxane
Column Length x ID x Film thickness: 30 m x 0.25 mm x 0.25 µm
GC temperature programme: 40 °C (3 min), 8 °C/min to 215 °C, 20 °C/min to 300 °C (1 min)

MS Instrument manufacturer and type: Thermo Fisher TSQ Quantum
Solvent delay time: 3 min
Electron energy: 100 eV
Reaction gas: Methane
Ionisation polarity: Negative
Scan range/time: 50-500 m/z in 1 second
Mass resolution: 0.7
Type of MS/MS scan: Product ion scan
Precursor ion(s): m/z 137
Collision gas: Argon
Collision Energy: 10

Ions/transitions	Sample		Standard		Criteria	Result
	Area*	Ion Ratio	Area*	Ion Ratio	Tolerance [#]	
79	6684	0.104	243159	0.109	±30%	4.4%
95	61457	1.0	2338092	1.0	N/A	N/A
137	12236	0.149	347521	0.199	±30%	33.0%

*Peak area of the ion, % intensity compared to the most abundant ion


[#] Compare with the relative ion intensity of the standard: >50%: ±20%; >20% to 50%: ±25%; >10% to 20%: ±30%; ≤10%: ±50%

Remarks


--

C:\OPCWBioMed2013\TSQ0357
CW-4-153-2-U-BLK4/24/2013 9:51:07 PM
C:\Xcalibur\methods\BioMed 3rd\CW-NCI-CH4_MS-MS.meth


RT: 16.3 - 20.3

C:\OPCWBioMed2013\TSQ0363
CW-4-151-2-U-3054/24/2013 10:31:31 PM
C:\Xcalibur\methods\BioMed 3rd\CW-NCI-CH4_MS-MS.meth

RT: 16.3 - 20.3

C:\OPCWBioMed2013\TSQ0363
CW-AV-1-134-4-U-STD4/25/2013 1:53:31 AM
C:\Xcalibur\methods\BioMed 3rd\CW-NCI-CH4_MS-MS.meth

RT: 16.3 - 20.3



GC-MS/MS (CI) chromatograms supporting identification of **Chemical E-1 (pentafluorobenzyl derivative)**; TIC

Top: Chromatogram of the blank.


Center: Chromatogram of the sample (aliquot code listed in header).

Bottom: Chromatogram of reference chemical of **Pentafluorobenzyl derivative of isopropyl methylphosphonate**.

C:\OPCWBioMed2013\TSQ0358
CW-4-151-2-U-305

4/24/2013 10:31:31 PM
C:\Xcalibur\methods\BioMed 3rd\CW-NCI-CH4_MS-MS.meth


TSQ0358 #893-896 RT: 18.28-18.33 AV: 4 SB: 51 18.07-18.19, 18.46-19.19 NL: 1.50E4
T: - c EI Full ms2 137.000 [50.000-550.000]



C:\OPCWBioMed2013\TSQ0363
CW-AV-1-134-4-U-STD

4/25/2013 1:53:31 AM
C:\Xcalibur\methods\BioMed 3rd\CW-NCI-CH4_MS-MS.meth

TSQ0363 #893-896 RT: 18.28-18.33 AV: 4 SB: 51 18.07-18.19, 18.46-19.19 NL: 5.60E5
T: - c EI Full ms2 137.000 [50.000-550.000]



CI product mass spectra of:

Top: **Chemical E-1 (pentafluorobenzyl derivative)** (aliquot code listed in header).

Bottom: Reference chemical of **Pentafluorobenzyl derivative of isopropyl methylphosphonate**.

GC-MS/MS (CI) TECHNIQUE METHOD AND ANALYSIS DESCRIPTION

Identification

Chemical: E-1
Aliquot code: CW-4-151-2-U-305
Datafile name: TSQ0379
Compound identified as: Pentafluorobenzyl derivative
Compound reference: Reference Chemical (own synthesis)
Match algorithm and match factor: NIST, 848/999

Analysis Method

GC Instrument manufacturer and type: Thermo Fisher Trace GC Ultra
Carrier gas: Helium
Flow control/rate: Constant flow, 1 mL/min
Injection mode: Splitless, 0.60 min
Injector temperature: 230 °C
Column brand/phase: Agilent DB-1: 100% Dimethylpolysiloxane
Column Length x ID x Film thickness: 30 m x 0.25 mm x 0.25 µm
GC temperature programme: 40 °C (3 min), 8 °C/min to 215 °C, 20 °C/min to 300 °C (1 min)

MS Instrument manufacturer and type: Thermo Fisher TSQ Quantum
Solvent delay time: 3 min
Electron energy: 100 eV
Reaction gas: Methane
Ionisation polarity: Negative
Scan range/time: 50-500 m/z in 1 second
Mass resolution: 0.7
Type of MS/MS scan: Product ion scan
Precursor ion(s): m/z 137
Collision gas: Argon
Collision Energy: 10

Ions/transitions	Sample		Standard		Criteria	Result
	Area*	Ion Ratio	Area*	Ion Ratio	Tolerance [#]	
79	2700	0.022	211094	0.113	±30%	78.1%
95	120425	1.000	1862500	1.000	N/A	N/A
137	7176	0.109	256792	0.138	±30%	21.0%

*Peak area of the ion, % intensity compared to the most abundant ion

[#] Compare with the relative ion intensity of the standard: >50%: ±20%; >20% to 50%: ±25%; >10% to 20%: ±30%; ≤10%: ±50%


Remarks

--

C:\OPCWBioMed2013\TSQ0378
CW-4-153-2-U-BLK

4/26/2013 7:43:31 PM
C:\Xcalibur\methods\BioMed 3rd\CW-NCl-CH4_MS-MS_DB1.meth


RT: 15.9 - 19.9



C:\OPCWBioMed2013\TSQ0379
CW-4-153-2-U-305

4/26/2013 8:23:55 PM
C:\Xcalibur\methods\BioMed 3rd\CW-NCl-CH4_MS-MS_DB1.meth


RT: 15.9 - 19.9 SM: 9G



C:\OPCWBioMed2013\TSQ0384
CW-AV-1-134-4-U-STD

4/26/2013 11:45:55 PM
C:\Xcalibur\methods\BioMed 3rd\CW-NCl-CH4_MS-MS_DB1.meth

RT: 15.9 - 19.9



GC-MS/MS (CI) chromatograms supporting identification of **Chemical E-1 (pentafluorobenzyl derivative)**; TIC

Top: Chromatogram of the blank.


Center: Chromatogram of the sample (aliquot code listed in header).

Bottom: Chromatogram of reference chemical of **Pentafluorobenzyl derivative of isopropyl methylphosphonate**.

C:\OPCWBioMed2013\TSQ0379
CW-4-151-2-U-305

4/26/2013 8:23:55 PM
C:\Xcalibur\methods\BioMed 3rd\CW-NCI-CH4_MS-MS_DB1.meth


TSQ0379 #870-872 RT: 17.88-17.92 AV: 3 SB: 51 18.07-18.19, 18.46-19.18 NL: 1.76E4
T: - c EI Full ms2 137.000 [50.000-550.000]



C:\OPCWBioMed2013\TSQ0384
CW-AV-1-134-4-U-STD

4/26/2013 11:45:55 PM
C:\Xcalibur\methods\BioMed 3rd\CW-NCI-CH4_MS-MS_DB1.meth

TSQ0384 #869 RT: 17.87 AV: 1 SB: 4 17.82-17.87 NL: 5.07E4
T: - c EI Full ms2 137.000 [50.000-550.000]




CI product mass spectra of:

Top: **Chemical E-1 (pentafluorobenzyl derivative)** (aliquot code listed in header).

Bottom: Reference chemical of **Pentafluorobenzyl derivative of isopropyl methylphosphonate**.

SAMPLE SUMMARY: F**Sample Code:** U-306/07**Laboratory Assigned Code:** CW-4-147-6**Description and condition of sample:** Approximately 10 mL of urine**Chemical:** F-1


Chemical name & Structure		CAS #	Schedule
Isopropyl methylphosphonate		1832-54-8	2.B.04
			
Aliquot(s)	Original/derivative	Analysis technique	
CW-4-152-2-U-306	Pentafluorobenzyl derivative	GC-MS/MS(CI)	
CW-4-152-2-U-306	Pentafluorobenzyl derivative	GC-MS/MS(CI)	
Comments:			

SAMPLE PREPARATION DESCRIPTION: F**1. Sample preparation**

Initial Aliquot Code	Type of Sample Preparation	Amount/ Volume	Sample/Blank Preparation Procedures	End Volume	Resulting Aliquot Code
CW-4-147-6	Urine Clean-Up and Pentafluorobenzyl Derivatization	3 mL	Silica SPE cartridge cleanup. Eluted using 3 mL 25% H ₂ O in acetonitrile. Reduced sample volume to complete dryness with nitrogen gas. Added 300 µL of acetone, 5 µL pentafluorobenzyl bromide and 5 µL diisopropylamine. Sample heated at 70°C for one hour.	310 µL	CW-4-152-1-U-306
CW-4-152-1-U-306	Sample Split	100 µL	Split of sample CW-4-152-1-U-306.	100 µL	CW-4-152-2-U-306
CW-4-152-1-U-306	Sample Split	100 µL	Split of sample CW-4-152-1-U-306.	100 µL	CW-4-152-3-U-306
CW-4-152-3-U-306	Preparation for LC-MS Analysis	100 µL	Reduced volume to ~25 µL with nitrogen gas. Added 75 µL of H ₂ O.	100 µL	CW-4-153-5-U-306

2. Additional information

SPE=solid phase extraction



Description of sample preparation and analysis methods

Sample preparation and analysis methods were developed using an in-house standard made by spiking a mixture of methylphosphonic acid, ethyl methylphosphonic acid, isopropyl methylphosphonic acid, and pinacolyl methylphosphonic acid into commercially procured human urine. The resultant acid-containing urine was used for subsequent method development.

Several SPE methods were attempted to isolate the phosphonic acids from the urine sample:

1. Alltech Silica (500 mg/4 mL)
Conditioning: 25% H₂O in acetonitrile; 3 mL of acetonitrile
Load 3x1 mL of urine sample.
Wash: 2mL acetonitrile; 2 mL 10% H₂O in acetonitrile.
Elution: 25% H₂O (in acetonitrile)
2. Agilent ABS Elut-NEXUS (200 mg/ 6 mL)
Conditioning: N-hexane; Ethyl Acetate; H₂O
Load 3x1 mL of urine sample
Elution: Ethyl Acetate
3. Phenomenex Strata-X 33µm (30 mg/ 3 mL)
Sample pretreatment: Dilute sample 1:1 with Acetate Buffer (pH ~3.5)
Conditioning: Methanol; Acetate Buffer (pH~3.5)
Load: 1 mL of urine sample
Wash: Acetate Buffer (pH ~3.5); methanol
Elution: 5% ammonium hydroxide in methanol

The following procedure describes, in detail, the method that was determined to be most successful and used for the sample (SPE preparation #1):

Spiked urine samples were processed following cleanup conditions found in Mawhinney (2007). An Alltech silica SPE was conditioned with 4mL 25% water in acetonitrile, followed by 3 mL of acetonitrile. Sample, 3x1 mL, was loaded onto cartridge and washed with 2 mL acetonitrile and 2 mL of 10% water in acetonitrile. Samples were eluted with using 25% water in acetonitrile.

Mawhinney took urine sample to dryness and reconstituted before introduction to SPE cartridge; however, in a side by side comparison, we found loading urine directly onto a conditioned cartridge resulted in the best recovery. Additionally, success has been reported with using polymeric SPE for cleanup and although we found success with one polymeric SPE, our best recovery was using silica SPE.

The urine sample eluted from the SPE cartridge was derivatized before analysis. The eluted sample was taken to dryness. Following the process outlined in Palit (2004), the dried samples were derivatized by adding 5 µL pentafluorobenzyl bromide, 5 µL diisopropylamine and 300 µL acetone and heating to 70°C for 1 hour.

References

- Mawhinney, D.B., Hameli, E.I., Fraser, R., Silva, S.S., Pavlopoulos, A.J. Kobelski, R.J. J.; The determination of organophosphate nerve agent metabolites in human urine by hydrophilic interaction liquid chromatography tandem mass spectrometry. J. Chromatogr. B. 852 (2007) 235-243.
- Palit, M., Gupta, A.K., Jain, R. and Raza, S.K.; Determination of pentafluorobenzyl derivatives of phosphonic and phosphonothionic acids by gas chromatography-mass spectrometry. J Chromatogr. A, 1043 (2004) 275-284.

GC-MS/MS (CI) TECHNIQUE

METHOD AND ANALYSIS DESCRIPTION

Identification

Chemical: F-1

Aliquot code: CW-4-152-2-U-306

Datafile name: TSQ0360

Compound identified as: Pentafluorobenzyl derivative

Compound reference: Reference Chemical (own synthesis)

Match algorithm and match factor: NIST, 920/999

Analysis Method

GC Instrument manufacturer and type: Thermo Fisher Trace GC Ultra

Carrier gas: Helium

Flow control/rate: Constant flow, 1 mL/min

Injection mode: Splitless, 0.60 min

Injector temperature: 230 °C

Column brand/phase: Agilent HP-5MS: (5%-Phenyl)-methylpolysiloxane

Column Length x ID x Film thickness: 30 m x 0.25 mm x 0.25 µm

GC temperature programme: 40 °C (3 min), 8 °C/min to 215 °C, 20 °C/min to 300 °C (1 min)

MS Instrument manufacturer and type: Thermo Fisher TSQ Quantum

Solvent delay time: 3 min

Electron energy: 100 eV

Reaction gas: Methane

Ionisation polarity: Negative

Scan range/time: 50-500 m/z in 1 second

Mass resolution: 0.7

Type of MS/MS scan: Product ion scan

Precursor ion(s): m/z 137

Collision gas: Argon

Collision Energy: 10


Ions/transitions	Sample		Standard		Criteria	Result
	Area*	Ion Ratio	Area*	Ion Ratio	Tolerance [#]	
79	3424	0.025	243296	0.104	±30%	73.2%
95	139394	1	2338092	1.0	N/A	N/A
137	13329	0.096	347521	0.199	±30%	27.2%

*Peak area of the ion, % intensity compared to the most abundant ion

[#] Compare with the relative ion intensity of the standard: >50%: ±20%; >20% to 50%: ±25%; >10% to 20%: ±30%; ≤10%: ±50%

Remarks

--



GC-MS/MS (CI) chromatograms supporting identification of **Chemical F-1 (pentafluorobenzyl derivative)**; TIC

Top: Chromatogram of the blank.


Center: Chromatogram of the sample (aliquot code listed in header).

Bottom: Chromatogram of reference chemical of **Pentafluorobenzyl derivative of isopropyl methylphosphonate**.

C:\OPCWBioMed2013\TSQ0360
CW-4-152-2-U-306

4/24/2013 11:52:19 PM
C:\Xcalibur\methods\BioMed 3rd\CW-NCI-CH4_MS-MS.meth


TSQ0360 #893-896 RT: 18.28-18.33 AV: 4 SB: 51 18.07-18.19, 18.46-19.18 NL: 3.23E4
T: - c EI Full ms2 137.000 [50.000-550.000]



C:\OPCWBioMed2013\TSQ0363
CW-AV-1-134-4-U-STD

4/25/2013 1:53:31 AM
C:\Xcalibur\methods\BioMed 3rd\CW-NCI-CH4_MS-MS.meth

TSQ0363 #893-896 RT: 18.28-18.33 AV: 4 SB: 51 18.07-18.19, 18.46-19.19 NL: 5.60E5
T: - c EI Full ms2 137.000 [50.000-550.000]



CI product mass spectra of:

Top: **Chemical F-1(pentafluorobenzyl derivative)** (aliquot code listed in header).

Bottom: Reference chemical of **Pentafluorobenzyl derivative of isopropyl methylphosphonate**.

GC-MS/MS (CI) TECHNIQUE METHOD AND ANALYSIS DESCRIPTION

Identification

Chemical: F-1

Aliquot code: CW-4-152-2-U-306

Datafile name: TSQ0381

Compound identified as: Pentafluorobenzyl derivative

Compound reference: Reference Chemical (own synthesis)

Match algorithm and match factor: NIST, 820/999

Analysis Method

GC Instrument manufacturer and type: Thermo Fisher Trace GC Ultra

Carrier gas: Helium

Flow control/rate: Constant flow, 1 mL/min

Injection mode: Splitless, 0.60 min

Injector temperature: 230 °C

Column brand/phase: Agilent DB-1: 100% Dimethylpolysiloxane

Column Length x ID x Film thickness: 30 m x 0.25 mm x 0.25 µm

GC temperature programme: 40 °C (3 min), 8 °C/min to 215 °C, 20 °C/min to 300 °C (1 min)

MS Instrument manufacturer and type: Thermo Fisher TSQ Quantum

Solvent delay time: 3 min

Electron energy: 100 eV

Reaction gas: Methane

Ionisation polarity: Negative

Scan range/time: 50-500 m/z in 1 second

Mass resolution: 0.7

Type of MS/MS scan: Product ion scan

Precursor ion(s): m/z 137

Collision gas: Argon

Collision Energy: 10

Ions/transitions	Sample		Standard		Criteria	Result
	Area*	Ion Ratio	Area*	Ion Ratio	Tolerance [#]	
79	2700	0.022	211094	0.113	±30%	80.2%
95	120425	1.000	1862500	1.000	N/A	N/A
137	5865	0.049	256792	0.138	±30%	64.7%

*Peak area of the ion, % intensity compared to the most abundant ion


[#] Compare with the relative ion intensity of the standard: >50%: ±20%; >20% to 50%: ±25%; >10% to 20%: ±30%; ≤10%: ±50%

Remarks


--

C:\OPCWBioMed2013\TSQ0380
CW-4-152-2-U-BLK4/26/2013 9:04:19 PM
C:\Xcalibur\methods\BioMed 3rd\CW-NCI-CH4_MS-MS_DB1.meth


RT: 15.9 - 19.9 SM: 9G

C:\OPCWBioMed2013\TSQ0381
CW-4-152-2-U-3064/26/2013 9:44:43 PM
C:\Xcalibur\methods\BioMed 3rd\CW-NCI-CH4_MS-MS_DB1.meth

RT: 15.9 - 19.9 SM: 9G

C:\OPCWBioMed2013\TSQ0384
CW-AV-1-134-4-U-STD4/26/2013 11:45:55 PM
C:\Xcalibur\methods\BioMed 3rd\CW-NCI-CH4_MS-MS_DB1.meth

RT: 15.9 - 19.9



GC-MS/MS (CI) chromatograms supporting identification of **Chemical F-1 (pentafluorobenzyl derivative)**; TIC

Top: Chromatogram of the blank.


Center: Chromatogram of the sample (aliquot code listed in header).

Bottom: Chromatogram of reference chemical of **Pentafluorobenzyl derivative of isopropyl methylphosphonate**.

C:\OPCWBioMed2013\TSQ0381
CW-4-152-2-U-306

4/26/2013 9:44:43 PM
C:\Xcalibur\methods\BioMed 3rd\CW-NCI-CH4_MS-MS_DB1.meth


TSQ0381 #870-872 RT: 17.88-17.92 AV: 3 SB: 51 18.07-18.19, 18.46-19.18 NL: 2.46E4
T: - c EI Full ms2 137.000 [50.000-550.000]



C:\OPCWBioMed2013\TSQ0384
CW-AV-1-134-4-U-STD

4/26/2013 11:45:55 PM
C:\Xcalibur\methods\BioMed 3rd\CW-NCI-CH4_MS-MS_DB1.meth

TSQ0384 #869 RT: 17.87 AV: 1 SB: 4 17.82-17.87 NL: 5.07E4
T: - c EI Full ms2 137.000 [50.000-550.000]



CI product mass spectra of:


Top: **Chemical F-1(pentafluorobenzyl derivative)** (aliquot code listed in header).

Bottom: Reference chemical of **Pentafluorobenzyl derivative of isopropyl methylphosphonate**.


SAMPLE SUMMARY: G

Sample Code: U-307/07	Laboratory Assigned Code: CW-4-147-7
Description and condition of sample: Approximately 10 mL of urine	

Chemical: G-1

Chemical name & Structure		CAS #	Schedule
Ethyl methylphosphonate		1832-53-7	2.B.04
			
Aliquot(s)	Original/derivative	Analysis technique	
CW-4-152-5-U-307	Pentafluorobenzyl derivative	GC-MS/MS(CI)	
CW-4-152-5-U-307	Pentafluorobenzyl derivative	GC-MS/MS(CI)	
Comments:			

Chemical: G-2

Chemical name & Structure		CAS #	Schedule
Pinacolyl methylphosphonate		616-52-4	2.B.04
			
Aliquot(s)	Original/derivative	Analysis technique	
CW-4-152-5-U-307	Pentafluorobenzyl derivative	GC-MS/MS(CI)	
CW-4-152-5-U-307	Pentafluorobenzyl derivative	GC-MS/MS(CI)	
CW-4-153-6-U-307	Pentafluorobenzyl derivative	LC-MS/MS	
Comments:			


SAMPLE PREPARATION DESCRIPTION: G

1. Sample preparation

Initial Aliquot Code	Type of Sample Preparation	Amount/ Volume	Sample/Blank Preparation Procedures	End Volume	Resulting Aliquot Code
CW-4-147-7	Urine Clean-Up and Pentafluorobenzyl Derivatization	3 mL	Silica SPE cartridge cleanup. Eluted using 3 mL 25% H ₂ O in acetonitrile. Reduced sample volume to complete dryness with nitrogen gas. Added 300 µL of acetone, 5 µL pentafluorobenzyl bromide and 5 µL diisopropylamine. Sample heated at 70°C for one hour.	310 µL	CW-4-152-4-U-307
CW-4-152-4-U-307	Sample Split	100 µL	Split of sample CW-4-152-4-U-307.	100 µL	CW-4-152-5-U-307
CW-4-152-4-U-307	Sample Split	100 µL	Split of sample CW-4-152-4-U-307.	100 µL	CW-4-152-6-U-307
CW-4-152-6-U-307	Preparation for LC-MS Analysis	100 µL	Reduced volume to ~25 µL with nitrogen gas. Added 75 µL of H ₂ O.	100 µL	CW-4-153-6-U-307

2. Additional information

SPE=solid phase extraction



Description of sample preparation and analysis methods

Sample preparation and analysis methods were developed using an in-house standard made by spiking a mixture of methylphosphonic acid, ethyl methylphosphonic acid, isopropyl methylphosphonic acid, and pinacolyl methylphosphonic acid into commercially procured human urine. The resultant acid-containing urine was used for subsequent method development.

Several SPE methods were attempted to isolate the phosphonic acids from the urine sample:

1. Alltech Silica (500 mg/4 mL)
Conditioning: 25% H₂O in acetonitrile; 3 mL of acetonitrile
Load 3x1 mL of urine sample.
Wash: 2mL acetonitrile; 2 mL 10% H₂O in acetonitrile.
Elution: 25% H₂O (in acetonitrile)
2. Agilent ABS Elut-NEXUS (200 mg/ 6 mL)
Conditioning: N-hexane; Ethyl Acetate; H₂O
Load 3x1 mL of urine sample
Elution: Ethyl Acetate
3. Phenomenex Strata-X 33µm (30 mg/ 3 mL)
Sample pretreatment: Dilute sample 1:1 with Acetate Buffer (pH ~3.5)
Conditioning: Methanol; Acetate Buffer (pH~3.5)
Load: 1 mL of urine sample
Wash: Acetate Buffer (pH ~3.5); methanol
Elution: 5% ammonium hydroxide in methanol

The following procedure describes, in detail, the method that was determined to be most successful and used for the sample (SPE preparation #1):

Spiked urine samples were processed following cleanup conditions found in Mawhinney (2007). An Alltech silica SPE was conditioned with 4mL 25% water in acetonitrile, followed by 3 mL of acetonitrile. Sample, 3x1 mL, was loaded onto cartridge and washed with 2 mL acetonitrile and 2 mL of 10% water in acetonitrile. Samples were eluted with using 25% water in acetonitrile.

Mawhinney took urine sample to dryness and reconstituted before introduction to SPE cartridge; however, in a side by side comparison, we found loading urine directly onto a conditioned cartridge resulted in the best recovery. Additionally, success has been reported with using polymeric SPE for cleanup and although we found success with one polymeric SPE, our best recovery was using silica SPE.

The urine sample eluted from the SPE cartridge was derivatized before analysis. The eluted sample was taken to dryness. Following the process outlined in Palit (2004), the dried samples were derivatized by adding 5 µL pentafluorobenzyl bromide, 5 µL diisopropylamine and 300 µL acetone and heating to 70°C for 1 hour.

References

- Mawhinney, D.B., Hameli, E.I., Fraser, R., Silva, S.S., Pavlopoulos, A.J. Kobelski, R.J. J.; The determination of organophosphate nerve agent metabolites in human urine by hydrophilic interaction liquid chromatography tandem mass spectrometry. *J. Chromatogr. B.* 852 (2007) 235-243.
- Palit, M., Gupta, A.K., Jain, R. and Raza, S.K.; Determination of pentafluorobenzyl derivatives of phosphonic and phosphonothionic acids by gas chromatography-mass spectrometry. *J Chromatogr. A*, 1043 (2004) 275-284.

GC-MS/MS (CI) TECHNIQUE

METHOD AND ANALYSIS DESCRIPTION

Identification

Chemical: G-1

Aliquot code: CW-4-152-5-U-307

Datafile name: TSQ0362

Compound identified as: Pentafluorobenzyl derivative

Compound reference: Reference Chemical (own synthesis)

Match algorithm and match factor: NIST, 804/999

Analysis Method

GC Instrument manufacturer and type: Thermo Fisher Trace GC Ultra

Carrier gas: Helium

Flow control/rate: Constant flow, 1 mL/min

Injection mode: Splitless, 0.60 min

Injector temperature: 230 °C

Column brand/phase: Agilent HP-5MS: (5%-Phenyl)-methylpolysiloxane

Column Length x ID x Film thickness: 30 m x 0.25 mm x 0.25 µm

GC temperature programme: 40 °C (3 min), 8 °C/min to 215 °C, 20 °C/min to 300 °C (1 min)

MS Instrument manufacturer and type: Thermo Fisher TSQ Quantum

Solvent delay time: 3 min

Electron energy: 100 eV

Reaction gas: Methane

Ionisation polarity: Negative

Scan range/time: 50-500 m/z in 1 second

Mass resolution: 0.7

Type of MS/MS scan: Product ion scan

Precursor ion(s): m/z 123

Collision gas: Argon

Collision Energy: 10

Ions/transitions	Sample		Standard		Criteria	Result
	Area*	Ion Ratio	Area*	Ion Ratio	Tolerance [#]	
79	2358	0.049	307914	0.178	±30%	72.4%
95	48048	1	1731474	1	N/A	
123	445	0.009	145233	0.084	±50%	89.0%

*Peak area of the ion, % intensity compared to the most abundant ion


[#] Compare with the relative ion intensity of the standard: >50%: ±20%; >20% to 50%: ±25%; >10% to 20%: ±30%; ≤10%: ±50%

Remarks


--

C:\OPCWBioMed2013\TSQ0361
CW-4-153-2-U-BLK4/25/2013 12:32:43 AM
C:\Xcalibur\methods\BioMed 3rd\CW-NCI-CH4_MS-MS.meth


RT: 16.0 - 20.0 SM: 9G

NL:
5.00E1
m/z=
122.50-
123.50
MS
TSQ0361C:\OPCWBioMed2013\TSQ0362
CW-4-152-5-U-3074/25/2013 1:13:07 AM
C:\Xcalibur\methods\BioMed 3rd\CW-NCI-CH4_MS-MS.meth

RT: 16.0 - 20.0 SM: 9G

NL:
3.60E2
m/z=
122.50-
123.50
MS
Genesis
TSQ0362C:\OPCWBioMed2013\TSQ0363
CW-AV-1-134-4-U-STD4/25/2013 1:53:31 AM
C:\Xcalibur\methods\BioMed 3rd\CW-NCI-CH4_MS-MS.meth

RT: 16.0 - 20.0 SM: 9G

NL:
3.06E4
m/z=
122.50-
123.50
MS
TSQ0363

GC-MS/MS (CI) chromatograms supporting identification of **Chemical G-1 (pentafluorobenzyl derivative)**; TIC

Top: Chromatogram of the blank.


Center: Chromatogram of the sample (aliquot code listed in header).

Bottom: Chromatogram of reference chemical of **Pentafluorobenzyl derivative of ethyl methylphosphonate**.

C:\OPCWBioMed2013\TSQ0362
CW-4-152-5-U-307

4/25/2013 1:13:07 AM
C:\Xcalibur\methods\BioMed 3rd\CW-NCI-CH4_MS-MS.meth


TSQ0362 #876 RT: 17.98 AV: 1 SB: 9 17.90-17.95, 18.02-18.09 NL: 1.96E4
T: - c EI Full ms2 123.000 [50.000-550.000]



C:\OPCWBioMed2013\TSQ0363
CW-AV-1-134-4-U-STD

4/25/2013 1:53:31 AM
C:\Xcalibur\methods\BioMed 3rd\CW-NCI-CH4_MS-MS.meth

TSQ0363 #875 RT: 17.97 AV: 1 SB: 3 17.97-17.98, 17.95 NL: 4.53E5
T: - c EI Full ms2 123.000 [50.000-550.000]



CI product mass spectra of:

Top: **Chemical G-1(pentafluorobenzyl derivative)** (aliquot code listed in header).

Bottom: Reference chemical of **Pentafluorobenzyl derivative of ethyl methylphosphonate**.

GC-MS/MS (CI) TECHNIQUE METHOD AND ANALYSIS DESCRIPTION

Identification

Chemical: G-1

Aliquot code: CW-4-152-5-U-307

Datafile name: TSQ0383

Compound identified as: Pentafluorobenzyl derivative

Compound reference: Reference Chemical (own synthesis)

Match algorithm and match factor: NIST, 851/999

Analysis Method

GC Instrument manufacturer and type: Thermo Fisher Trace GC Ultra

Carrier gas: Helium

Flow control/rate: Constant flow, 1 mL/min

Injection mode: Splitless, 0.60 min

Injector temperature: 230 °C

Column brand/phase: Agilent DB-1: 100% Dimethylpolysiloxane

Column Length x ID x Film thickness: 30 m x 0.25 mm x 0.25 µm

GC temperature programme: 40 °C (3 min), 8 °C/min to 215 °C, 20 °C/min to 300 °C (1 min)

MS Instrument manufacturer and type: Thermo Fisher TSQ Quantum

Solvent delay time: 3 min

Electron energy: 100 eV

Reaction gas: Methane

Ionisation polarity: Negative

Scan range/time: 50-500 m/z in 1 second

Mass resolution: 0.7

Type of MS/MS scan: Product ion scan

Precursor ion(s): m/z 123

Collision gas: Argon




Collision Energy: 10

Ions/transitions	Sample		Standard		Criteria	Result
	Area*	Ion Ratio	Area*	Ion Ratio	Tolerance [#]	
79	2141	0.038	293708	0.182	±30%	79.4%
95	57029	1.000	1610485	1.000	N/A	N/A
123	2209	0.039	127164	0.079	±50%	40.9%

*Peak area of the ion, % intensity compared to the most abundant ion

[#] Compare with the relative ion intensity of the standard: >50%: ±20%; >20% to 50%: ±25%; >10% to 20%: ±30%; ≤10%: ±50%

Remarks



GC-MS/MS (CI) chromatograms supporting identification of **Chemical G-1 (pentafluorobenzyl derivative)**; TIC

Top: Chromatogram of the blank.


Center: Chromatogram of the sample (aliquot code listed in header).

Bottom: Chromatogram of reference chemical of **Pentafluorobenzyl derivative of ethyl methylphosphonate**.

C:\OPCWBioMed2013\TSQ0383
CW-4-152-5-U-307

4/26/2013 11:05:31 PM
C:\Xcalibur\methods\BioMed 3rd\CW-NCI-CH4_MS-MS_DB1.meth


TSQ0383 #848 RT: 17.50 AV: 1 SB: 40 20.57-20.93, 21.59-21.89 NL: 1.58E4
T: - c EI Full ms2 123.000 [50.000-550.000]



C:\OPCWBioMed2013\TSQ0384
CW-AV-1-134-4-U-STD

4/26/2013 11:45:55 PM
C:\Xcalibur\methods\BioMed 3rd\CW-NCI-CH4_MS-MS_DB1.meth

TSQ0384 #850 RT: 17.54 AV: 1 SB: 40 20.57-20.93, 21.60-21.89 NL: 4.26E4
T: - c EI Full ms2 123.000 [50.000-550.000]



CI product mass spectra of:

Top: **Chemical G-1(pentafluorobenzyl derivative)** (aliquot code listed in header).

Bottom: Reference chemical of **Pentafluorobenzyl derivative of ethyl methylphosphonate**.

GC-MS/MS (CI) TECHNIQUE METHOD AND ANALYSIS DESCRIPTION

Identification

Chemical: G-2

Aliquot code: CW-4-152-5-U-307

Datafile name: TSQ0362

Compound identified as: Pentafluorobenzyl derivative

Compound reference: Reference Chemical (own synthesis)

Match algorithm and match factor: NIST, 815/999

Analysis Method

GC Instrument manufacturer and type: Thermo Fisher Trace GC Ultra

Carrier gas: Helium

Flow control/rate: Constant flow, 1 mL/min

Injection mode: Splitless, 0.60 min

Injector temperature: 230 °C

Column brand/phase: Agilent HP-5MS: (5%-Phenyl)-methylpolysiloxane

Column Length x ID x Film thickness: 30 m x 0.25 mm x 0.25 µm

GC temperature programme: 40 °C (3 min), 8 °C/min to 215 °C, 20 °C/min to 300 °C (1 min)

MS Instrument manufacturer and type: Thermo Fisher TSQ Quantum

Solvent delay time: 3 min

Electron energy: 100 eV

Reaction gas: Methane

Ionisation polarity: Negative

Scan range/time: 50-500 m/z in 1 second

Mass resolution: 0.7

Type of MS/MS scan: Product ion scan

Precursor ion(s): m/z 179

Collision gas: Argon

Collision Energy: 10

Ions/transitions	Sample		Standard		Criteria	Result
	Area*	Ion Ratio	Area*	Ion Ratio	Tolerance [#]	
79	14388	0.175	128937	0.155	±30%	13.1%
95	82225	1.0	833281	1.0	N/A	N/A
179	5910	0.067	60141	0.072	±50%	7.0%

*Peak area of the ion, % intensity compared to the most abundant ion


[#] Compare with the relative ion intensity of the standard: >50%: ±20%; >20% to 50%: ±25%; >10% to 20%: ±30%; ≤10%: ±50%

Remarks


--

C:\OPCWBioMed2013\TSQ0361
CW-4-152-5-U-BLK4/25/2013 12:32:43 AM
C:\Xcalibur\methods\BioMed 3rd\CW-NCI-CH4_MS-MS.meth


RT: 19.5 - 23.4 SM: 9G

C:\OPCWBioMed2013\TSQ0362
CW-4-152-5-U-3074/25/2013 1:13:07 AM
C:\Xcalibur\methods\BioMed 3rd\CW-NCI-CH4_MS-MS.meth

RT: 19.5 - 23.5 SM: 9G

C:\OPCWBioMed2013\TSQ0363
CW-AV-1-134-4-U-STD4/25/2013 1:53:31 AM
C:\Xcalibur\methods\BioMed 3rd\CW-NCI-CH4_MS-MS.meth

RT: 19.5 - 23.5



GC-MS/MS (CI) chromatograms supporting identification of **Chemical G-2 (pentafluorobenzyl derivative)**; TIC

Top: Chromatogram of the blank.


Center: Chromatogram of the sample (aliquot code listed in header).

Bottom: Chromatogram of reference chemical of **Pentafluorobenzyl derivative of pinacolyl methylphosphonate**.

C:\OPCWBioMed2013\TSQ0362
CW-4-152-5-U-307

4/25/2013 1:13:07 AM
C:\Xcalibur\methods\BioMed 3rd\CW-NCI-CH4_MS-MS.meth


TSQ0362 #1078-1081 RT: 21.44-21.50 AV: 4 SB: 9 21.27-21.41 NL: 2.00E4
T: - c EI Full ms2 179.000 [50.000-550.000]



C:\OPCWBioMed2013\TSQ0363
CW-AV-1-134-4-U-STD

4/25/2013 1:53:31 AM
C:\Xcalibur\methods\BioMed 3rd\CW-NCI-CH4_MS-MS.meth

TSQ0363 #1078-1080 RT: 21.44-21.48 AV: 3 SB: 40 20.57-20.93, 21.60-21.89 NL: 2.60E5
T: - c EI Full ms2 179.000 [50.000-550.000]



CI product mass spectra of:

Top: **Chemical G-2(pentafluorobenzyl derivative)** (aliquot code listed in header).

Bottom: Reference chemical of **Pentafluorobenzyl derivative of pinacolyl methylphosphonate**.

GC-MS/MS (CI) TECHNIQUE METHOD AND ANALYSIS DESCRIPTION

Identification

Chemical: G-2

Aliquot code: CW-4-152-5-U-307

Datafile name: TSQ0383

Compound identified as: Pentafluorobenzyl derivative

Compound reference: Reference Chemical (own synthesis)

Match algorithm and match factor: NIST, 931/999

Analysis Method

GC Instrument manufacturer and type: Thermo Fisher Trace GC Ultra

Carrier gas: Helium

Flow control/rate: Constant flow, 1 mL/min

Injection mode: Splitless, 0.60 min

Injector temperature: 230 °C

Column brand/phase: Agilent DB-1: 100% Dimethylpolysiloxane

Column Length x ID x Film thickness: 30 m x 0.25 mm x 0.25 µm

GC temperature programme: 40 °C (3 min), 8 °C/min to 215 °C, 20 °C/min to 300 °C (1 min)

MS Instrument manufacturer and type: Thermo Fisher TSQ Quantum

Solvent delay time: 3 min

Electron energy: 100 eV

Reaction gas: Methane

Ionisation polarity: Negative

Scan range/time: 50-500 m/z in 1 second

Mass resolution: 0.7

Type of MS/MS scan: Product ion scan

Precursor ion(s): m/z 179

Collision gas: Argon

Collision Energy: 10

Ions/transitions	Sample		Standard		Criteria	Result
	Area*	Ion Ratio	Area*	Ion Ratio	Tolerance [#]	
79	6607	0.098	79743	0.153	±30	35.8%
95	67100	1.000	520164	1.000	N/A	N/A
179	3430	0.051	37110	0.071	±50%	28.3%

*Peak area of the ion, % intensity compared to the most abundant ion

[#] Compare with the relative ion intensity of the standard: >50%: ±20%; >20% to 50%: ±25%; >10% to 20%: ±30%; ≤10%: ±50%


Remarks

--

C:\OPCWBioMed2013\TSQ0382
CW-4-153-2-U-BLK

4/26/2013 10:25:07 PM
C:\Xcalibur\methods\BioMed 3rd\CW-NCl-CH4_MS-MS_DB1.meth


RT: 19.1 - 23.1



C:\OPCWBioMed2013\TSQ0384
CW-4-153-2-U-BLK

4/26/2013 11:05:31 PM
C:\Xcalibur\methods\BioMed 3rd\CW-NCl-CH4_MS-MS_DB1.meth


RT: 19.1 - 23.1



C:\OPCWBioMed2013\TSQ0384
CW-AV-1-134-4-U-STD

4/26/2013 11:45:55 PM
C:\Xcalibur\methods\BioMed 3rd\CW-NCl-CH4_MS-MS_DB1.meth

RT: 19.1 - 23.1



GC-MS/MS (CI) chromatograms supporting identification of **Chemical G-2 (pentafluorobenzyl derivative)**; EIC

Top: Chromatogram of the blank.


Center: Chromatogram of the sample (aliquot code listed in header).

Bottom: Chromatogram of reference chemical of **Pentafluorobenzyl derivative of pinacolyl methylphosphonate**.

C:\OPCWBioMed2013\TSQ0383
CW-4-152-5-U-307

4/26/2013 11:05:31 PM
C:\Xcalibur\methods\BioMed 3rd\CW-NCI-CH4_MS-MS_DB1.meth


TSQ0383 #1056-1058 RT: 21.07-21.10 AV: 3 SB: 10 21.75-21.90 NL: 8.12E3
T: - c EI Full ms2 179.000 [50.000-550.000]



C:\OPCWBioMed2013\TSQ0384
CW-AV-1-134-4-U-STD

4/26/2013 11:45:55 PM
C:\Xcalibur\methods\BioMed 3rd\CW-NCI-CH4_MS-MS_DB1.meth

TSQ0384 #1059-1060 RT: 21.12-21.14 AV: 2 SB: 10 21.75-21.91 NL: 3.10E4
T: - c EI Full ms2 179.000 [50.000-550.000]



CI product mass spectra of:

Top: **Chemical G-2(pentafluorobenzyl derivative)** (aliquot code listed in header).

Bottom: Reference chemical of **Pentafluorobenzyl derivative of pinacolyl methylphosphonate**.

LC-API-MS/MS TECHNIQUE METHOD AND ANALYSIS DESCRIPTION

Identification

Chemical: G-2**Aliquot code:** CW-4-153-6-U-307**Datafile name:** ORBI004531**Compound identified as:** Pentafluorobenzyl Derivative**Compound reference:** Reference Chemical (own synthesis)

Analysis Method


LC Instrument manufacturer and type: Thermo Surveyor Plus**Injection volume:** 10 μ L**Eluent composition:** A = H₂O with 0.1% formic acid
B = acetonitrile with 0.1% formic acid**Elution programme:** 95% A for 5 min, linear gradient to 20% A in 10 min, hold at 20% A for 10 min, regenerate column by returning to 95% A in 3 min, hold at 95% A for 10 min.**Flow rate:** 200 μ L/min**Column brand/phase:** Waters Atlantis T3/C18**Column Length x ID x Particle size:** 150 mm x 2.1 mm x 3 μ m**Column temperature:** 30 °C**MS Instrument manufacturer and type:** Thermo LTQ XL**Ionization type:** Electrospray**Ionisation polarity:** Positive**Electrospray/APCI voltage:** 4.0 kV**Type of MS/MS scan:** Product ion scan**Scan range:** m/z 75-750 u**Scan time:** 60 microsec/u**Collision energy:** 13 eV**Collision gas:** Helium**Mass resolution:** Unit**Precursor ion(s):** m/z 277

Ions/transitions	Sample		Standard		Criteria	Result
	Area*	Ion Ratio	Area*	Ion Ratio	Tolerance [#]	
181	9489	0.449	128565	0.313	±25%	43.7%
257	1179	0.056	14673	0.036	±50%	56.5%
277	21115	1.000	411161	1.000	N/A	N/A

*Peak area of the ion, % intensity compared to the most abundant ion

[#] Compare with the relative ion intensity of the standard: >50%: ±20%; >20% to 50%: ±25%; >10% to 20%: ±30%; ≤10%: ±50%


Remarks



C:\OPCWBioMed2013\ORBI004531
CW-4-153-6-U-307

4/25/2013 5:17:18 PM
C:\Xcalibur\methods\OPCW-MSMS-U307-T3-pos-022713-dsa.meth


ORBI004531 #3275-3285 RT: 19.96-20.08 AV: 11 SB: 14 21.75-21.91 NL: 1.41E3
T: ITMS + c ESI Full ms2 277.00@cid13.00 [75.00-750.00]



C:\OPCWBioMed2013\ORBI004533
B-DA-111-2

4/25/2013 6:35:20 PM
C:\Xcalibur\methods\OPCW-MSMS-U307-T3-pos-022713-dsa.meth

ORBI004533 #3144-3161 RT: 19.96-20.14 AV: 18 SB: 15 21.75-21.91 NL: 3.87E4
T: ITMS + c ESI Full ms2 277.00@cid13.00 [75.00-750.00]



ESI product mass spectra of:

Top: **Chemical G-2 (Pentafluorobenzyl derivative)** (aliquot code listed in header).

Bottom: Reference chemical of **Pentafluorobenzyl derivative of pinacolyl methylphosphonate**.