

The Optical Properties of Uranium Carbide and its Rate of Oxidation in Air Determined by Spectroscopic Ellipsometry

W. J. Siekhaus, S, B. Donald, C. Saw, A. Nelson

October 20, 2015

AVS 62 San Jose, CA, United States October 18, 2015 through October 23, 2015

Disclaimer

This document was prepared as an account of work sponsored by an agency of the United States government. Neither the United States government nor Lawrence Livermore National Security, LLC, nor any of their employees makes any warranty, expressed or implied, or assumes any legal liability or responsibility for the accuracy, completeness, or usefulness of any information, apparatus, product, or process disclosed, or represents that its use would not infringe privately owned rights. Reference herein to any specific commercial product, process, or service by trade name, trademark, manufacturer, or otherwise does not necessarily constitute or imply its endorsement, recommendation, or favoring by the United States government or Lawrence Livermore National Security, LLC. The views and opinions of authors expressed herein do not necessarily state or reflect those of the United States government or Lawrence Livermore National Security, LLC, and shall not be used for advertising or product endorsement purposes.

The optical properties of uranium carbide and its rate of oxidation in air determined by spectroscopic ellipsometry.

AVS 62nd International Symposium and Exhibition

18-23 October 2015

Wigbert J. Siekhaus, Scott B. Donald, Cheng K. Saw, Art Nelson







LLNL-Conf-678401

This work was performed under the auspices of the U.S. Department of Energy by Lawrence Livermore National Laboratory under contract DE-AC52-07NA27344. Lawrence Livermore National Security, LLC

Motivation


- Uranium carbide is used
 - as fuel in fast breeder reactors
 - as fuel (UC₂) in high temperature gas cooled reactors
- Uranium carbide
 - Appears as inclusions in uranium metal
 - Up to 300 weight ppm, 0.009 volume fraction






Uranium Carbide Oxidation

■ UC(s) +
$$2H_2$$
 O(g) \rightarrow UO₂(s) + C(s) + $2H_2$ (g) (1) UC(s) + $2H_2$ O(g) \rightarrow UO₂(s) + CH_4 (g) (2)




24 h at 20 mbar water vapour pressure (NO O_2) and 20° C. T.B. Scott et al. / Journal of Hazardous Materials 195 (2011) 115–123




• UC(s) +
$$2O_2(g) \rightarrow UO_2(s) + CO_2(g)$$
 (3)
2UC(s) + $3O_2(g) \rightarrow 2UO_2(s) + 2CO(g)$ (4)


UC Oxidation in UHV: XPS of O and U



Sputter-cleaned UC crystal



X-ray diffraction of UC single crystal




By scanning χ and ϕ , a two dimensional map of the (200) peak is acquired that indicates that

[100] planes are parallel to the sample's surface



Ellipsometric spectroscopy measures

the complex reflection coefficient ρ as f(photon energy).



 $\rho = \frac{R_P}{R_S} = \tan(\psi)e^{i\Delta}$

 R_{p} , R_{s} are the complex Fresnel

reflection coefficients

- Ellipsometry: measures Δ and Ψ with 1mm Ø beam as a function of photon energy at several angles of incidence
- Constants (n,k,ε₁,ε₂): are derived by fitting optical materials models (Cauchy, Urbach, Oscillator ..model)
- Quality of fit is determined by calculating the mean square error (MSE) value.

$$MSE = \sqrt{\frac{1}{2N - M} \sum_{n=1}^{N} \left[\left(\frac{\Psi_i^{\text{mod}} - \Psi_i^{\text{exp}}}{\sigma_{\Psi,i}^{\text{exp}}} \right)^2 + \left(\frac{\Delta_i^{\text{mod}} - \Delta_i^{\text{exp}}}{\sigma_{\Delta,i}^{\text{exp}}} \right)^2 \right]}$$

Exp. Data Measurement Gen. Data Model n,k Compare Fit Fit Parameters Thickness Results Roughness Uniformity

M = Number of fit parameter

Ellipsometry Work Outline




Remove by sputtering before measuring UC.
Grows as a function of air exposure time


Several mm thick. Thicker than the depth of penetration of 1-> 4 eV light

- Measure Δ,Ψ of a UO₂ crystal
 - Create optical model of UO₂ and a file of n,k (or ε₁, ε₂)
- Measure UC's Δ and Ψ after sputter-removing UO₂
 - Create optical model of UC and a file of n,k (or ε₁, ε₂)
- Measure Δ and Ψ changes as UO₂ grows on UC
 - Determine the thickness of UO₂ that minimizes MSE


<111> UO₂: Δ , Ψ & optical model

Siekhaus, Crowhurst OP Conf. Series: Mat. Sci. and Eng.9 (2010) 012055






UC's optical constants ε_1 , ε_2 & n, k.



Ellipsometry of UO₂ growth on UC




Growing UO_2 changes Δ and Ψ


UC, contribution to Δ and Ψ is fixed

Thickness grows as a function of exposure time to air


Fixed thickness and ellipsometry material model

- Measure ∆ and Ψ as UO2 grows on UC
 - "Fit" the thickness of UO₂ that minimizes UO₂ MŞE






 Δ and Ψ change due to 1 week air exposure



Oxide growth (~\time) on UC in air & fit quality




- The "Fit quality"
 = 1/MSE
 decreases with
 increasing UO2
 thickness
 - >>The "oxide" is NOT pure UO₂

$$MSE = \sqrt{\frac{1}{2N - M} \sum_{n=1}^{N} \left[\left(\frac{\Psi_i^{\text{mod}} - \Psi_i^{\text{exp}}}{\sigma_{\Psi,i}^{\text{exp}}} \right)^2 + \left(\frac{\Delta_i^{\text{mod}} - \Delta_i^{\text{exp}}}{\sigma_{\Delta,i}^{\text{exp}}} \right)^2 \right]}$$


$$M = \text{Number of fit parameter}$$

Model Fit and experimental Δ and Ψ at Final Spot 1 (~ 7 nm UO₂)



- Good fit to ∆ below 3.5 eV
- Bad fit to Ψ
 - Particularly the 'upward' movement of Ψ at 75° is completely missing
- >>The "oxide" is NOT pure UO₂

XPS: A graphitic carbon peak grows in/on the oxide of UC as it oxidizes in UHV



1) The optical properties of uranium carbide 2) its rate of oxidation in air determined by spectroscopic ellipsometry.

- 1)The optical constants (ε₁, ε₂) and (n, k) of a sputter-cleaned surface of single crystal UC have been determined by spectroscopic ellipsometry.
- 2) Changes in ∆ and Ψ during exposure to air of UC were measured as a function time.
- 2) A growth curve (~ √time) of oxide on UC was developed using an optical materials file of UO₂ to fit those changes in Δ and Ψ.
 - The quality of fit decreased with increasing oxide thickness, indicating that the oxide is most likely NOT pure UO₂
 - The 'fate' of carbon during oxidation is a possible reason.
 - XPS sees only the surface.
 - SIMS depth profiling will determine C distribution throughout the layer



